



AERO Z



CHOICE MATTERS AND THE NEW TILITE **AERO Z** DELIVERS.

TILITE INNOVATION.

- **Tru-Fit System:**
 - Reduces the number of holes in the frame, improving durability.
 - Allows for easy adjustment and reduces set-up time.
- **Adjustable Camber Mount** incorporates TiLite's Tru-Fit System with up to 3.5" of Rear Seat Height adjustment.
- **Speedloader** Adjustable Front Caster Housing (patent pending) quickly and accurately adjusts fork angle in 1° increments.
- **SlipStream** Single-Sided Fork with Integrated Bumper (patent pending).
- **Performance 5-Spoke Soft Roll Wheel** weighs just 1.6 ounces, 40% lighter than the original 3-Spoke Wheel.
- **TiShaft** back release bar makes folding and unfolding the backrest easy from any angle.

TECHNICAL SPECIFICATIONS.

Frame Type	Aluminum Mono-Tube Rigid - Full Adjustability
Frame Styles	SuperLite Edition Active Edition Rugged Edition Armed Forces Edition Power Adaptable Frame
Total Chair Width*	23"
Weight Limit	265 lbs.
Seat Widths	10" - 20"
Front Seat Widths (Seat Taper)	3 Options
Seat Depths	10" - 20"
Front Seat Heights	16" - 21"
Rear Seat Heights	12.5" - 20"
Front Frame Angles	75°, 80°, 85° or 90°
Footrest Width Options	Standard Taper or V-Front End
Cambers	0°, 2°, 4°, 6°, 8°, 12°, 15° or 18°
Front Wheels	22 Options
Rear Wheels	62 Options
Handrims	12 Options
Center of Gravity Adjustment	6"
Transport Weight*	12.4 lbs.
Convenience	Transport Bracket Option Folding Angle Adjustable Seat Back Depth Adjustable Seat Back
Aesthetics	6 Color Anodize Packages 19 Frame Colors 7 Ultrasuede Colors 8 Spinergy Wheel Colors
HCPCS Code	K0005
Warranty	Lifetime on Frame; Limited Lifetime on Components

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lbs.

CUSTOM CONFIGURATION THE TILITE WAY.

Custom configuration is when you are able to select options and components on a chair to meet your needs. While other companies offer some custom configuration, TiLite offers the widest range of options from which to choose. More choices mean more ability to truly create a custom configuration of options that best meets your needs. Choice matters and that is why TiLite offers more options than any other manufacturer! The choice is yours.



TiLite | 800.545.2266 | tilite.com
© 2011 TiSport, LLC.

All TiLite wheelchairs are shown with
available options. CF_AEROZS2_0911RevB