

Seat Colors  
Gray Vinyl

Body Color  
Red


Roll

New!

AlantéDX

The all-new Alanté DX features a fresh new compact base design with front-wheel drive technology to maneuver indoors and out. With a 300 lb. weight capacity and a standard Captain's seat, the Alanté DX offers end users an economical power base choice with many value-added features included.


GOLDEN  
technologies

YOUR LIFE IN MOTION

# AlantéDX

P O W E R C H A I R

<b>Model Number</b>	<b>GP204</b>
Medicare Code	K0823
Weight Capacity	300 lbs
Drive Wheels	Front
Maximum Speed	3.5 mph
Operating Range <sup>1</sup>	8 Miles
Ground Clearance to mid frame	2.5"
Turning Radius	24.5"
Type Batteries	2-22AH
Colors	Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	PG VR-2 50 Amp
<b>Dimensions:</b>	
Length	40"
Width	24"
Height (ground to top of back of seat)	34.5" - 36.5"
Height (ground to top of back of headrest)	41.5" - 43.5"
Ground to Top of Seat	20.5" - 22.5"
Footplate to Top of Seat	19"-21"
Number of Seat Height Positions	3
Size of Increments	1"
<b>Weight of Unit: (assembled)</b>	133 lbs
Base	67 lbs
Seat w/Arms	36 lbs
Batteries	30 lbs ea
<b>Tires:</b>	Flat Free
Front	9.5"
Casters	6"
Anti-Tippers Front	4"
<b>Standard Seat:</b>	Captain's Highback
Slide Seat	N/A
Back Height - no headrest	18"
Back Height with headrest	24.5" - 26.5"
Width x Depth	18" X 16"
Color	Gray Vinyl
<b>Optional Seat:</b>	Stadium Style
Back Height - no headrest	16"
Width x Depth	20" X 18"
Color	Charcoal Vinyl
<b>Warranty:</b>	
Frame	5 years
Drive Train	13 months
Electronics <sup>2</sup>	13 months


<sup>1</sup>. Battery range will vary due to rider weight, drive surface and terrain. Range listed tested at maximum weight.

<sup>2</sup>. Electronics warranty excludes batteries.

Disassembles easily  
for transport!



Disconnect the controller wire



Pull the captain's seat off the seat post.



Pull the battery pack out of the base

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Golden Technologies • 401 Bridge Street, Old Forge PA 18518 • Tel: (800) 624-6374 • Fax: (800) 628-5165 • [www.goldentech.com](http://www.goldentech.com)

SHT-GP204-011911