

Gemino selection guide

Choose the Gemino rollator that fits you 100%



GEMINO 20

Lightweight design rollator.
Great functionality at an attractive price



Indoor and outdoors use

Everyday use



331 lbs
287 lbs (M)
276 lbs (S)



4'11" - 6'7"
4'5" - 5'7" (M)
4'1" - 5'5" (S)
Also suitable for children



GEMINO 30

Handicare's 100% rollator.
Award-winning design combining high indoor maneuverability with outdoor performance. Lots of color options and accessories for easy customization



Use indoors and out, mainly in the city

Everyday use



331 lbs
287 lbs (M)
276 lbs (S)



4'11" - 6'7"
4'5" - 5'7" (M)
4'1" - 5'5" (S)
Also suitable for children



Small turning radius for great indoor maneuverability



Easy personalization with frame colors and bag colors. Wide range of accessories



GEMINO 60

Designed for any kind of terrain and heavy-duty use.
Lightweight, yet rock-solid.
All the benefits of the award-winning Gemino 30 and more



Mostly outdoor use, also on rugged terrain, occasionally indoors

Active / heavy-duty use



331 lbs
331 lbs (M)



4'11" - 6'7"
4'5" - 5'7" (M)




Extra wide and long frame for ultimate stability



Extra-large and soft wheels ensure ultimate comfort and easy negotiation of any obstacle



Extra wide seat



Extra strong bag (22 lbs load capacity)



GEMINO 30 walker

If extra support is needed.
Also great for gait training



Indoor use

Everyday use / balance training



331 lbs
287 lbs (M)



4'11" - 6'7"
4'5" - 5'7" (M)



Height-adjustable forearm supports for extra support and stability




handicare

www.handicare.com/us

Gemino range - specifications



		GEMINO 20			GEMINO 30			GEMINO 60		GEMINO 30 walker	
		20	20M	20S	30	30M	30S	60	60M	30	30M
DIMENSIONS											
Maximum user weight	lbs	330.7	286.6	275.6	330.7	286.6	275.6	330.7	330.7	330.7	286.6
Recommended user height	inch	4'11" - 6'7"	4'5" - 5'7"	4'1" - 5'5"	4'11" - 6'7"	4'5" - 5'7"	4'1" - 5'5"	4'11" - 6'7"	4'5" - 5'7"	4'11" - 6'7"	4'5" - 5'7"
Maximum load capacity in basket	lbs	11.0	11.0	11.0	11.0	11.0	11.0	22.0	22.0	11.0	11.0
Total length A	inch	25.6"	25.6"	26.8"	25.6"	25.6"	26.8"	29.1"	29.1"	25.6"	25.6"
Total width B	inch	23.6"	23.6"	21.7"	23.6"	23.6"	21.7"	25.2"	25.2"	24.0"	24.0"
Total weight without basket (with basket)	lbs	15.7 (16.1)	15.4 (15.9)	14.6 (15.0)	15.7 (16.1)	15.4 (15.9)	14.6 (15.0)	18.7 (19.4)	18.1 (18.7)	23.4 (10.9)	22.9 (10.7)
Width between center of push handles C	inch	18.5"	18.5"	16.5"	18.5"	18.5"	16.5"	20.1"	20.1"	10"	10"
Seat height D	inch	24.4"	21.7"	18.9"	24.4"	21.7"	18.9"	24.4"	20.9"	24.4"	21.7"
Total height E	inch	30.7 - 39.4	27.2 - 34.6	25.6 - 30.3	30.7 - 39.4	27.2 - 34.6	25.6 - 30.3	29.9 - 37.4	28.7 - 32.7	38.6 - 45.3	33.1 - 39.4
Total height folded	inch	31.5"	29.9"	27.2"	31.5"	29.9"	27.2"	33.1"	31.5"	38.2"	35.0"
Total length folded	inch	25.6"	25.6"	26.8"	25.6"	25.6"	26.8"	29.1"	29.1"	25.6"	25.6"
Total width folded	inch	9.1"	9.1"	9.1"	9.1"	9.1"	9.1"	9.1"	9.1"	14.0"	14.0"
Diameter front wheels F	inch	7.9" x 1.4"	7.9" x 1.4"	7.9" x 1.4"	7.9" x 1.4"	7.9" x 1.4"	7.9" x 1.4"	11 x 1.7	11 x 1.7	7.9" x 1.4"	7.9" x 1.4"
Diameter rear wheels G	inch	7.9" x 1.4"	7.9" x 1.4"	7.9" x 1.4"	7.9" x 1.4"	7.9" x 1.4"	7.9" x 1.4"	9.1 x 1.7	9.1 x 1.7	7.9" x 1.4"	7.9" x 1.4"
PERFORMANCE											
Turning radius (ISO 7176-5)	inch	33.1"	33.1"	31.3"	33.1"	33.1"	31.3"	37.4"	37.4"	33.1"	33.1"
COLOR OPTIONS											
Silver grey		√	√	√	√	√	√	√	√	√	√
Champagne					√	√					
Midnight blue					√	√	√				
Pink					√		√				
Bronze					√						
Black					√						
CERTIFICATES											
CE		√	√	√	√	√	√	√	√	√	√
EN 12182		√	√	√	√	√	√	√	√	√	√
ISO 11199-2		√	√	√	√	√	√	√	√		
ISO 11199-3										√	√



handicare

Handicare
2201 Hangar Place, Suite 200
Allentown, PA 18109
Tel: 610 266 5260
Toll Free: 866 276 5438
customerservice.usa@handicare.com
www.handicare.com/us