

Backpacker[®] SERIES

Backpacker[®] Plus Backpacker[®] MV Backpacker[®] AVP

Base weight lift	155 lbs. ¹	191 lbs. ¹	193 lbs. ¹
Mobility device overall weight	350 lbs.	350 lbs.	350 lbs.
Maximum product width	27"	27"	27"
Maximum travel distance from the ground	26"	37"	37"
Motor	Sealed 12 volt DC	Sealed 12 volt DC	Sealed 12 volt DC
Overall platform length	44" <small>(43" usable surface)</small>	44" <small>(43" usable surface)</small>	42" - 34" <small>(41" - 33" usable surface)</small>
Maximum product wheelbase	43"	43"	41" - 33"
Warranty	3-year transferable	3-year transferable	3-year transferable
Optional battery pack	12 volt, on-board	12 volt, on-board	12 volt, on-board
Optional charger	12 volt or 110 volt	12 volt or 110 volt	12 volt or 110 volt

Backpacker[®] SERIES



For the latest vehicle and product compatibility information, visit **"Find My Compatible Lift"** at prideliftsandramps.com or call for more information.

¹ Weight does not include any accessories

Battery Pack

NEW, OPTIONAL MULTI-PURPOSE BATTERY BOX



- Industry first
- Industry exclusive
- Charges the mobility device from the lift as you drive

PRIDE TOUGH

Backpacker[®] SERIES



**Backpacker[®]
PLUS**



**Backpacker[®]
MV**



**Backpacker[®]
AVP**



For more applications, see **"Find My Compatible Lift"** at prideliftsandramps.com.

**Pride[®]
Lifts & Ramps**

**PRIDE[®]
PERFORMANCE MOBILITY
— LIVE YOUR BEST[™] —**

Pride[®] Lifts & Ramps

www.prideliftsandramps.com • www.pridemobility.com
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

www.prideliftsandramps.com • www.pridemobility.com
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

Backpacker[®]

SERIES

Backpacker[®] Interior Lifts from Pride[™] Lifts and Ramps make it easier than ever to take your mobility product on the road. Backpacker Interior Lifts are the ultimate solution to lifting your powered mobility device into your vehicle, engineered specifically to be a perfect fit for a wide variety of popular vehicles including SUVs and minivans. Like all Pride Lifts, Backpacker Interior Lifts feature all-steel construction and a protective finish, along with an industrial-rated 12 volt motor and battery pack, and easy-to-use controls to combine rugged durability with safe, simple and reliable operation.

- Removable front drive panel for easier serviceability
- Adjustable wheel chocks on the platform surface
- Angled, low-profile platform
- Easy-to-use two-button hand control
- Manual override
- 3-year transferable warranty
- 350 lbs. lifting capacity



Charges the mobility device from the lift as you drive



One piece, fully assembled drive-on/drive-off platform



Angled, low-profile platform



Removable, front drive panel for easier serviceability



Adjustable wheel chocks on the platform surface



Easy-to-use two-button hand control