



Seat Color

Black Vinyl



Body Colors

Blue



Red



3 wheel GB-117 XLHD
4 wheel GB-147 XLHD

The Buzzaround XL HD features compact portability with a heavy-duty weight capacity! The Buzzaround XL HD features a robust 325 lb. weight capacity yet disassembles quickly and easily for transport. Standard features include a new LED angle-adjustable headlight, LED battery gauge, and more foot room than other scooters in its class. It's tight turning radius makes it ideal for indoor use, while it's ample ground clearance makes it great for outdoor use, as well.





LIGHTWEIGHT PORTABLE SCOOTER

Wireless disassembly in 4 easy steps!

Model	GB117 XLHD	GB147 XLHD
Medicare Code (pending)	K0801	K0801
Weight Capacity	325 lbs.	325 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	4 mph	4 mph
Operating Range up to ¹	14 miles	14 miles
Under Deck Clearance	3"	3"
Turning Radius	37"	47"
Type Batteries	2-22 AH	2-22 AH
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off-board	Off-board
Controller	50 Amp Dynamic R-Series	50 Amp Dynamic R-Series
Dimensions		
Length	42.5"	41.5"
Width	21.5"	21.5"
Height (ground to top of back of seat)	33"-35"	33"-35"
Ground to Top of Seat	21"-23"	21"-23"
Ground to Top of Deck	4.75"	4.75"
Deck to Top of Seat	16"-18"	16"-18"
Number of Seat Height Positions	3	3
Size of Increments	1"	1"
Front Axle to Rear Axle	33.5"	32"
Weight of Unit: (assembled)³		
Front Half	110 lbs.	118 lbs.
Rear Half	35 lbs.	43 lbs.
Seat w/Arms	25 lbs.	25 lbs.
Batteries	29 lbs pack	29 lbs pack
Tires:		
Front	Solid 8"	Solid 8"
Rear	Solid 8"	Solid 8"
Rear Anti-Tip Wheels	2"	2"
Standard Seat		
Seat Back Height - no headrest	Stadium Style 16"	Stadium Style 16"
Width x Depth	17" x 16"	17" x 16"
Color	Black Vinyl	Black Vinyl
Optional Seat:		
Back Height - no headrest	Stadium Style 16"	Stadium Style 16"
Width x Depth	18"x 16", 20"x 18"	18"x 16", 20"x 18"
Color	Black Vinyl	Black Vinyl
Warranty:		
Frame	Lifetime Warranty	Lifetime Warranty
Drive Train	2 Year Warranty	2 Year Warranty
Electronics ²	2 Year Warranty	2 Year Warranty
Batteries	1 Year Warranty	1 Year Warranty

- ¹ Battery range tested at 200 lbs. and will vary due to rider weight, drive surface and terrain.
² Electronics warranty excludes batteries.
³ Assembled weight calculated with the standard battery pack 2-22 AH batteries.



Step 1: Pop off the seat!



Step 2: Pull out the battery pack!



Step 3: Fold down tiller and lock in place!



Step 4: Pull forward on the drivetrain release lever to separate the front and rear frame sections!



Disassembles into 5 pieces!

4 wheel model shown



A battery charging port is conveniently located on the tiller under the control panel.

Features

- 325 lb. weight capacity
- Tiller mounted basket
- Angle-adjustable LED headlight
- 22 amp battery box

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Golden Technologies • 401 Bridge Street, Old Forge PA 18518 • Tel: (800) 624-6374 • Fax: (800) 628-5165 • www.goldentech.com

SHT-GB117XLHD-GB147XLHD-092115