



HCPCS
E1390


Pure Oxygen Concentrator

Item #	Product Description	UOM
CH5000S	Pure Oxygen Concentrator, Oxygen Sensing	1/ea
CH5000	Pure Oxygen Concentrator, Non-Oxygen Sensing	1/ea

- The Drive Pure is an easy to use, service, and maintain 5L concentrator that will have patients, referral sources and delivery technicians thanking you.
- Simple user interface with easy to view and operate power switch, LED indicator lights, and flow selector.
- Simple humidifier bottle connection within the compact, sleek, inconspicuous design.
- Available with both oxygen sensing and non-oxygen sensing configurations.
- Superior mobility with top handle, easy roll casters, and minimal weight.
- Easy to connect transfill port on oxygen sensing units.
- Simple single door access to both the foam and the HEPA filter.
- Exhaust system design directs airflow away from the user and reduces noise.

specifications	
Liter Flow Range	0.5 to 5 LPM
Dimensions	16.75" (W) x 10.25" (D) x 24.5" (H)
Weight	35 lbs.
Power	360 W at 3 LPM
Max Outlet Pressure	5.5 PSI +/- 0.5 PSI
Oxygen Purity	96% to 87%
Sound Level	45 dBA
Duty Cycle	Continuous Use
Operating Altitude	0 to 13,123 ft (4000 m) above sea level
Warranty	Limited 3 Year

accessories		UOM
HUM 001	Humidifier Bottle	50/cs
HUM-TUB	Humidifier Bottle Connection Adapter, 16"	50/cs
SOFT 225	Non-Kinking Cozy Nasal Cannula, 25'	25/cs
SOFT 250	Non-Kinking Cozy Nasal Cannula, 50'	20/cs
TUB NK 25	Non-Kinking Oxygen Tubing, 25'	25/cs
TUB NK 50	Non-Kinking Oxygen Tubing, 50'	20/cs
CON 550	Ridged Tubing Extension Connector	50/pkg
PL-2001844	Swivel Caster	1/ea
PL-3000706	Intake Filter Door	1/ea
PL-2001853	HEPA Intake Filter	1/ea
PL-2001854	Foam Intake Filter	1/ea