

Gait Trainers for patients who require ankle prompts to prevent scissoring.

Star Posterior Gait Trainer

SR 3000 Adult Star Posterior Gait Trainer, 1/cs

SR 3100 Pediatric Star Posterior Gait Trainer, 1/cs

SR 3200 Tyke Star Posterior Gait Trainer, 1/cs

Comet Anterior Gait Trainer

CO 2000 Adult Comet Anterior Gait Trainer, 1/cs

CO 2100 Pediatric Comet Anterior Gait Trainer, 1/cs

CO 2200 Tyke Comet Anterior Gait Trainer, 1/cs

FEATURES:

- Ankle Prompts – prevents legs from scissoring.
- Seat Harness – supports patients who cannot stand completely on their own.
- Pelvic Stabilizer – supports and stabilizes the pelvis.
- Positioning Bar – provides adjustable depth.
- Vertical handgrips – provides alternate handgrip position, posterior models only
- Variable Resistance Tabs – controls the speed of the gait trainer and aligns an asymmetrical gait pattern.
- One-Directional Rear Wheels – allows only forward mobility, posterior models only.
- Swivel Wheel Locking Bracket – provides option to lock swivel wheels.
- Limited Lifetime Warranty on the frame.

ACCESSORIES:

CE 1035 Set of adjustable forearm platforms, 1pr/bx

CE 1035 S Set of mounting brackets for forearm platform, 1pr/bx

CE 1060 Set of ankle prompts, 1pr/bx

CE 1080 Trunk support for anterior and posterior gait trainers, 1/bx

CE 1080 S Trunk support for pediatric & tyke anterior & posterior gate trainers, 1/bx

CE 1080 L Trunk support for adult anterior & posterior gait trainers, 1/bx

CE 1260 Set of extended uprights for the LT 3100, LT 3200, LT 2100 and LT 2200, 1pr/bx

CE 1270 Set of extended uprights for the LT 3000 and LT 2000, 1pr/bx



Posterior



Anterior

ITEM #	ADJUSTABLE HEIGHT	INSIDE WIDTH BETWEEN HANDLEBAR	WIDTH
SR 3000	29" - 36"	21.5"	26.5"
SR 3100	22" - 28"	16"	23.5"
SR 3200	17" - 20"	14.5"	22"
CO 2000	29" - 36"	21.5"	26.5"
CO 2100	22" - 28"	16"	23.5"
CO 2200	17" - 20"	14.5"	22"



mobility
gait trainers

Wenzel Re/hab

To order call toll free:
New York: 877.224.0946 • California: 800.371.2266

drive
MEDICAL DESIGN & MANUFACTURING