



Yes, you can.®

CS9™ FX600 Adjustable-Width Bed

Bed featured with ThinkSoft® Positioning
Device, Kendall bed ends in Biltmore
Cherry and Solace® mattress



The innovative CS9 FX600 adjusts from 36" to 39" or 42" to accommodate a variety of residents, offering facilities the utmost in flexibility. The frame can easily be adjusted without any tools and supports up to 600 pounds for a truly versatile long-term care bed.

CS9™ FX600



Expandable deck
Easily adjusts
without tools to
meet your needs



Expandable headboard
Helps ensure resident
safety at all widths



Optional under-bed light
Helps the resident
to safely get out of
bed at night



Hand Control Pendant
A rugged, water-resistant
pendant with fewer
recessed areas is easy to
clean and very intuitive for
the resident to use

CS9 FX600 Adjustable-Width Bed

- Deck quickly adjusts from 36" to 39" or 42" without additional parts or tools
- Optional Safety Package includes under-bed lighting and integrated nurse call, all through the hand pendant
- Includes tool-less adjustable headboard and a fixed 39" footboard

Accessories

- Trapeze, 200 lbs.
- Trapeze, 300 lbs.
- Battery Backup
- Optional Trendelenburg/
Reverse Trendelenburg
- 4" Bed Extender
- Optional Attendant Control Panel

Specifications

- 9"– 30" Height Range
- 600 lbs. Weight Capacity
- Auto Contour
- Ergonomically Designed, Rugged,
Fluid Resistant 8 Button Pendant
- Mattress Keepers

Invacare Continuing Care, Inc.

www.invacare-cc.com

info_icc@invacare.com

USA

One Invacare Way
Elyria, Ohio 44035
USA
800-668-2337

Canada

994 Hargrieve Road
London, Ontario
N6E 1P5 Canada
800-668-2337

Trademarks are identified by the symbol™ and ®.
All trademarks are owned or licensed to Invacare
Corporation or its subsidiaries unless otherwise noted.

© 2010 Invacare Corporation. All rights reserved.
Form 08.70.1/10



Yes, you can.®