



# THE COMFORT COMPANY

*Quietly making lives more comfortable®*

## **PRODUCT CARE INSTRUCTIONS** (COVER WASHING INSTRUCTIONS)

---

1. All Comfort Company covers, should be washed on cold or warm cycle using any common laundry detergents. Once cover has been washed, please hang dry damp cover.
2. If the cover cannot be removed: Spray and wipe with any household or commercial disinfectant, cleaner, or anti bacterial spray. It is recommended that any cleaners, or disinfectant be wiped away with a "water only" damp cloth afterwards.