


The Cozzia MC-520



Cranberry



Briarwood



Saddle

The MC-520 is the only lift chair in the world that helps you get out, *but also helps you get in.*

The MC-520 features a lift function, retractable arms, an adjustable headrest, air massage, and heat.

cozzia®

The Cozzia MC-520

Model #	MC-520
Chair Upright Dimensions (LxWxH)	37.4" x 29.9" x 59.4"
Chair Full Recline Dimensions (LxWxH)	73.2" x 29.9" x 22.8"
Weight Capacity	300 LBS
Headrest Adjustable Height	3"
Seat Width Between Arms	21"
Seat Height	20"
Seat Depth	20"
Back Height from Seat	28"
Distance from Wall	15"
Recline Position #	Infinite Recline Positions
Fabric Type	Bonded Fabric
Trendelenburg Position	✓
Memory Foam Seating	✓
Heat Therapy	✓
Air Massage	✓
Air Lumbar Support	✓
Retractable Arm	✓
Adjustable Headrest	✓
Footrest Extension	✓
Casters for Easy Moving	✓
Battery Backup	✓

Easy in. Easy out.

The MC-520 has an electric lift function to give you easy access to the chair. Both armrests of the MC-520 are power adjustable to make getting in to and out of the chair seamless.

Relax with heat and massage.

The MC-520 has many features built into it that make it our most advanced lift chair. These include: Stretch air massage, Relax air massage, and heating.

Adjustable headrest.

The headrest of the MC-520 is adjustable in two ways. You can alter the headrest's height, and you're able to move the wings of the headrest in and out.

WWW.COZZIAUSA.COM

cozzia®

TEL: 887-977-0656 / FAX: 800-521-4712

861 S. OAK PARK ROAD, COVINA, CA 91724