



Invacare® Top End® Crossfire™ All Terrain Wheelchair

The **Invacare® Top End® Crossfire™ All Terrain Wheelchair** is the monster truck of wheelchairs and is designed to leave some serious tracks. Do you like to hike, fish, camp, hunt, kayak, garden or go to the beach? You can do it all with the Crossfire All Terrain Wheelchair. Built of 6061 T6 Aluminum, it is equipped with enormous 25" x 2.75" knobby wheels and huge 8" x 2" front casters for extreme traction. The Crossfire All Terrain Wheelchair is engineered with mega ground clearance, plus forward and side-to-side stability so that you can conquer the toughest terrain.

Features

- 6061 T6 Aluminum frame
- Adjustable center-of-gravity
- Camber 0°, 3°, 6°, 9°
- Lock down folding back
- Adjustable back height
- Adjustable tension seat upholstery
- Adjustable tension back upholstery
- 24" or 25" high performance Sun® wheels
- Anodized handrims
- 24" or 25" x 2.75" knobby tires
- 8" x 2" flat free front casters
- Adjustable height and angle footplate
- Powder coat paint
- 250 lb. weight capacity

Model #: CRFAT



1. 24" or 25" x 2.75" knobby wheels and tires*
2. 8" x 2" flat free front casters
3. Engineered for front and side-to-side stability
4. Center-of-gravity adjustment
5. Anti-tip option
6. Optional push handles, wheel locks and anti-tip bars

*24" tires are 24" x 2 gray, 25" tires are 25 x 2.5 black



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(800) 532-8677
Fax: (727) 522-1007

©2018

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols™ and®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted. Natural
Fit is a registered
trademark of Three Rivers
Out Front. Spinerger is
a registered trademark
of Spinerger, Inc.
Specifications are subject
to change without
notification.
Form No. 18-321
181123

Options

- Spinerger® Spox wheels (narrow)
- Quad release axles
- Coated handrims
- Natural Fit® handrims
- Fixed back
- Push handles
- Carbon fiber sideguards
- Removeable sideguards
- Folding sideguards
- Wheel locks
- Positioning straps
- Alignment gauge
- Footplate risers
- 2" cushion
- Anti-tip bars