



 Breezy® *Elegance*

Starting at \$499



Elegance

SILVER


- Ultra lightweight, brilliant silver frame
- Accommodates foot propulsion
- Easily removable, lightweight, and low maintenance spoke wheels
- Secure transfers
- Comfortable, contoured cushion
- Improved back support



Elegance GOLD

Includes Silver package options PLUS:

- Increased control and maneuverability
- Increased comfort for the user
- Increased comfort for the caregiver
- Reduced shoulder tension



Height Adjustable
Push Handles with
Locking Brake Lever



Height Adjustable
Armrest



Height-adjustable
push handles

Height-adjustable,
flip-back armrest

Contoured plush
cushion

Integrated drum
braking system

Elegance PLATINUM

Includes Gold package options

PLUS:

- Ultimate comfort
- Ultimate back support



7-Strap Tension
Adjustable Back



7-strap tension
adjustable backrest

Contoured ULTRA
plush cushion

Accessorize for even more convenience

Anti-tips

Inserted into the rear frame of the chair, anti-tips will give you confidence and extra stability when climbing steep hills or curbs.

Elevating Legrests

For clients requiring the elevation of their legs for various conditions.

Step Tube

This easily attachable insert can be "stepped" on by a caregiver to help boost the chair and user up over a curb.

Extra Cushion

For times when your cushion becomes soiled or damaged, having an extra cushion prevents any delays in adventure!

Crutch Holder

An attachment that allows crutches to be carried on the back of the Elegance.

Get your loved one an Elegance

The Breezy Elegance is stocked and available for pick up at a Premium Retailer near you.

Visit the Elegance website for additional information or to find a Premium Retailer:



Specifications	
Weight	32 lbs. (16 kgs)
Frame Material	Aluminum
Frame Color	Brilliant Silver
Weight Capacity	275 lbs. (125 kgs)
Seat Widths	16", 18" (41 cm, 46 cm)
Seat-to-Floor Height	17" (43 cm) expandable to 20" (51 cm) <i>consult your trained Breezy dealer</i>
Legrest & Footplates	Angle adjustable footplate
Rear Wheels	24" x 1" (61 cm x 2.5 cm) solid type
Seat Upholstery	Standard padded nylon