

# Escape!



## Freedom to go anywhere...

**Escape** is available in 3 seat heights to fit almost any user and has a compact folding system that is easy to fold with a simple lift of the release handle. This design makes it easier to lift into the trunk of a car but, as an added benefit, it also stands when folded for easy storage when not in use.

**Escape** features a large padded seat for added comfort, pressure reducing anatomical handgrips to maintain your hand strength and the easiest to apply hand brakes available on any walker.

In addition **Escape** offers a depth adjustable padded backrest, spring loaded levers for easy handle height adjustment and innovative Gentle Touch braking system with no cables to catch on objects or get in your way. **Escape** turns a world of obstacles into a future of possibilities.

## Escape Rollator Features

Escape!



Available Colours:



- Light weight aluminum construction
- Large padded seat for added comfort
- Adjustable padded backrest
- Folding shopping bag
- Easy to apply, cable-free brakes
- Easy to navigate with large 8" wheels



Easy-to-adjust  
Handle Height



Anatomical  
Handles



Curb Climber



Spring-assist  
Folding

Folds with a single lift  
of the seat handle

**Escape** stands  
independently  
when folded

## Optional Accessories

*Sold separately*



Cane Holder *(cane not included)*



Cup Holder



Slow Down Brake



Plastic Tray

# Specifications

## Model Numbers

Charcoal	500-10191	500-10211	500-10241
Red	500-10195	500-10215	500-10245

Description	Super Low	Low	Standard
Seat Height	19" 48cm	21" 54cm	24" 61cm
Handle height range	28 - 36" 71 - 91.5cm	29.5 - 37.5" 75 - 95.3cm	32 - 40" 81.3 - 101.6cm
Distance between handles	18" 45.7cm	18" 45.7cm	18" 45.7cm
Distance between wheels	19" 48.25cm	19" 48.25cm	19" 48.25cm
Seat size	18.5 x 9.5" 47 x 24cm	18.5 x 9.5" 47 x 24cm	18.5 x 9.5" 47 x 24cm
Overall width	25" 63.5cm	25" 63.5cm	25" 63.5cm
Overall length	25.6" 65cm	25.6" 65cm	25.6" 65cm
Depth folded	11" 28cm	11" 28cm	11" 28cm
Wheel diameter	8" 200mm	8" 200mm	8" 200mm
Weight capacity of rollator	250lbs 113Kgs	250lbs 113Kgs	250lbs 113Kgs
Weight capacity of shopping bag	10lbs 5Kgs	10lbs 5Kgs	10lbs 5Kgs
Product weight	15.4lbs 7Kgs	15.6lbs 7.1Kgs	15.8lbs 7.2Kgs

Escape™

Escape™

Limited lifetime warranty.  
Please visit [triumphmobility.com](http://triumphmobility.com)  
for additional information.

Authorized Dealer



Triumph Mobility, Inc.

151 Bentley Ave. Unit 4  
Nepean, Ontario, K2E 6T7

Tel: 1-855-546-0711  
Fax: 1-866-225-3059  
e-mail: [info@triumphmobility.com](mailto:info@triumphmobility.com)

[www.triumphmobility.com](http://www.triumphmobility.com)

R-031015\_Escape-Brochure | Revision # 4 | Printed in Canada

Escape™



THE NEXT STEP  
IN MOBILITY

where will it take you?

