


## **SPECIFICATIONS**

- Frame Material: 1 1/4" Schedule 40 Medical Grade PVC
- Weight Capacity: 300lbs [136kg]
- Tires & Wheels: (4ea) 19" Dia (49cm) gray non-corrosive Wheel-eez(tm)
- Seat: 22" wide cushioned (approx. 3"-4" thick) (blush mesh)
- Umbrella: 5' diameter (color may vary)
- Footrest: Fixed
- Dual action handle parking brake
- Rear storage bag included (blue mesh)
- Easily breaks down into smaller pieces for transport



PRODUCT:

**F-6001**