

The Aftermarket Group Introduces



FIT GRIPS™

- Ergonomic, Affordable Handrim Grips
- More Control & Comfort
- Latex Free, Dishwasher Safe, Made in USA
- One Size Fits All



FIT GRIPS PRO

High Endurance for Active Use

Power • Efficiency • Control



- High endurance rubber grips designed for active indoor/outdoor use
- Installs in seconds and fits any 24"-26" push rim
- Latex Free
- Great for use on ultralight wheelchairs
- Price includes a set of 2 Grips in Black, Blue, or Red

FIT Grips PRO Part #:

BLACK: FTGRPPRO-BLK

BLUE: FTGRPPRO-BLUE

RED: FTGRPPRO-RED

FIT GRIPS COMFORT

Softer for Extra Comfort in the Home

Comfort • Ergonomic • Ease



- Extra Soft rubber grips designed for comfort in the home
- Installs in seconds and fits any 24"-26" push rim
- Latex Free
- Great for use with standard wheelchairs
- Price includes a set of 2 Grips in Black, Blue, or Red

FIT Grips COMFORT Part #:

BLACK: FTGRPC-BLK

BLUE: FTGRPC-BLUE

RED: FTGRPC-RED

Latex Free



Dishwasher Safe



Made In USA



Call Today to order your FIT GRIPS!

39400 Taylor Parkway • North Ridgeville, OH 44039 • Phone: 1-888-TAG-8200 • Fax: 1-800-468-0033 • www.AftermarketGroup.com

Fit Grips is a registered trademark of Intelliwheels.

FIT GRIPS F.A.Q

What size Fit Grips do I need to fit on my wheelchair?

Fit Grips are one size fits all. Any Fit Grip will fit on any 24", 25", or 26" wheelchair wheel.

What is the difference between Fit Grips Comfort and Fit Grips Pro?


Fit Grips Pro are a great solution for active users looking for durable, all-weather grip to improve power performance in both indoor and outdoor environments. Fit Grips Comfort provide additional control, comfort and improved ergonomic grip for wheelchair users pushing mostly in the home.

How long do Fit Grips last?


The Fit Grips are made of high quality specialized rubber that will endure at least 8-12 months of daily use. Unlike rubberized or coated handrims, should the grips wear over time they can be easily and affordably replaced.

How do you install Fit Grips?

1. With the teeth facing in, line up the Fit Grips over one section of the hand rim.



2. Pull the Fit Grips over the rest of the rim and align as preferred.



Are the Fit Grips Latex Free?

Yes. Fit Grips are made of soft, dishwasher safe, flexible rubber. The soft material has enough stretch for easy installation yet enough grip to stay put once it is installed.

How do Fit Grips work in the cold?

The Fit Grips maintain grip in the cold just as they do at room temperature. The Fit Grips feel much warmer than a typical metal handrim in cold weather, giving you more comfort and better control with or without gloves.

How do Fit Grips work in the heat?

The Fit Grips will retain their shape even in hot weather. However, we advise that Fit Grips should not be left in extreme heat over 150 degrees, or a hot car in the middle of the summer.

Can you take them on and off?

With the patent pending design, Fit Grips can be taken on and off your handrims as many times as you desire. They are simple to install and do not require any glue or tape, and they don't leave any type of residue on your handrims.

Where can I buy Fit Grips?

Fit Grips are available for purchase at The Aftermarket Group. Contact your sales representative or our friendly customer service reps to place your order.



Call Today to order your FIT GRIPS!

39400 Taylor Parkway • North Ridgeville, OH 44039 • Phone: 1-888-TAG-8200 • Fax: 1-800-468-0033 • www.AftermarketGroup.com

**Fit Grips is a registered trademark of Intelliwheels.*