

# FLEXION FOOTBOX™ - SINGLE



**HCPCS CODING:** K0108  
WHEELCHAIR ACCESSORY

## PRODUCT FEATURES

- + Available in four sizes: 5"W, 6"W, 7"W, and 8"W.
- + Easy angle adjustment from 90° to 160° with a single screw.
- + Lightweight, durable aluminum shell.
- + After drilling mounting holes in footbox, bolt on to chair for permanent placement.
- + All surfaces and edges are foam padded for added comfort.
- + Comes standard with Comfort-Tek cover.



### STANDARD COVER



**COMFORT-TEK FABRIC:** the fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multi-dimensional stretch and helps alleviate pressure on bony prominences.

### ANGLE ADJUSTMENT

Easy angle adjustment from 90° to 160° with a single screw.



### FLEXION FOOTBOX™ SINGLE

PART #	WIDTH X HEIGHT X DEPTH
<b>FBS-5</b>	5.25"W x 11"H x 11"D
<b>FBS-6</b>	6.25"W x 11"H x 11"D
<b>FBS-7</b>	7.25"W x 11"H x 11"D
<b>FBS-8</b>	8.25"W x 11"H x 11"D

NOTES

# FLEXION FOOTBOX™ - DOUBLE

**HCPCS CODING:** K0108  
WHEELCHAIR ACCESSORY



## PRODUCT FEATURES

- + Available in four sizes: 14"W, 16"W, 18"W, and 20" chairs.
- + Easy angle adjustment from 90° to 160° with a single screw.
- + Lightweight, durable aluminum shell.
- + After drilling mounting holes in footbox, bolt on to chair for permanent placement.
- + All surfaces and edges are foam padded for added comfort.
- + Comes standard with Comfort-Tek cover.



### FLEXION FOOTBOX™ DOUBLE

PART #	CHAIR SIZE	WIDTH X HEIGHT X DEPTH
<b>FBD-14</b>	14"	12"W x 11"H x 11"D
<b>FBD-16</b>	16"	14"W x 11"H x 11"D
<b>FBD-18</b>	18"	16"W x 11"H x 11"D
<b>FBD-20</b>	20"	18"W x 11"H x 11"D

### STANDARD COVER



**COMFORT-TEK FABRIC:** the fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multi-dimensional stretch and helps alleviate pressure on bony prominences.

### ANGLE ADJUSTMENT

Easy angle adjustment from 90° to 160° with a single screw.



NOTES