

# GO-GO<sup>®</sup> ES

## Outstanding Value. Exceptional Durability.


The new Go-Go<sup>®</sup> ES Scooter from Pride<sup>®</sup> delivers high-performance at an affordable price. The lightweight Go-Go ES features easy disassembly and Pride's exclusive black, non-scuffing tires. Designed for exceptional durability, the Go-Go ES provides reliability and outstanding value.



### Easy disassembly.



**Pride<sup>®</sup>**  
Mobility Products Corp.  
— LIVE YOUR BEST<sup>®</sup> —



## Specifications

Model numbers	S83 (3-Wheel)
Max. weight capacity	275 lbs.
Maximum speed <sup>1</sup>	Up to 4 mph
Ground clearance <sup>2</sup>	2.5" at mid-frame
Turning radius <sup>2</sup>	32.75"
Overall length <sup>2</sup>	38"
Overall width <sup>2</sup>	19.5"
<b>Tires</b>	
Front	8" solid
Rear	8" solid
Anti-tips	Rear
Standard seating	Foldable, compact, vinyl-covered molded plastic
Drivetrain	24-volt DC motor, sealed mini-transaxle
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	Off-board, 1.5 amp
Per-charge range <sup>1,5</sup>	Up to 6 miles
Battery requirements <sup>3</sup>	(2) 12AH, 12V deep-cycle
Battery pack weight <sup>4</sup>	16 lbs.
Weight w/o battery pack <sup>2</sup>	72.5 lbs.
Heaviest piece <sup>2</sup>	27.5 lbs. (front section)
Warranty	5-year limited on frame; 2-year limited on electronics; 1-year limited on motors

Standard colors



Red


## Features

- Pride's exclusive black, non-scutting tires
- Front frame-mounted seat post offers maximum stability
- Easily disassembles into 6 super lightweight pieces for convenient transport and storage
- Modular design for easy serviceability
- Convenient off-board dual voltage charger can charge battery pack on-board or off-board


## Pride® Lifts & Ramps

These Pride® Lifts are compatible with the Go-Go® ES.

- Backpacker® Plus
- Backpacker® MV
- Power Tote Mini
- Outlander



Backpacker AVP



SilverBoom 250

### Preliminary specifications; subject to change

<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of + or - 10%.

<sup>2</sup>Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup>AGM or gel cell type required.

<sup>4</sup>Battery weight may vary based on manufacturer.

<sup>5</sup>This range was obtained with the 4-wheel model.

*The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.*