

Silver Ion
Antimicrobial CoverNEW
*****Multiple Therapeutic
Benefits in 1 System1P
ACK

Harmony True Low Air Loss Tri-Therapy Mattress Replacement System

Item #	Product Description	Product Dimensions	UOM
14200	Harmony True Low Air Loss Tri-Therapy Mattress Replacement System	80" (W) x 36" (L) x 8" (D)	1/ea
<ul style="list-style-type: none"> Feature filled True Low Air Loss System for microclimate and patient comfort. Easy-to-use, intuitive 1300 LPM Blower System with Pulsation and Alternating Pressure for the prevention and treatment of all stages of pressure ulcers. Comfort setting control based on patient's weight for optimal comfort and therapy. Seat Inflation feature for added support in the Fowler position. Fully inflates in under two minutes with Autofirm feature. Digital micro-processing continuously monitors pressure for optimal therapeutic pressure redistribution. Visual and audio low pressure, power failure and system failure alarms. 2" foam base included for power outage protection. Silver Stretch Cover: permanent silver-ion kills 99.5% of germs; fluid resistant, vapor permeable, soft and supple. Cable Management: 5 side mattress loops keep cables off the floor. Five minute panel lockout feature protects from unintentional setting changes. Handy, durable carrying bag included with each system. 			

specifications

Airflow	1300 LPM
Cycle Time	5, 10, 15, 20 Minutes
Type	Blower
Alarms	Visual/Audio
Mattress Cover Material	Silver Ion 4 way stretch
Weight Capacity	400 lbs
Warranty	1 Year Limited