

INSIDE LIFT SYSTEMS



EASY AS 1, 2, 3!

Harmar's Inside Lifts come standard with many useful features and provide robust heavy-duty lifting capacities for today's heavier scooters and power chairs. Harmar offers a solution for virtually every chair and vehicle on the market. The Universal design of the AL400 Series and the AL800 allows installation into either side of most vehicles.

STEP 1.

Attach the lifting strap to the power-chair or scooter.



STEP 2.

Raise and swing into the vehicle.



STEP 3.

Close the hatch and go!



* Typical AL400 Series application.

Harmar Mobility produces lift and ramp products that enhance people's mobility, independence and quality of life. Products include Inside, Outside and Hybrid Auto Lifts, along with Single-fold and Multi-fold Portable Ramps.

HYBRID LIFT SYSTEMS

Hybrid Lifts provide the ultimate in convenience and simplicity of use by combining the best of both inside and outside lifts into unique solutions.



OUTSIDE LIFT SYSTEMS

Harmar is the leader and benchmark all others look to for ease of use, quality and performance. Models capable of fitting Class I, II or III hitches.



MICRO LIFT SYSTEMS

Designed to combine quality with budget, these user-friendly models feature manual, powered or combination capabilities.



ACCESS RAMP SYSTEMS

Designed for ease of use at home or while traveling. These commercial quality ramps are built to access vehicles, landings and other unkind elevation changes.



Harmar Mobility
2075 47th St., Sarasota, FL 34234
Phone: 800.833.0478
Toll-Free Fax: 866.234.5680

www.harmar.com

INSIDE LIFT SYSTEMS

MODELS AL410, AL420 & AL800



INTRODUCING...

America's Leader in Lifts & Ramps

Harmar Mobility designs, engineers and produces American made lift and ramp products that enhance people's mobility, independence and quality of life. Our products include a complete line of Auto Lifts and Access Ramps.

Harmar Mobility products are available exclusively through a nationwide network of authorized dealers. Customers receive care from the nation's top mobility dealers. Harmar's dealers will assist in determining the best product to fit your needs and budget. It is imperative to have an authorized Harmar dealer demonstrate, install and service your Harmar lift.

Harmar's lifts and ramps are designed and built with the "Harmar difference" - they are simpler, stronger, lighter and easier to use.

Chad C. Williams
President, Harmar Mobility



Model AL420 Shown



www.harmar.com

AL410 Inside Lift Specifications:

Lifting Capacity:	400 lbs.
Manual Rotation:	Yes
Adjustable Boom Height:	30 - 40"
Adjustable Boom Length:	23.5 - 33.5"
Adjustable Boom Pitch:	0 - 30 deg.
Adjustable Rotation Limits:	Yes
Manual Crank Backup:	Yes

AL420 Inside Lift Specifications:

Lifting Capacity:	400 lbs.
Powered Rotation:	Yes
Adjustable Boom Height:	30 - 40"
Adjustable Boom Length:	23.5 - 33.5"
Adjustable Boom Pitch:	0 - 30 deg.
Adjustable Rotation Limits:	Yes
Manual Crank Backup:	Yes

AL800 Inside Lift Specifications:

Lifting Capacity:	350 lbs.
Powered Rotation:	Yes
Adjusts for Cab Size:	Yes
Adjusts for Bed Length:	Yes
Adjustable Rotation Limits:	Yes
Manual Crank Backup:	Yes

Eight years of customer input and design innovation have combined to create Harmar's AL400 Series and AL800 Inside Lifts. With an appearance that integrates with today's vehicle interiors, Harmar's latest Inside Lifts include Best-in-Class features and lifting capacities to offer high-quality solutions for virtually every chair and vehicle. No other products on the market have the level of safety, quality and adjustability found with Harmar Inside Lifts.



Model
AL420
FULLY POWERED INSIDE LIFT

- Fully-powered rotation for effortless loading
- Only 30" of headroom required
- Extremely adjustable base and boom
- Mounts on either side of vehicle
- Weather resistant
- Remote hand control
- Removable hook with swivel

Model
AL410
UNIVERSAL INSIDE LIFT

- Manual rotation with adjustable limits and lock
- Only 30" of headroom required
- Extremely adjustable base and boom
- Mounts on either side of vehicle
- Weather resistant
- Remote hand control
- Removable hook with swivel



Model
AL800
SIDE-LOADING INSIDE LIFT

- Brings chair to driver or passenger door
- Fully-powered for effortless loading
- Extremely adjustable base and boom
- Weather resistant
- Only 2 motors to control with remote hand control
- Removable hook with swivel



A totally new and innovative design from Harmar that delivers a wheelchair to the side door of today's most popular vehicle, the pickup truck. The AL800 provides ample clearance and reach with one less motor than the competition, making operation much simpler. With 350 lbs. of capacity, the AL800 will load virtually any mobility device available with ease into the bed.