


Invacare® Personal Mobility Vehicles



Yes, you can.™

Invacare three- and four-wheel Personal Mobility Vehicles (PMVs) offer the best in aesthetics, reliability, performance and features. Now, the line has progressed even further with innovative features like the ergonomically-enhanced delta tiller, with its sleek and stylish good looks and optimum functionality. These PMVs also feature an expanded line of fully-functional accessories. Best of all, every PMV offers versatility, from compact indoor maneuverability to rugged outdoor performance and power.



Easy-to-use controls feature a speed control, LED battery gauge, and a canopy over the gauge to improve readability.



The beautifully-contoured rear shroud gives the Zoom-3 an elegant look.



Invacare® Zoom™-3

With its compact design and stylish good looks, the new Invacare® Zoom™-3 offers a 29" turning radius for maximum maneuverability and performance indoors. The fold-down tiller and removable seat make transporting simple. Personal items fit conveniently into the front basket.

FEATURES

- 29" turning radius
- Three-wheel frame design
- Candy Red color (as shown)
- Simple fold-down tiller design
- Quiet, powerful transaxle
- 38.5" overall length, 22" overall width
- Contoured seat with cloth upholstery
- Foam-filled tires
- Approximately 10-15 mile battery range
- 0-3.5 mph speed range
- 250 lb. weight capacity
- Remote battery charger
- Front basket



Components disassemblable for easy transport

Invacare® Lynx™ SX-3

Known for its superior looks, the Invacare Lynx SX-3 is a small-sized three-wheel PMV designed to go anywhere. It allows for maximum maneuverability indoors, and is great for outdoor use, too.

The Lynx SX-3 features beautifully contoured shrouds available in two color options. The tiller adjusts easily and securely to a comfortable position, and the urethane armrests flip back for easy entry or exit.

FEATURES

- Three-wheel frame design
- Choice of two colors (Candy Red or Deep Blue)
- Quiet, powerful transaxle
- 43.5" overall length, 23" overall width
- Contoured seat with swivel and cloth upholstery
- Foam-filled tires
- Angle- and width-adjustable arms
- Approximately 10-15 mile battery range
- 0-4.5 mph speed range
- 250 lb. weight capacity
- On-board battery charger
- Front and floor baskets



The 3-amp charger is mounted on-board to ensure that it is always available to recharge the batteries.



Easy-to-use controls feature a speed control, LED battery gauge, and a canopy over the gauge to improve readability.

The Lynx SX-3 is not intended for street or highway use.



Tiller includes a removable key, speed control dial, an LED battery gauge and a battery gauge canopy to improve readability.



The beautifully-contoured rear shroud gives the LX-3 and LX-3Plus an elegant look and features two visible reflectors. (LX-3Plus shown)



The LX-3 model comes standard with a contoured seat with cloth upholstery and a traditional tiller.



Invacare® Lynx™ LX-3 and LX-3Plus


The Invacare Lynx LX-3 and LX-3Plus are mid-sized three-wheel PMVs that offer superior looks, performance, and functionality. They are versatile enough, both indoors and outdoors, to suit a wide range of owners.

The LX-3 features a full tiller shroud that adds to the overall good looks of the PMV, as well as a contoured seat with cloth upholstery. The Lynx LX-3Plus features the ergonomically-enhanced delta tiller, as well as a non-sliding deluxe vinyl seat. Front and floor baskets on both models allow room to carry personal items. A headlight is also standard on both models.

FEATURES

- Three-wheel frame design
- Choice of three colors (Candy Red, Deep Blue or Magic Black)
- Quiet, powerful transaxle
- 46" overall length, 24.5" overall width
- LX-3 standard seat is a contoured seat with cloth upholstery, swivel, and width- and angle-adjustable arms
- LX-3Plus standard seat is a deluxe no slide, vinyl seat with swivel; and height-, width-, angle- and depth-adjustable arms
- Foam-filled tires
- Approximately 18-23 mile battery range
- 0-5 mph speed range
- 300 lb. weight capacity
- On-board battery charger

The Lynx LX-3 and LX-3Plus are not intended for street or highway use.



Ergonomically enhanced delta tiller for optimum function and style. Features include a removable key, LED battery gauge and speed control.*



The beautifully contoured rear shroud gives the LX-4 an elegant look and features two visible reflectors.



The tiller adjustment is both easy-to-use and secure; five positions for maximum comfort.

Invacare® Panther™ LX-4

Perfect for a wide range of drivers, the Invacare Panther LX-4 is a mid- to large-sized four-wheel PMV that looks as good as it performs.

The Panther LX-4 features four 3" x 10" tires for outdoor use, and a deluxe padded vinyl seat to provide comfort and adjustability. The ergonomically-enhanced delta tiller, with its sleek and stylish good looks and optimum functionality, is standard on this model. A headlight and front bumper are also standard.

FEATURES

- Four-wheel frame design
- Choice of three colors (Candy Red, Deep Blue or Magic Black)
- Quiet, powerful transaxle
- 48" overall length, 24.5" overall width
- Deluxe vinyl seat with swivel and slide
- Height-, width-, depth-, and angle-adjustable arms
- Foam-filled tires
- Approximately 18-23 mile battery range
- 0-5 mph speed range
- 300 lb. weight capacity
- Remote battery charger
- Headlight
- Front basket

*Standard on the Lynx LX-3^{AS} and Panther LX-4

The Panther LX-4 is not intended for street or highway use.



Ergonomically enhanced wood grain delta tiller for optimum function and style. Features include a removable key, LED battery gauge and speed control.*



The beautifully-contoured rear shroud gives the Panther MX-4 an elegant look and features turn signals and taillights.



Designed with four-wheel suspension to ensure a comfortable ride.

*Standard on the Panther MX-4



Invacare® Panther™ MX-4

The Invacare Panther MX-4 is a large-sized, four-wheel PMV designed for outdoor use and maximum performance. Full suspension allows the Panther MX-4 to be rated for 400 lb. weight capacity. The beautiful and comfortable ergonomically-enhanced delta tiller comes standard on the MX-4.

The Panther MX-4 has four 4" x 12" tires. A Hi/Lo speed switch allows a range of 0-7 mph. Front and floor baskets allow the user to take personal items along for the ride.

FEATURES

- Four-wheel frame with four-wheel suspension
- Champagne color (as shown)
- Quiet, powerful transaxle
- 56" overall length, 26" overall width
- Extra wide (22") deluxe vinyl seat with swivel and slide
- Height-, width-, angle- and depth-adjustable arms
- Foam-filled tires
- Approximately 20-25 mile battery range
- Hi/Lo speed switch
- 0-7 mph speed range
- 400 lb. weight capacity
- Full lighting package
- On-board battery charger
- Front and floor baskets

The Panther MX-4 is not intended for street or highway use.


Seats and Accessories



Contour Seat with 360° swivel comes standard with cloth upholstery on the Lynx SX-3 and LX-3 models only. The seat features flip-back width- and angle-adjustable armrests. A vinyl cover is available by ordering part #1088758.



Deluxe Vinyl Seat with 360° swivel features flip-back height-, width-, depth-, and angle-adjustable arms. Available on the Lynx LX-3^{Plus} and Panther LX-4 models only. The LX-4 Deluxe Seat also includes a slide feature.



Extra Wide (22") Deluxe Vinyl Seat includes swivel and slide features and flip-back height-, width-, depth-, and angle-adjustable arms. Available on the Panther MX-4 model only.



Walker Holder helps user safely and conveniently transport walker.




Crutch and Cane Holder keeps user's walking aids close by.




Oxygen Holder carries a "D" size cylinder.



Rear Basket provides additional room for user's personal items.



Extended Throttle Levers provide added length for easier access to the throttle control. Rear View Mirror is available on all models.



Safety Flag is available on all models.

Seat Positioning Strap is available on all models. *(not shown)*

Available Colors



Candy Red
Zoom-3, SX-3,
LX-3, LX-3^{Plus},
LX-4



Deep Blue
SX-3, LX-3,
LX-3^{Plus}, LX-4



Magic Black
LX-3^{Plus}, LX-4,



Champagne
MX-4



Specifications

	Zoom-3	Lynx SX-3	Lynx LX-3	Lynx LX3Plus	Panther LX-4	Panther MX-4
Style	3 wheel	3 wheel	3 wheel	3 wheel	4 wheel	4 wheel w/supension
Candy Red	Candy Red	Candy Red Deep Blue	Candy Red Deep Blue	Candy Red Deep Blue Magic Black	Deep Blue Magic Black	Champagne
Dimensions						
Overall Width	22"	23"	24.5"	24.5"	24.5"	26"
Overall Length*	38.5"	43.5"	46"	46"	48"	56"
Seat Style	contour, cloth	contour, cloth	contour, cloth	deluxe vinyl, no slide	deluxe vinyl, slide	deluxe, vinyl slide
Seat Depth	14.5"	16.5"	16.5"	18"	18"	18"
Seat Width	16.5"	16.5"	16.5"	18"	18"	22"
Seat Height (adjusts in 1" increments)	13.5"-15.5"	14.5"-16.5"	14.5"-17.5"	15"-18"	16.5"-19.5"	17"-19"
Arm Width	23"	22"-25"	22"-25"	22"-25"	22"-25"	24"-27"
Arm Height	7"	7"	7"	7" (adjustable)	7" (adjustable)	6" (adjustable)
Front Tire Size	2.5" x 7"	2.5" x 7"	2.5" x 9"	2.5" x 9"	3" x 10"	4" x 12"
Rear Tire Size	3" x 8"	3" x 8"	3" x 10"	3" x 10"	3" x 10"	4" x 12"
Weight (of heaviest components)						
Base Section	56 lb.	39 lb.	47 lb.	52 lb.	53 lb.	84 lb.
Rear Section	n/a.	35 lb.	38 lb.	38 lb.	38 lb.	58 lb.
Seat and Arms	20 lb.	22 lb.	22 lb.	31 lb.	35 lb.	43 lb.
Battery (each)	10 lb.	14 lb.	23 lb.	23 lb.	23 lb.	31 lb.
Total Weight	96 lb.	126 lb.	156 lb.	170 lb.	175 lb.	250 lb.
Capabilities						
Speed Range	0-3.5 mph	0-4.5 mph	0-5 mph	0-5 mph	0-5 mph	0-7 mph
Maximum Turning Radius	29"	33"	40"	40"	58"	70"
Weight Limit	250 lb.	250 lb.	300 lb.	300 lb.	300 lb.	400 lb.
Incline Capability	8°	8°	8°	8°	8°	8°
Charger	remote	on-board	on-board	on-board	remote	on-board
Wood Grain Decals	no	yes	yes	yes	yes	yes
Battery Size	12 ah	17 ah**	31 ah**	31 ah**	31 ah**	38 ah**
Battery Range***	10-15 miles	10-15 miles	18-23 miles	18-23 miles	18-23 miles	20-25 miles
Warranty						
Frame	3 years	3 years	3 years	3 years	3 years	3 years
Electronics	18 months	18 months	18 months	18 months	18 months	18 months
Transaxle	18 months	18 months	18 months	18 months	18 months	18 months

Standard Features and Accessories

Baskets						
Floor	n/a	standard	standard	standard	n/a	standard
Front	standard	standard	standard	standard	standard	standard
Delta Tiller	n/a	n/a	n/a	standard	standard	standard
Headlight	n/a	n/a	standard	standard	standard	standard
Taillight	n/a	n/a	n/a	n/a	n/a	standard
Contour Seat/Cloth	standard	standard	standard	n/a	n/a	n/a
Deluxe Seat/Vinyl	n/a	n/a	n/a	standard no slide	standard slide	n/a
Extra Wide Deluxe Seat (22")	n/a	n/a	n/a	n/a	n/a	standard slide
Walker Holder	optional	optional	optional	optional	optional	optional
Crutch & Cane Holder	optional	optional	optional	optional	optional	optional
Oxygen Holder	optional	optional	optional	optional	optional	optional
Rear Basket	optional	optional	optional	optional	optional	optional
Extended Throttle Levers	optional	optional	optional	optional	optional	optional
Seat Positioning Strap	optional	optional	optional	optional	optional	optional
Safety Flag	optional	optional	optional	optional	optional	optional
Rear View Mirror	optional	optional	optional	optional	optional	optional

* Measurement includes anti-tippers for Zoom-3, SX-3, LX-3, LX-3Plus and Panther LX-4.

** Batteries not included with the PMV.

*** Range varies depending on terrain, battery condition, and operator weight.

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Invacare, Lynx, Panther,
Zoom and Yes, you can.
are trademarks of
Invacare Corporation.
Specifications are subject
to change without
notification.

© 2001
Invacare Corporation

Form No. 00-255
Rev. 10/01