



Transport Chairs

Yes, you can:





Yes, you can.®

The Invacare® Lightweight Aluminum Transport Wheelchairs are designed for safety and comfortable patient transport. The lightweight frames weigh only 19 lb. and incorporate a fold-down back, making these chairs easy to transport and store.



Features

- Black nylon upholstery is flame retardant and bacteria resistant
- Full-length arms are padded for comfort
- Carry straps on seat allow for ease of folding and lifting
- Footrest painted to match frame finish
- Fold down back makes transport and storage easy
- Weighs 19 lb. (w/o riggings)



LTTB19HBFR, Lightweight Aluminum Transport chair in Blue.



Lightweight Aluminum Transport chair folded is compact for ease of transport and storage.

ALR19HBFR, Lightweight Aluminum Transport chair in Red.

Specifications

Model Numbers	ALR19HBFR ALB19HBFR	LTR19FR LTTB19FR
Chair Dimensions		
Seat width	19"	19"
Seat depth	18"	15.5"
Back height	17.5"	17.5"
Seat-to-floor height	19.25"	19.25"
Overall Dimensions		
Overall width	23.5"	23.5"
Overall depth w/o riggings	24.25"	24.25"
Overall height	37.5"	36.5"
Weight w/o riggings	19 lb.	19 lb.
Weight limitation	250 lb.	250 lb.

Note: All specifications are approximate

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Invacare, the Blue Medallion Design and Yes, you can., are registered trademarks of Invacare Corporation. Specifications are subject to change without notification.

© 2005
Invacare Corporation

Form No. 04-022
Rev. 10/05



"Invacare says 'Yes, you can', and I couldn't have said it better myself."