



Jazzy[®]

600

The Jazzy[®] 600 provides excellent performance and outstanding stability. Its compact design and standard Active-Trac[®] ATX Suspension allow it to excel indoors and out, making it a superb, all-around power chair choice.



Pride[®]
Mobility Products Corp.
LIVE YOUR BEST[®]

SPECIFICATIONS

Weight capacity	300 lbs.
Maximum speed	Up to 5.3 mph
Ground clearance	2.8"
Turning radius	22.5"
Overall length	35" without front riggings
Overall width	23"
Drive wheels	14"
Front casters	5" solid OMNI-casters
Rear casters	6" solid OMNI-casters
Suspension	Active-Trac® ATX Suspension
Drivetrain	Two-motor, Mid-Wheel 6®
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	60 amp, PG VR2
Specialty controls	No
Battery charger	5 amp, off-board
Per-charge range	Up to 15.9 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	NF-22
Battery weight	38 lbs. each
Base weight	127 lbs.
Standard seat weight	37 lbs. (medium-back)
Max. seat dimensions	22" x 20"
Warranty	K0822/K0823 Models: Lifetime limited warranty on frame; 13-month limited warranty on electronics; 13-month limited warranty on drive motors K0848/K0849/K0856/K0857 Models: Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-month limited warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

MODELS

Jazzy 600 2S-SS	(Code: K0822)
Jazzy 600 2S-C	(Code: K0823)
Jazzy 600 3S-SS	(Code: K0848)
Jazzy 600 3S-C	(Code: K0849)
Jazzy 600 3SP-SS	(Code: K0856)
Jazzy 600 3SP-C	(Code: K0857)

FEATURES

- OMNI-Casters (nylon, spherical-shaped casters) on front and rear to prevent wheel hang-ups
- Side-mounted, easily accessible freewheel levers
- Active-Trac® ATX Suspension for enhanced performance over more varied terrain
- Easy rear access to batteries

OPTIONS

- Solid seat pan
- Elevating leg rests
- Cup holder
- Oxygen tank holder
- Swing-away joystick
- Cane/crutch holder
- Swing-away footrests
- Walker holder
- TRU-Balance® Power Tilt (Group 3 only)
- Synergy® Seat

RECOMMENDED LIFTS

Pride Lifts & Ramps

- Backpacker® AVP
- Backpacker® MV
- Backpacker® Plus
- Outlander
- PrideBoom 250




COLOR CHOICES



Candy
Apple Red



Viper
Blue




Onyx
Black



Silver



Emerald
Green



Sapphire
Blue

(US) 800-800-8586
(Canada) 888-570-1113
www.pridemobility.com

PRIDE®
PERFORMANCE MOBILITY
— LIVE YOUR BEST® —