

# Jazzy<sup>®</sup>

## SPORT 2

The Jazzy<sup>®</sup> Sport 2 delivers a reliable blend of power, performance and style. In-line, front-wheel drive technology gives the Jazzy Sport 2 stability and maneuverability for solid performance indoors and out.

**Jazzy Sport<sup>®</sup> 2**  
shown with Jazzy<sup>®</sup> Armor<sup>®</sup> (standard)



**Pride<sup>®</sup>**  
Mobility Products Corp.

## Specifications

<b>Weight capacity</b>	300 lbs.
<b>Maximum speed</b>	Up to 4 mph
<b>Ground clearance</b>	1.5"
<b>Turning radius</b>	24.75"
<b>Overall length</b>	34.75"
<b>Overall width</b>	22.75"
<b>Drive wheels</b>	9"
<b>Front wheels</b>	3" solid casters
<b>Rear wheels</b>	6" solid casters
<b>Drivetrain</b>	Two-motor, in-line, front-wheel drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	40 amp, PG GC 3 controller
<b>Specialty controls</b>	N/A
<b>Battery charger</b>	Off-board, 3 amp
<b>Per-charge range</b>	Up to 12.89 miles
<b>Battery requirements</b>	(2) 12 volt, deep cycle
<b>Battery size</b>	U-1
<b>Battery weight</b>	24.5 lbs. each
<b>Base weight</b>	68.5 lbs.
<b>Standard seat weight</b>	50 lbs. (comfort, high-back)
<b>Max. seat dimensions</b>	20" x 20"
<b>Warranty</b>	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

## Features

- In-line, front-wheel drive technology
- 40 amp, PG GC3 controller
- Easy access to batteries from rear
- Simple main frame design for easy serviceability
- Comfort high-back seat with headrest
- Jazzy® Armor® includes shroud and controller guards to protect from daily wear and tear
- 60" lap belt

## Options

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Saddle bag
- Cane/crutch holder
- Walker holder
- Height adjustable foot platform bracket

## Recommended Lifts

- Backpacker® Plus
- Outlander



## Color Choices



Jazzy Red



Jazzy Blue