

Voyager[®] Portable Easytrack[®] redefining patient handling



Voyager[®] Portable

The Voyager Portable lift and Easytrack[®] systems combine to form a unique portable overhead lifting solution for the home and hospital environment. The Voyager Portable provides the perfect answer for an environment in which a fixed ceiling lifter is not suitable in a homecare, emergency or temporary solution.

- High lift capacity of 440 lbs SWL
- Lightweight design, enables the caregiver to transfer the Voyager Portable lift easily from one room to another
- Soft start and stop ensures the patient is transferred safely and comfortably
- **Dual controls** coupled with a removable handset allow total flexibility
- Simple trolley allows the Voyager Portable to be used with either the Easytrack or the standard Voyager ceiling track
- Unique quick release strap allows instant repositioning of the lift, which prevents the need to lift the unit onto the rail, while saving time and battery power
- Emergency stop/lower and strap twist prevention help ensure safety never becomes an issue
- The removable, chargeable battery allows 24 hour use and ensures the lift is never out of action
- Audible low battery indicator helps to ensure the battery will not run out at a critical moment



Easytrack

The Easytrack system is a unique overhead portable track designed for the Voyager Portable lift. It does not require permanent fixings or structural alterations, and can be assembled in as little as five minutes by just one person.

• Simple Installation: The unique floor-to-ceiling design eliminates the need for specialist fixings, alterations or electrical installations. This allows the system to be used permanently, semi-permanently or temporarily. The Easytrack is the ideal solution for an environment in which a fixed ceiling track cannot be installed.



• Multi-Post Options: The 2, 3 or 4 post options allow the Easytrack to solve patient handling problems in nearly every room layout. Each option is ideal for a different situation: the 2 post for a simple lift in or out of bed, the 3 post for lifting in two directions, and the 4 post for complete room coverage.

- **Portable Kit:** The lightweight aluminum Easytrack comes in portable kit form that can be easily transported along with the Voyager Portable lift. A useful carrying bag is available for both the Easytrack and Voyager Portable.
- Footprint: The small feet of the Easytrack allow use even in the smallest bedroom or bathroom where a conventional floor lifter cannot be maneuvered.
- Easy-Glide Rail: The unique design of the telescopic rail ensures that nearly any patient can be transferred without strain or risk of back injury.
- Rotating Rail Connection: This feature allows the rail or post to be positioned at any angle, no matter how close to a wall or obstruction.

- Extending Rail: Unlike regular gantry tracks, the Easytrack can be extended to suit different locations as and when they are required. The track can be extended when in position to allow adjustment, and then retracted to enable easier storage or transportation.
- **Telescopic Post:** This feature allows the Easytrack to be used with a variety of ceiling heights.
- Safety Systems: A multitude of unique safety systems have been incorporated into the Easytrack to ensure safety. These systems prevent the post losing grip if the floor or ceiling flexes, preventing the post being lowered if the rail is still connected, or the rail disconnecting from the post.

Easytrack Assembly in 4 Easy Steps



Extend posts and rotate to lock in place.



Fix one end of the rail onto the post.



Extend the rail to required length and fix the rail onto the second post.



Fix the end cap to lock the rail in place.

Easytrack Options

There are many variations of the Easytrack system's with 2, 3 and 4 post configurations, thereby providing a solution for almost any setting.

2 Post System:

The standard 2 post Easytrack is suited for transfer over a bed or bath in a peninsular setting:



The 2 post Easytrack



The 2 post Easytrack with wall bracket



The 2 post Easytrack with bath bracket



The 2 post Easytrack without posts

3 Post System:

The 3 post Easytrack[®] is suited for transfer over a bed or bath and down a room. This system negates the need for a traditional curved track:



The 3 post Easytrack

4 Post System:



The 3 post Easytrack with wall bracket



The 3 post Easytrack with bath bracket



The 3 post Easytrack without posts



The 4 post Easytrack operates like a traditional x-y system covering all 4 corners of a room. The 4 post system is not telescopic



An x-y system without posts and with 4 wall brackets.

| Voyager Portable Overhead Lift | | |
|-----------------------------------|--------------------|--|
| Safe Working Load | 440lbs | |
| Soft Start and Stop | \checkmark | |
| Emergency Lower | Manual | |
| Emergency Stop | \checkmark | |
| Low Battery Indicator | Audible | |
| Battery Capacity (load of 75kg) | Up to 30 transfers | |
| Battery Charger | Freestanding | |
| Weight (including battery) | 12lbs | |
| Quick Release Tape Mechanism | \checkmark | |
| Compatible with Full Oxford Sling | | |
| Range | ¥ | |
| Vertical Speed | 4 sec | |

| Easytrack (Portable Track System) | 2 Post System | 3 Post System | 4 Post System |
|-----------------------------------|---------------|---------------|---------------|
| Safe Working Load | 440lbs | 440lbs | 440lbs |
| Ceiling Height Range | 84 - 108" | 84 - 108" | 84 - 108" |
| Rail Length Range | 74 - 124.5" | 74 - 124.5" | 118" |
| Rail Weight (each) | 26lbs | 26lbs | 20lbs |
| Post Weight (each) | 9lbs | 9lbs | 9lbs |
| Total Weight | 44lbs | 79lbs | 114lbs |
| Standard Kit | \checkmark | \checkmark | \checkmark |
| Wall Bracket Option | \checkmark | \checkmark | - |
| Bath Bracket Option | ✓ | \checkmark | ✓ |
| Multi-Wall Bracket Option | ✓ | ✓ | - |

Specifications may vary due to manufacturing tolerances

Warranty

Patient Lifter including Hydraulics: 1 year Easytrack System: 1 year Voyager Portable Battery Unit: 3 months



Post Acute, Acute, HomeCare 5001 Joerns Drive Stevens Point WI 54481 (P) 800.826.0270 (F) 800.457.8827 VA/Government 7027 Hayvenhurst Avenue Van Nuys, CA 91406 (P) 800.966.6662 (F) 800.232.9796 Canadian Office 1000 Clarke Road, Suite 6 London, ON Canada N5V 3A9 (P) 866.546.1151 (F) 519.451.8662