

NEW! LC-525 S E R I E S



The new LC-525 Lift Chair Series features superior comfort, performance and style. A new “zero-gravity” position, back style, infinite-position hand control and super soft fabrics selections make the LC-525 the ultimate in lift chair comfort.

Features and Benefits:

- New! “Zero-gravity” position with active seat box for maximum comfort
- New! Back style for exceptional support
- New! Infinite-position hand control
- Available in Petite Wide, Small, Medium and Large



	LC-525PW	LC-525S	LC-525M	LC-525L
	Infinite-Position “Zero-gravity” Chaise Lounger 5’2” & Below	Infinite-Position “Zero-gravity” Chaise Lounger 5’3” & Below	Infinite-Position “Zero-gravity” Chaise Lounger 5’4” - 6’0”	Infinite-Position “Zero-gravity” Chaise Lounger 5’9” - 6’2”
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	19”	19”	20”	21”
Seat Depth	19”	18”	20”	22”
Seat Width	24”	18”	20”	22”
Top of Back to Seat	25.5”	25.5”	26.5”	27.5”
Back Style	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back
Fully Padded Chaise	Yes	Yes	Yes	Yes
Distance From Wall	17”	17”	17”	17.5”
Overall Chair Width	39”	35”	35”	37”
Chair Weight	118 lbs.	118 lbs.	117 lbs.	120 lbs.
Warranty	2-Year Limited	2-Year Limited	2-Year Limited	2-Year Limited

AVAILABLE FABRICS

NEW CLOUD 9 FABRICS



**BLACK
CHERRY**

PACIFIC


SPRUCE



STONE

WALNUT

LEXIS VINYL



BLACK

**NEW
CHESTNUT**

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

LC525SERIES-SS-04-10-12

1-800-800-8586 (U.S.) • 1-888-570-1113 (Canada)
www.pridemobility.com • www.prideliftchairs.com

Pride
Mobility Products Corp.

PRIDE
PERFORMANCE MOBILITY
— LIVE YOUR BEST! —