

# Mobile Activity Tray

2 year warranty



Adjustable, Maneuverable and Versatile

- Two sizes: Short & Tall
- Height adjustable tray surface
- Angle adjustable tray surface
- Easy to move around a room
- Four tray surface color options



3" (75mm) front swivel casters lock



4" (110mm) fixed rear casters allow for easy access and keep the tray secure





Mobile Activity Tray	
Mobile Activity Tray Short - Lilac	65000102
Mobile Activity Tray Short - Aqua	65000104
Mobile Activity Tray Short - Chocolate	65000105
Mobile Activity Tray Short - Grey	65000109
Mobile Activity Tray Tall - Lilac	65000302
Mobile Activity Tray Tall - Aqua	65000304
Mobile Activity Tray Tall - Chocolate	65000305
Mobile Activity Tray Tall - Grey	65000309

Mobile Activity Tray - Short Specifications	
Lowest Floor-to-tray surface	16" H
Highest Floor-to-tray surface	21" H
Weight	24lbs.
Tray Dimensions	28" W x 20" D

Mobile Activity Tray - Tall Specifications	
Lowest Floor-to-tray surface	20" H
Highest Floor-to-tray surface	28" H
Weight	29lbs.
Tray Dimensions	28" W x 20" D



## Color Choices



## Additional Features



- 1 Tray Tilt Angle**  
0 - 80 degrees

- 2 Tray Height Adjustment**  
Short: 16-21"  
Tall: 20-28"

- 3 Use with Wheelchair**  
28" wide base opening gives many wheelchair users access to the Mobile Activity Table.

- 4 Use with Soft-Touch® Sitters**  
Designed to integrate with the Soft-Touch Sitters for a variety of activities.

- 5 Use with Walker**  
28" wide base allows children with walkers access to the Mobile Activity Table from a standing position.

