

SIMPLIFIED WHEELCHAIR BY QUICKIE

The QUICKIE® QX/QXi® is a durable, low-maintenance wheelchair frame that is perfect for the person looking for simplicity, longevity, and QUICKIE's most popular options at the right price with quick adjustments and multiple options.



QX & QXi



Quick and Intuitive Adjustments

Most adjustments can be made in less than one minute, using only one tool.

- Center of Gravity
- Seat-to-Floor Height
- Camber
- Axle Spacing



QUICKIE QX Frame Option

Love QUICKIE's style, but don't qualify for an ultra lighweight wheelchair? The QX may be the option for you. The simplified QX gives you the customization, options and frame colors of a QUICKIE wheelchair in a high strength, wheelchair package.



Simple Option Mix

- · Options include the most frequently used armrests, wheels and tires
- New Swing-in/Swing-out Hangers with simple footrest adjusments and intuitive release levers
- Transit option



eDurable Frame Design

- Ultra lightweight and durable aluminum frame with a 300 lb. weight capacity
- Three-arm cross brace sits inside frame for a more rigid ride

SPECIFICATIONS	QUICKIE QX	QUICKIE QXi
HCPCS Code	K0004	K0005
Frame Style	Folding – Traditional	Folding – Traditional
Transport Weight	16.8 lbs	17.6 lbs
Overall Weight	28.0	28.8 lbs
User Weight Capacity	300 lbs	300 lbs
Seat Width	14" to 20"	14" to 20"
Seat Depth	14" to 20"	14" to 20"
Back Angle Adjustment	-15° to 3°	-15° to 3°
Front Seat-to-Floor Height	12.5" to 19.5"	12.5" to 19.5"
Rear Seat-to-Floor Height	12" to 19"	12" to 19"
Hanger Angle	60°, 70°, 80° & Elevating	60°, 70°, 80° & Elevating
Warranty	Limited Lifetime	Limited Lifetime



