



# JAY® X2 Cushion

E2622

February 2023

## Mark For:

Date: \_\_\_\_\_

Dealer Acct #: \_\_\_\_\_

Dealer: \_\_\_\_\_

Dealer Contact: \_\_\_\_\_

Dealer Address: \_\_\_\_\_

Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Dealer Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_

Confirmation Email: \_\_\_\_\_

Confirm Via: ☐ Fax ☐ Email

## Submitting for:

☐ Quote

☐ Order

PO#: \_\_\_\_\_

## ADDITIONAL SHIPPING INFORMATION

Ship To: \_\_\_\_\_

Attention: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Ship To Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

## JAY X2 Cushion - E2622

### JAY X2 CUSHION Less than 22" in Width (300lb weight capacity)

✓	Option #	HCPCS	Definition	Actual Width	Actual Depth	Price
<input type="checkbox"/>	JX21414	E2622	JAY X2 Cushion 14" X 14"	13.75"	14"	\$480.00
<input type="checkbox"/>	JX21416	E2622	JAY X2 Cushion 14" X 16"	13.75"	16"	\$480.00
<input type="checkbox"/>	JX21515	E2622	JAY X2 Cushion 15" X 15"	15.25"	15"	\$480.00
<input type="checkbox"/>	JX21516	E2622	JAY X2 Cushion 15" X 16"	15.25"	16"	\$480.00
<input type="checkbox"/>	JX21616	E2622	JAY X2 Cushion 16" X 16"	15.75"	16"	\$480.00
<input type="checkbox"/>	JX21617	E2622	JAY X2 Cushion 16" X 17"	15.75"	17"	\$480.00
<input type="checkbox"/>	JX21618	E2622	JAY X2 Cushion 16" X 18"	15.75"	18"	\$480.00
<input type="checkbox"/>	JX21620	E2622	JAY X2 Cushion 16" X 20"	15.75"	20"	\$480.00
<input type="checkbox"/>	JX21716	E2622	JAY X2 Cushion 17" X 16"	17.25"	16"	\$480.00
<input type="checkbox"/>	JX21717	E2622	JAY X2 Cushion 17" X 17"	17.25"	17"	\$480.00
<input type="checkbox"/>	JX21816	E2622	JAY X2 Cushion 18" X 16"	17.75"	16"	\$480.00
<input type="checkbox"/>	JX21818	E2622	JAY X2 Cushion 18" X 18"	17.75"	18"	\$480.00
<input type="checkbox"/>	JX21820	E2622	JAY X2 Cushion 18" X 20"	17.75"	20.25"	\$480.00
<input type="checkbox"/>	JX22020	E2622	JAY X2 Cushion 20" X 20"	19.75"	20"	\$480.00

### JAY X2 CUSHION with Extra Cover - Less than 22" in Width (300lb weight capacity)

<input type="checkbox"/>	JX21414-2	E2622	JAY X2 Cushion 14" X 14"	13.75"	14"	\$592.00
<input type="checkbox"/>	JX21416-2	E2622	JAY X2 Cushion 14" X 16"	13.75"	16"	\$592.00
<input type="checkbox"/>	JX21515-2	E2622	JAY X2 Cushion 15" X 15"	15.25"	15"	\$592.00
<input type="checkbox"/>	JX21516-2	E2622	JAY X2 Cushion 15" X 16"	15.25"	16"	\$592.00
<input type="checkbox"/>	JX21616-2	E2622	JAY X2 Cushion 16" X 16"	15.75"	16"	\$592.00
<input type="checkbox"/>	JX21617-2	E2622	JAY X2 Cushion 16" X 17"	15.75"	17"	\$592.00
<input type="checkbox"/>	JX21618-2	E2622	JAY X2 Cushion 16" X 18"	15.75"	18"	\$592.00
<input type="checkbox"/>	JX21620-2	E2622	JAY X2 Cushion 16" X 20"	15.75"	20"	\$592.00
<input type="checkbox"/>	JX21716-2	E2622	JAY X2 Cushion 17" X 16"	17.25"	16"	\$592.00
<input type="checkbox"/>	JX21717-2	E2622	JAY X2 Cushion 17" X 17"	17.25"	17"	\$592.00
<input type="checkbox"/>	JX21816-2	E2622	JAY X2 Cushion 18" X 16"	17.75"	16"	\$592.00
<input type="checkbox"/>	JX21818-2	E2622	JAY X2 Cushion 18" X 18"	17.75"	18"	\$592.00
<input type="checkbox"/>	JX21820-2	E2622	JAY X2 Cushion 18" X 20"	17.75"	20.25"	\$592.00
<input type="checkbox"/>	JX22020-2	E2622	JAY X2 Cushion 20" X 20"	19.75"	20"	\$592.00

Please visit [www.sunparts.us](http://www.sunparts.us) for spare parts information.

## JAY YOUR WAY MODIFICATIONS

Note: To order sizes offered via JAY Your Way add a "M" to the front of the part number and follow the prompts.

### Active Fluid Bladder

☐ AFB Active Fluid Bladder \$141.00

Notes: 40% Larger fluid bladder that provides fluid support to the outside edges of the cushion.

### Fluid Volume Change

☐ OF Fluid Volume Overfill NC

☐ 5% ☐ 10% ☐ 15% ☐ 20% ☐ 25% ☐ 30% ☐ 35% ☐ 40%

☐ UF Fluid Volume Underfill NC

☐ -5% ☐ -10% ☐ -15% ☐ -20%

### Obliquity Fluid Volume Change

☐ OOFR Obliquity Overfill - Right Side NC

☐ 5% ☐ 10% ☐ 15% ☐ 20% ☐ 25% ☐ 30% ☐ 35% ☐ 40%

☐ OOFI Obliquity Overfill - Left Side NC

☐ 5% ☐ 10% ☐ 15% ☐ 20% ☐ 25% ☐ 30% ☐ 35% ☐ 40%

☐ OUFU Obliquity Underfill - Right Side NC

☐ -5% ☐ -10% ☐ -15% ☐ -20%

☐ OULI Obliquity Underfill - Left Side NC

☐ -5% ☐ -10% ☐ -15% ☐ -20%

### Cover Options

☐ CVR Cover Options NC

☐ X-static w/ Stretch ☐ Incontinent ☐ Microclimatic


### JAY X2 Size Matrix

		WIDTH										
		12	13	14	15	16	17	18	19	20	21	22
DEPTH	12	N/A	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A
	13	N/A	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A
	14	N/A	JYW		JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A
	15	N/A	JYW	JYW		JYW	JYW	JYW	N/A	N/A	N/A	N/A
	16	N/A	JYW					JYW	JYW	N/A	N/A	N/A
	17	N/A	N/A	N/A	JYW			JYW	JYW	JYW	N/A	N/A
	18	N/A	N/A	N/A	JYW		JYW		JYW	JYW	N/A	N/A
	19	N/A	N/A	N/A	JYW	JYW	JYW	JYW	JYW	N/A	N/A	N/A
	20	N/A	N/A	N/A	JYW		JYW		JYW		N/A	N/A
	21	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	22	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A

	=	Standard Offering
JYW	=	Offered Via Jay Your Way
N/A	=	Not available

### JAY X2 Measuring Guide

Measurements made from completed cushion with covers.



#### Width

X = Width

Note: Measurement made from rear bottom edge to bottom edge of the cushion.  
Ex: JX21618 X=16"

#### Depth

Y = Depth

Note: Measurement made from front edge of rear cutout and front lower of cushion.  
Ex: JX21618 Y = 18"

### Size Change

☐ RL Reduce Length NC

☐ -1" ☐ -2" ☐ -3" ☐ -4"

☐ RW Reduce Width NC

☐ -1"

### Positioning Components

☐ APC Add Positioning Components NC

Notes: Add 1" tall Lateral, Medial, and Pelvic Positioning Pieces.

Don't see the option you want? NO PROBLEM! Write your request below and we can evaluate it and provide a quote.

Notes:



Sunrise Medical (US) LLC

2842 Business Park Ave. • Fresno, CA 93727 • USA

Customer Service: 800-333-4000 Fax: 800-300-7502 Email: [orders@sunmed.com](mailto:orders@sunmed.com) [www.sunrisemedical.com](http://www.sunrisemedical.com)

Sunrise Medical (US) LLC. © 2/2023 127757 REV. E