


Valiant Series



The New Valiant (Model 311N)



The Companion HD (Model 316)
America's Most Durable Adjustable Bed



The Valiant HD (Model 312) shown with the head and foot elevated, in the Trendelenburg mode

"The bed is great! Thank you so much for helping me get my life back to normal after Katrina."

Pamela Ortega
Baton Rouge, LA

Every year we receive thousands of phone calls seeking a variety of special requirements. Despite the tremendous success of the *Original*, we could not meet all of the diverse requests with one model. In 1995 we added the Transfer Master Plus. The more affordable Valiant Series later replaced the Plus in 1998.

The Valiant Adjustable Hi-Low Bed has the advantage of offering three versions for any weight requirement. The standard model handles up to **300 pounds**, the enhanced HD model can support **600 pounds**, and the super heavy duty SHD has a weight capacity of **750 pounds**.

These beds fit attractively into any home yet meet the strict requirements for many assisted living facilities. While the *Original* features an upholstered foundation, the Valiant's foundation is a powder coated MDF Board, which has been treated with an anti-microbial agent (easy to clean and affordable).

The Valiant operates with four independent CSA-US approved motors allowing the bed to rise in multiple positions such as: head and foot articulation, Hi-Low, Trendelenburg and Reverse Trendelenburg.

The Valiant's standard minimum height is 22 inches from the top of the mattress to the floor (it will raise 10 inches from this point). It is compatible with most "Hoyer" type systems. Each bed also comes with an orthopedic mattress that can be finished in cloth or vinyl.

The Companion, a head and foot adjustable bed is identical in style and a perfect compliment to the Valiant. This bed is often ordered in a Dual King with the Valiant.

Call 800-475-8122 or visit us at www.transfermaster.com

Valiant Series

The New Valiant (Model 311N)



Standard Features:

- Head, Foot and Hi-Low Adjustable
- Trendelenburg and Reverse Trendelenburg
- Weight capacity up to 300 pounds
- Deck boards are powder coated MDF board which has been treated with an antimicrobial agent
- Safety thermal overload protected motors
- Orthopedic mattress, choice of cloth or vinyl
- Lowest setting about 22 inches from top of mat to floor
- Raises up 12 inches
- Lighter, easier for delivery and set-up



Dual King with the Valiant HD on the left and the Companion on the right

The Valiant HD (Model 312)

- Enhanced version with increased weight capacity of **600 pounds**.
- Lowest setting about 22 inches from top of mat to floor
- Raises up 10 inches



Same Dual King demonstrating the Hi-Low motion of the Valiant

The Valiant SHD (Model 313)

- Our super heavy duty version, with a weight capacity of **750 pounds**.
- Lowest setting about 22 inches from top of mat to floor
- Raises up 10 inches

The Valiant Clone (Model 305)

Standard Features:


- Standard, non-adjustable bed built to match the height of the Valiant
- Orthopedic mattress, choice of cloth or vinyl

The New Companion (Model 314N)



Standard Features:

- Head and Foot Adjustable
- Looks identical to the Valiant
- Level surface with the Valiant when lowered
- Safety thermal overload protected motors
- Orthopedic mattress, choice of cloth or vinyl
- Weight capacity up to 300 pounds



The New Valiant demonstrating the Reverse Trendelenburg motion


The Companion HD (Model 316)

Enhanced version with increased weight capacity of **600 pounds**.


The Companion SHD (Model 317)

Our super heavy duty version, with a weight capacity of **750 pounds**.


The Valiant HD (Model 312) shown with optional Fracture Frame System and optional Half Rails.



Valiant and Companion Options:



Half Rails



Big Button Hand Control



Joystick Control



Hand Control Extension Cord



Magic Massage

Headboard, Footboard, Memory Foam Overlay, Hand Control Clip & Battery Back-up

Available on Valiant HD/SHD and Companion HD/SHD only:

Fracture Frame System (All the above options may be used for HD models)

DUAL KING OPTIONS:

Valiant (Model 311/312) with a Valiant (Model 311/312)

Valiant (Model 311/312) with a Companion (Model 314/316)

Valiant (Model 311/312) with a Clone (Model 305)

Companion (Model 314/316) with a Companion (Model 314/316)

Call 800-475-8122 or visit us at www.transfermaster.com