

Soft-Touch® Booster Car Seat Order Form

Online Dealer Resources: Get Quotes Instantly & Get Immediate Order Confirmations!

Please mark clearly with an **X** the items you wish to purchase.

Fax form to Special Tomato at (866) 901- 8707 or call in your order to our Pediatric Helpline at (866) 529-8407.

Account #: _____	Ship to: _____
Date: _____	Name: _____
PO #: _____	Address: _____
Buyer: _____	City: _____
Attention: _____	Phone: () _____

Soft-Touch Booster Car Seat Features:

- FAA Compliant
- Portable, lightweight and durable
- FMVSS213 compliant
- Vehicle LATCH system required
- Easy to clean
- Small Weight Capacity: 51 lbs. to 90 lbs.
- Large Weight Capacity: 81 lbs. to 130 lbs.



Gianna in the Small:
Height (45")
Weight (51 lbs)
Age (6 y)



Seat Specifications	Small	Large
Inside Seat Depth	11.5" (29 cm)	13.5" (34 cm)
Inside Seat Width	12" (30.5 cm)	14" (35.5 cm)
Seat Back Height	29" (74 cm)	35" (89 cm)
Harness Shoulder Height	14.5-18.5" (37- 47 cm)	17.75-22.25" (45 – 56.5 cm)
Overall Seat Depth	16" (41 cm)	19" (48 cm)
Overall Seat Height*	34" (86 cm)	40" (102 cm)
Overall Seat Width	16" (41 cm)	18" (46 cm)
Weight Capacity (Min)	51 lbs (23 kg)	81 lbs (37 kg)
Weight Capacity (Max)	90 lbs (41 kg)	130 lbs (59 kg)
Weight of Seat	13 lbs (6 kg)	15 lbs (7 kg)

Overall seat height: Customer **MUST measure from the seat surface to the headliner to be sure the large will fit into your vehicle.*





Online Dealer Resources: Get Quotes Instantly & Get Immediate Order Confirmations!

Please mark clearly with an X the items you wish to purchase.

Fax form to Special Tomato at (866) 901- 8707 or call in your order to our Pediatric Helpline at (866) 529-8407.

Color Options and Features



Item #	Item Description	MSRP
71010302	Size Small Special Tomato® Booster Car Seat™ - LILAC	\$1,619.34
71010304	Size Small Special Tomato® Booster Car Seat™ - AQUA	\$1,619.34
71010309	Size Small Special Tomato® Booster Car Seat™ - DARK GRAY	\$1,619.34
71110302	Small ST Booster Car Seat™ w/ Stationary Base - LILAC	\$2,448.65
71110304	Small ST Booster Car Seat™ w/ Stationary Base - AQUA	\$2,448.65
71110309	Small ST Booster Car Seat™ w/ Stationary Base - DARK GRAY	\$2,448.65
71210302	Small ST Booster Car Seat™ w/ Mobile Base - LILAC	\$2,978.98
71210304	Small ST Booster Car Seat™ w/ Mobile Base - AQUA	\$2,978.98
71210309	Small ST Booster Car Seat™ w/ Mobile Base - DARK GRAY	\$2,978.98
71310302	Small ST Booster Car Seat™ w/ Floor Sitter Base - LILAC	\$1,803.29
71310304	Small ST Booster Car Seat™ w/ Floor Sitter Base - AQUA	\$1,803.29
71310309	Small ST Booster Car Seat™ w/ Floor Sitter Base - DARK GRAY	\$1,803.29
71010402	Size Large Special Tomato® Booster Car Seat™ - LILAC	\$2,159.17
71010404	Size Large Special Tomato® Booster Car Seat™ - AQUA	\$2,159.17
71010409	Size Large Special Tomato® Booster Car Seat™ - DARK GRAY	\$2,159.17
71110402	Large ST Booster Car Seat™ w/ Stationary Base - LILAC	\$3,133.84
71110404	Large ST Booster Car Seat™ w/ Stationary Base - AQUA	\$3,133.84
71110409	Large ST Booster Car Seat™ w/ Stationary Base - DARK GRAY	\$3,133.84
71210402	Large ST Booster Car Seat™ w/ Mobile Base - LILAC	\$3,776.45
71210404	Large ST Booster Car Seat™ w/ Mobile Base - AQUA	\$3,776.45
71210409	Large ST Booster Car Seat™ w/ Mobile Base - DARK GRAY	\$3,776.45



Soft-Touch® Booster Car Seat™ Small ONLY w/ Floor Base



Soft-Touch® Booster Car Seat™ shown w/ Stationary Base



Soft-Touch® Booster Car Seat™ shown w/ Mobile Base

*Overall seat height: Customer MUST measure from the seat surface to the headliner to be sure the large will fit into your vehicle.