



Quickie GP Swing Away

January 2014

EIGPSA

Mark For:

Submitting for:

☐ Quote

☐ Order

Date: _____

PO#: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.

The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME STYLE

☐	025WM26	K0005	Quickie GP Swing-Away	\$1,995.00
--------------------------	---------	-------	-----------------------	------------

2. FRAME STYLE & SET UP

Standard Frame Colors

☐	025C140	Aztec Gold	NEW!	NC
☐	025C2	Black (Not avail as base w/Zebra)		NC
☐	025C141	Black Cherry	NEW!	NC
☐	025C10	Black Opal		NC
☐	025C38	Blue Opal		NC
☐	025C36	Candy Blue		NC
☐	025C14	Candy Purple		NC
☐	025C35	Candy Red		NC
☐	025C142	Electric Blue	NEW!	NC
☐	025C143	Evergreen	NEW!	NC
☐	025C126	Glow		NC
☐	025C144	Green Applie	NEW!	NC
☐	025C147	Hot Sparkle Pink	NEW!	NC
☐	025C148	Mauve Pink	NEW!	NC
☐	025C103	Rootbeer		NC
☐	025C114	Sparkle Silver		NC
☐	025C145	Sunrise Orange	NEW!	NC
☐	025C113	Titanium Color Paint		NC
☐	025C8	Yellow		NC

Matte Frame Colors

☐	025C116	Matte Black		NC
☐	025C151	Matte Black Cherry	NEW!	NC
☐	025C150	Matte Electric Blue	NEW!	NC
☐	025C149	Matte Evergreen	NEW!	NC
☐	025C152	Matte Purple	NEW!	NC

Custom Frame Color

☐	025C112	Desert Camouflage	\$295.00
☐	025C43	Mossy Oak Camouflage	\$295.00
☐	025C42	Stars & Stripes	\$295.00
☐	025C108	Zebra	\$295.00

Frame Width

Seat Width 1/2" Narrower

☐	025W3	12"	NA w/ adj. upholstery	NC
☐	025W4	13"	NA w/ adj. upholstery	NC
☐	025W5	14"		NC
☐	025W6	15"		NC
☐	025W7	16"		NC
☐	025W8	17"		NC
☐	025W9	18"		NC

Frame Width - Continued

☐	025W41	19"		\$125.00
☐	025W11	E2201	20"	\$450.00

Sling Depth

☐	025SD5	12"		NC
☐	025SD6	13"		NC
☐	025SD7	14"		NC
☐	025SD8	15"		NC
☐	025SD9	16"		NC
☐	025SD10	17"		NC
☐	025SD11	18"		NC
☐	025SD74	19"		\$125.00
☐	025SD13	E2203	20"	\$475.00
☐	025SD14	E2203	21"	\$475.00
☐	025SD33	E2204	22"	\$795.00

3. SEATING

Cushion

☐	025CU1	K0669	2" Cushion	\$95.00
☐	025CU2		Omit Cushion	NC
☐	025CU4		Omit Cushion and Seat Sling	-\$15.00

JAY Bundled Cushions - Optional

To order any JAY cushions on-chair, please see JAY Price List. Cushion name or part number must be filled out with desired width & depth.

☐	025JC20	JAY Bundled Cushion		Price List
		Cushion Name or Part #:		
		Width:	Depth:	

4. BACKREST

Back Type			
☐	025BT1	Folding Backrest (Adjustable)	
☐	025BA1	11" - 14" Back Height	\$65.00
☐	025BA2	12" - 16" Back Height	\$65.00
☐	025BA3	14" - 18" Back Height	\$65.00
☐	025BA4	16" - 20" Back Height	\$65.00
☐	025BT3	Non-Folding Backrest (Fixed)	
		<i>NA with push handles or adj. upholstery</i>	
☐	025BB1	8" Back Height	NC
☐	025BB2	9" Back Height	NC
☐	025BB3	10" Back Height	NC
☐	025BB4	11" Back Height	NC
☐	025BB5	12" Back Height	NC
☐	025BB6	13" Back Height	NC
☐	025BB7	14" Back Height	NC
☐	025BB8	15" Back Height	NC
☐	025BB9	16" Back Height	NC

Back Upholstery			
☐	025BU1	Standard	NC
☐	025BU8	Adjustable	\$175.00
☐	025BU16	Reinforced (NA w/ adj. upholstery)	\$30.00
☐	025BU14	Omit Back Upholstery	-\$15.00

Backrest Options			
☐	025BT33	Folding w/ Integral Push Handles	NC
☐	025A85	Bolt-On Push Handle (Folding only) <i>NA with Xtender Power Assist</i>	\$115.00
☐	025SR2	Regular Stroller Handle <i>NA with amputee axle plate or Xtender Power Assist</i>	\$125.00
☐	025SR3	Tall Stroller Handle <i>NA with amputee axle plate or Xtender Power Assist</i>	\$125.00

JAY Back (Optional)			
☐	025JB01	J3 Back (See J3 Order Form)	Price List

5. FOOTREST

Hangers			
☐	025H1	60° Swing Away	NC
☐	025H2	70° Swing Away	NC
☐	025H6	90° Swing Away <i>Only available with kids angle adjustable, adult angle adj. NA with extension tubes</i>	NC
☐	025H14	E0990 Elevating Hanger <i>Only available w/ alum or angle adj footplates.</i>	\$295.00
☐	025H24	K0053 Articulating <i>Only available w/ alum or angle adj footplates.</i>	\$350.00

Footrest Extension Tubes			
<i>See Chart on the bottom of the page for lengths</i>			
☐	025E1	Short NA w/ ELR or ALR	NC
☐	025E2	Medium	NC
☐	025E5	Long	NC
☐	025E16	K0037 High Mount <i>Only avail. with adult and kids angle adjustable footplates</i>	\$55.00
☐	025E12	Short ELR	NC
☐	025E15	Long ELR	NC

Footplates			
<i>Heel loop or leg strap not included, See Footplate Options Section</i>			
☐	025F1	Composite Footplate 14"-20" <i>NA w/ 90° hangers and ELRs</i>	NC
☐	025F3	K0040 Adult Angle Adjustable	\$185.00
☐	025F4	K0040 Kids Angle Adjustable <i>Only available on 12"-16" width</i>	\$185.00
☐	025F48	K0040 Platform Flip Up <i>Available w/ 60° and 70° hangers</i>	\$190.00
☐	025F5	Aluminum <i>Only available w/ ELRs</i>	NC
☐	025F7	Omit Footplate & Ext Tube	-\$15.00

Footrest Options			
☐	025HL1	E0951 Heel Loop	\$30.00
☐	025HL5	K0038 Leg Strap	\$35.00
☐	025HL4	Omit Heel Loop	NC
☐	025IG2	Impact Guard Neoprene (<i>Only avail. On 60° or 70° hangers</i>)	\$85.00

6. CASTERS

Caster Wheels			
☐	025CW1	3" Micro	NC
☐	025CW3	4" Low Profile Polyurethane	NC
☐	025CW31	4" Micro Lighted	\$70.00
☐	025CW5	5" Low Profile Polyurethane	NC
☐	025CW7	5" x 1.5" Semi-Pneumatic <i>only avail. for use in bottom fork holes</i>	\$85.00
☐	025CW6	5" Soft Roll Semi-Pneumatic	\$85.00
☐	025CW33	E2214 E2215 6" Pneumatic	\$75.00
☐	025CW34	6" Polyurethane	NC
☐	025CW35	6" x 1.5" Semi-Pneumatic <i>only avail. for use in bottom fork holes</i>	\$85.00
☐	025CW21	E2214 E2215 8" Pneumatic	\$75.00
☐	025CW32	8" Polyurethane	NC

Caster Wheels - Continued			
☐	025CW36	E2214 E2215 8" x 2" Pneumatic	\$145.00
☐	025CW37	E2214 K0108 8" x 2" Pneumatic w/Airless Insert	\$225.00

Caster Fork Options			
☐	025CS5	3/4" Longer Fork Stem Bolt N/A w/ Xtender Power Assist Option	NC
☐	025CS9	1-1/2" Longer Fork Stem Bolt N/A w/ Xtender Power Assist Option	NC

7. ARMREST

Single Post

☐	025AA9	E0973	Single Post Ht. Adjust. Low <i>Must be removed to fold back down</i>	\$265.00
☐	025AM1	E0973	Desk Length - Classic (10")	NC
☐	025AM2	E0973	Full Length - Classic (14")	NC
☐	025AA8	E0973	Single Post Ht. Adjust. Std <i>Must be removed to fold back down</i>	\$265.00
☐	025AM1	E0973	Desk Length - Classic (10")	NC
☐	025AM2	E0973	Full Length - Classic (14")	NC

Padded Swing Away Armrest

All Padded Swing-Away NA with amputee axle plate or Xtender Power Assist

☐	025AA19	E0973	Padded S-A Medium	\$195.00
☐	025AA16	E0973	Padded S-A Adult	\$195.00

8. AXLE & WHEEL OPTIONS

Axle Plate

N/A w/ Power Assist

☐	025AL2		Standard Adjustable	NC
☐	025AL6	E0959	Amputee Adjustable NA with padded swing-away armrests or stroller handles	\$175.00
☐	025AL13		Double Slotted Adjustable	\$60.00

Axle

☐	025AX1		Quick Release Steel	NC
☐	025AX11		Quad Quick Release	\$90.00

Rear Anti-Tip

☐	025AT1	E0971	Rear Anti-Tip	\$135.00
--------------------------	--------	-------	---------------	----------

Rear Wheel Type

☐	025RW8		22" Mag NA with Xtender or 8" x 2" casters	NC
☐	025RW9		22" Spoke NA with Xtender or 8" x 2" casters	\$200.00
☐	025RW1		24" Mag	NC
☐	025RW2		24" Spoke	NC
☐	025RW11		26" Spoke NA with Xtender or 8" x 2" casters	\$65.00

Rear Wheel Tire

☐	025RT1	E2211 E2212	Pneumatic	NC
☐	025RT2	E2211 E2213	Pneumatic Airless Insert NA with 26" wheels	\$80.00
☐	025RT3	E2218	Full Poly NA with 26" wheels	NC
☐	025RT6	E2211 E2212	High-Pressure Clincher <i>Only avail. On 24" or 26" spoke rims</i>	\$95.00

Hand Rim

NA with Xtender Power Assist

☐	025HR1		Aluminum Anodized	NC
☐	025HR4		Plastic Coated -NA on all rims	\$95.00
☐	025HR7		Omit Hand Rim	NC

Projection Hand Rim

NA with plastic coated hand rims.

☐	025HR8		8 Oblique	\$240.00
--------------------------	--------	--	-----------	----------

Wheel Lock

☐	025WL1		High Mount Push	NC
☐	025WL2		High Mount Pull	NC
☐	025WL5		High Mount Scissor NA w/ grade aids or extension handles	\$125.00
☐	025WL12		Do Not Mount	NC
☐	025WL10		Omit	NC

Extension Handle

☐	025EH1	E0961	6" Pair	\$75.00
--------------------------	--------	-------	---------	---------

Grade Aid

☐	025GA1	E0974	Grade Aid NA w/ poly, high-pressure clincher tires	\$150.00
--------------------------	--------	-------	--	----------

Seat Sling to Top of Footplate Range

60° Swing- Away	Short Ext.	12 - 15 ^{3/4}
60° Swing- Away	Medium Ext.	14 ^{1/2} - 18 ^{1/4}
60° Swing- Away	Long Ext.	17 - 20 ^{3/4}
60° Platform Flip-Up	Short Ext.	12 ^{3/8} - 15 ^{5/8}
60° Platform Flip-Up	Medium Ext.	14 ^{1/2} - 17 ^{3/4}
60° Platform Flip-Up	Long Ext.	17 - 20 ^{1/4}
70° Swing- Away	Short Ext.	11 - 14 ^{3/4}
70° Swing- Away	Medium Ext.	13 ^{1/2} - 17 ^{1/4}
70° Swing- Away	Long Ext.	16 - 19 ^{3/4}
70° Platform Flip-Up	Short Ext.	11 ^{5/8} - 15 ^{3/8}
70° Platform Flip-Up	Medium Ext.	13 ^{3/4} - 17 ^{1/2}
70° Platform Flip-Up	Long Ext.	16 ^{1/4} - 20

*All adjustment ranges may vary due to seat height restrictions. Recommended 2^{1/2}" footrest clearance to ground surface.

9. Xtender Power Assist Options

Xtender Power Assist

Rear seat height of 16.5" only. Higher seat heights available through B4M.

☐	025XT2	E0986	Xtender Power Assist 4X Multiplier NICKEL	\$6,595.00
☐	025XT3	E0986	Xtender Power Assist 4X Multiplier LITHIUM	\$7,595.00

Please select one of the following handrim options

☐	025XH3		Plastic Coated	NC
☐	025XH4		Stainless Steel	NC

Xtender - Tire Option

☐	025XW1	E2211 E2212	Pneumatic	NC
--------------------------	--------	----------------	-----------	----

Xtender - Other Options

☐	025XA1		Spare NiMH Power Pack - NICKEL	\$1,395.00
☐	025XA7		Spare NiMH Power Pack - LITHIUM	\$1,795.00

10. KEY OPTIONS & ACCESSORIES

Backpack

☐	025A1	Quickie Backpack Black	\$55.00
☐	025A2	Quickie Kids Backpack Black	\$50.00

Seat Pouch

☐	025A8	Quickie Seat Pouch Black	\$45.00
--------------------------	-------	--------------------------	---------

Crutch Holder

Can only mount on left side from rear

☐	025A53	E2207	Left Side Mount	\$195.00
--------------------------	--------	-------	-----------------	----------

Extra Leg Strap

☐	025A55	K0038	Extra Leg Strap Black	\$30.00
--------------------------	--------	-------	-----------------------	---------

Miscellaneous

☐	025A83		Tool Kit	\$45.00
--------------------------	--------	--	----------	---------

Positioning Belts

☐	025A142	E0978	Belt 1.5" Auto Buckle	\$70.00
☐	025A12	E0978	Belt 2" Aircraft Buckle	\$110.00
☐	025A13	E0978	Belt 2" Aircraft Buckle Padded	\$120.00

Spoke Guards

NA w/ Xtender Power Assist

☐	025A23	K0065	Black Spoke Guards	\$115.00
☐	025A24	K0065	Clear Spoke Guards	\$115.00

Tray Table

☐	025A31	E0950	Tray Table	\$250.00
--------------------------	--------	-------	------------	----------

Side Guard

NA w/ height adjustable arms

☐	025SG3	K0108	Plastic Regular	\$125.00
☐	025SG4	K0108	Plastic Kids	\$125.00

GP Series Rear STF Chart

*All measurements given are an estimate and have been rounded to the nearest .5"

SHORT Frame (9 1/4" Front / 8 1/4" Rear)			Rear Axle Plate Hole Position					
			Standard			Inverted		
Wheel	Rim		Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	14.5	15	15.5	16	16.5	17
22"	Mag / Spoke	Pneum / Insert	14.5	15	16	16	17	17.5
24"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17	17.5
24"	Spoke / Spinergy	Clincher	15	16	16.5	16.5	17	18
24"	Mag / Spoke / Spinergy	Pneum / Insert	15.5	16	17	17	17.5	18
26"	Spoke / Spinergy	Clincher	16	17	17.5	17.5	18	19
26"	Spoke / Spinergy	Pneumatic	16.5	17	17.5	18	18.5	19

TALL Frame (11" Front / 10" Rear)			Rear Axle Plate Hole Position					
			Standard			Inverted		
Wheel	Rim	Tire	Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	16	16.5	17.5	17.5	18	19
22"	Mag / Spoke	Pneum / Insert	16.5	17	17.5	18	18.5	19
24"	Mag / Spoke	Full Poly	17	17.5	18	18.5	18.5	19.5
24"	Spoke / Spinergy	Clincher	17	17.5	18	18.5	19	19.5
24"	Mag / Spoke / Spinergy	Pneum / Insert	17	17.5	18.5	19	19	20
26"	Spoke / Spinergy	Clincher	18	18.5	19	19.5	20	20.5
26"	Spoke / Spinergy	Pneumatic	18	19	19.5	19.5	20	21

STANDARD Frame (10" Front / 9" Rear)			Rear Axle Plate Hole Position					
			Standard			Inverted		
Wheel	Rim	Tire	Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17.5	18
22"	Mag / Spoke	Pneum / Insert	15.5	16	17	17	17.5	18
24"	Mag / Spoke	Full Poly	16	16.5	17	17	18	18.5
24"	Spoke / Spinergy	Clincher	16	16.5	17	17.5	18	18.5
24"	Mag / Spoke / Spinergy	Pneum / Insert	16	16.5	17.5	17.5	18	18.5
26"	Spoke / Spinergy	Clincher	17	17.5	18	18.5	19	19.5
26"	Spoke / Spinergy	Pneumatic	17	18	18.5	18.5	19.5	20

Maximum Dump 2 3/4"

*** Custom Frame Height Through Configurator - REAR 6 1/2 - 13" (Max frame squeeze = 3 1/2") Calculate Accordingly

GP Series Front STF Chart

*All measurements given are an estimate and have been rounded to the nearest .5"

SHORT Frame (9 1/4" Front / 8 1/4" Rear)									
Fork	Spacer	Fork Hole	3"	4" Poly	5" Poly	6" Poly	6" Pne	8" Poly	8" Pne
			4" Lighted	5" Soft Roll	6" x 1.5	6" x 1.5			8x2
4 1/4	STD	Bottom	15.50	16.00	16.50	17.00	17.00	x	x
4 1/4	3/4	Bottom	16.50	17.00	17.00	17.50	18.00	x	x
4 1/4	1 1/2	Bottom	17.00	17.50	18.00	18.50	18.50	x	x
5 1/4	STD	Top	x	x	16.50	17.00	17.00	x	x
5 1/4	STD	Bottom	16.50	17.00	17.5	18.00	18.00	19.00	20.5
5 1/4	3/4	Top	x	x	17.50	18.00	18.00	x	x
5 1/4	3/4	Bottom	17.50	18.00	18.50	18.50	19.00	19.50	20.00
5 1/4	1 1/2	Top	x	x	18.00	18.50	18.50	x	x
5 1/4	1 1/2	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50
6	STD	Top	x	x	17.50	18.00	18.00	18.50	19.00
6	STD	Bottom	17.50	18.00	18.50	19.00	19.00	19.50	20.00
6	3/4	Top	x	x	18.00	18.50	18.50	19.50	19.50
6	3/4	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50
6	1 1/2	Top	x	x	19.00	19.50	19.50	20.00	20.50
6	1 1/2	Bottom	19.00	19.50	20.00	20.50	20.50	21.00	21.50
7	STD	Top	x	x	18.50	19.00	19.00	19.50	20.00
7	STD	Bottom	18.50	19.00	19.50	20.00	20.00	20.50	21.00
7	3/4	Top	x	x	19.00	19.50	19.50	20.50	20.50
7	3/4	Bottom	19.00	19.50	20.00	20.50	20.50	21.50	21.50
7	1 1/2	Top	x	x	20.00	20.50	20.50	21.00	21.50
7	1 1/2	Bottom	20.00	20.50	21.00	21.00	21.50	22.00	22.50
8x2	STD	Top	x	x	x	x	x	x	19.00
8x2	STD	Bottom	x	x	x	x	x	x	20.00
8x2	3/4	Top	x	x	x	x	x	x	19.50
8x2	3/4	Bottom	x	x	x	x	x	x	20.50
8x2	1 1/2	Top	x	x	x	x	x	x	20.50
8x2	1 1/2	Bottom	x	x	x	x	x	x	21.00

STANDARD Frame (10" Front / 9" Rear) ADD 3/4" to each measurment above

TALL Frame (11" Front / 10" Rear) ADD 1 3/4" to each measurment above

*** Custom Frame Height through Configurator - Front 9 - 13" - Calculate Accordingly (Minimum Front STF = 15 1/2" / Max = 26")