

TiLite XSERIES Order Form

Effective July 5, 2016 (New Items in Red)

FOLDING FRAME

Page 1 of 6

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:		Name:	Ship Via:
TiLite Sales Rep:		Address:	
TiLite Information:	Promo Code:	City:	State: Zip:
Ship Date:	Order #:	Phone:	Fax:

Please Email or Fax to the TiLite Customer Service Team at orders@tilite.com or 866.586.2413 / 509.586.2413

FRAME STYLE

- AXFS1 ☐ **Aero X2 Aluminum Frame** (Full Adjustability. 265 Pound Weight Limit.).....\$2,250 K0005
- G2FS1 ☐ **2GX2 Titanium Frame Upgrade** (Full Adjustability. 265 Pound Weight Limit. Light Weight, High Strength and Vibration Dampening. Not Available with Frame Transit Tie-Down.).....\$1,145 K0108
- XFS1 ☐ **Swing Away**STD K0005
- XFS3 ☐ **Hemi Height Frame** (Available with a Front Seat Height of 13.5" - 20" (34-51cm). Not Available with Depth Adjustable Back/Frame, 12" or 13" Seat Depth or Heavy Duty Frame.)NCO
- XFS2 ☐ **Fixed Front**NCO K0005
- XHDR1 ☐ **Heavy Duty Frame** (350 Pound Weight Limit. Includes Tension Adjustable by Straps Back Upholstery. Available with 18" to 22" Seat Widths and 16" to 20" Seat Depths. Not Available with 90° front end).....\$410
- XPWR1 ☐ **Power Adaptable Frame for Frank Mobility System** (No charge when 20" Rear Seat Width **A** and 20" Seat Depth **B** are selected. Requires LiteSpeed Billet Aluminum Front Wheels and 0 degree Camber. Not available with Single-Sided Forks).....\$320
- XRFR1 ☐ **Reinforced Frame** (No charge when 20" Rear Seat Width and 20" Seat Depth are selected.)\$320

TRANSIT OPTIONS

- (Please refer to the Aero X2 Transit Tie-Down Option Brochure for important information regarding the Transit Tie-Down Options.)
- XTTDO1 ☐ **Aluminum Frame Transit Tie-Down** (Consists of four Tie-Down Brackets. Not Available with 2GX2, Heavy Duty Frame, Hemi Height Frame, Swing Away Tubular Armrests or Swing Away- Flip Back Tubular Armrests. Requires a Folding Stabilizer Bar unless a Depth Adjustable Backrest is selected. Requires minimum seat depth of 14").....\$325
- XWTB ☐ **Transport Bracket Option** (Consists of four Tie-Down Brackets. Unoccupied Use Only.) (Not Available when Swing Away frame and Heavy Duty Frame are both selected. Not available with Swing Away Tubular Armrest or Swing Away-Flip Back Tubular Armrests. Requires a minimum seat depth of 14" or 15" seat Depth with Hemi Height Frame.).....\$325

FRAME FINISHES

Natural Finishes

- XFF26 ☐ **Light Satin** (Available with Aluminum Frame only. Finished with clear powder coat.).....STD
- XFF1 ☐ **Satin** (Available with Titanium Frame only.).....STD
- XSCBT1 ☐ **Satin Finish Buffing Pads**.....\$15
- XFF2 ☐ **Polished** (Available with Titanium Frame only.).....\$695
- XFF3 ☐ **Tattoos** (Tattoos on front frame tubes. Finished in Satin (Titanium) or Light Satin (Aluminum). Click on the images below to view the Tattoos online.).....\$315

XFF31 ☐ **The King of Rock**



XFF32 ☐ **Aloha**



XFF33 ☐ **Snow Leopard**



XFF34 ☐ **Wood Grain**



High Gloss Paint

XFF8 ☐ **Super Wet Black**

XFF10 ☐ **TiLite Red**

XFF29 ☐ **Neon Pink**

Metallic Paint

XFF35 ☐ **Burgundy**

XFF15 ☐ **Tangerine Metallic**

XFF16 ☐ **Salsa Red Metallic**

XFF14 ☐ **Ocean Blue Metallic**

XFF30 ☐ **Acid Green Metallic**

XFF18 ☐ **White Pearl**

XFF27 ☐ **Titanium Grey**

Pearlescent Paint

XFF17 ☐ **Midnight Blue Pearl**

XFF20 ☐ **Deep Blue Ice**

XFF25 ☐ **Electric Plum**.....STD

Matte Paint

XFF19 ☐ **Black Lava**

XFF23 ☐ **Azure Blue**

Candy Paint

XFF22 ☐ **Ruby Red**



TiLite XSERIES Order Form

Effective July 5, 2016 (New Items in Red)

FOLDING FRAME

Page 2 of 6

FRAME MEASUREMENTS



SEAT WIDTH A Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

XSW1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19".....STD
XSW2 ☐ 20" ☐ 21" ☐ 22" (21" and 22" Seat Widths A are only available with Heavy Duty Frame.).....\$465 E2201

SEAT DEPTH B Measure from 1/2" in front of back post to the front edge of seat sling. (12"-14" Seat Depths B not available with 25" (559) and 26" (560) Rear Wheels.)

XSD1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19".....STD
XSD2 ☐ 20".....\$490 E2203
XSD3 ☐ Custom Frame Depth: _____" Frame Depth Longer Than Seat Depth (Use 1" increments with a maximum of 2".
A Custom Frame Depth with a maximum of 20" is \$475.).....NCO

FRONT SEAT HEIGHT C Measure from floor to top of seat tube at front edge of seat sling. Adjustable in 1/2" increments. (13.5" - 14" Front Seat Height C requires Hemi Height Frame Style. 20.5" - 21.5" Front Seat Height C not available with Hemi Height Frame Style.)

FRONT SEAT HEIGHT RANGES

Front Wheel Size	Front Seat Heights C	Hemi Front Seat Height C	
3" Front Wheel	14.5" - 19.0"	13.5" - 17.5"	
3" Front Wheel w/Glide™	16.5" - 18.0"	15" - 17.0"	
4" Front Wheel	15.5" - 19.5"	14.5" - 18.0"	
4" Front Wheel w/Glide™	17.0" - 18.5"	15.5" - 17"	
5" Front Wheel	16.5" - 20.0"	16" - 18.5"	
5" Front Wheel w/Glide™	17.5" - 19.0"	16.5" - 18.0"	
6" Front Wheel	17.5" - 20.5"	16.5" - 19.0"	
6" Front Wheel w/Glide™	18.5" - 19.5"	17.5" - 18"	
8" Front Wheel	19.5" - 21.5"	18.5" - 20.0"	

XFSH1 ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15"
☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17"
☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19"
☐ 19.5" ☐ 20" ☐ 20.5" ☐ 21"
☐ 21.5".....STD

REAR SEAT HEIGHT D Measure from floor to top of seat tube at back post. Adjustable in 3/8" increments.

REAR SEAT HEIGHT RANGES

Rear Wheel Size	Rear Seat Height D	
22"	13" - 18"	
24"	14" - 19"	
25"	14" - 19"	
26"	15" - 20"	

XRSH1 ☐ 13" ☐ 13.5" ☐ 14"
☐ 14.5" ☐ 15" ☐ 15.5"
☐ 16" ☐ 16.5" ☐ 17"
☐ 17.5" ☐ 18" ☐ 18.5"
☐ 19" ☐ 19.5" ☐ 20".....STD

SEAT BACK HEIGHT E/SEAT BACK/PUSH HANDLES Measure from top of back post to top of seat tube at rear of frame.

SEAT BACK TYPE

XSB13 ☐ Fixed Angle, Adjustable-Height Seat Back.....STD
XLB1 ☐ Folding Stabilizer Bar (Required with Transit Tie-Down Option unless Depth Adjustable Backrest is selected.
Includes Tension Adjustable by Strap Back Upholstery.).....\$365
XSB14 ☐ Fixed Angle, Adjustable Height Seat Back w/8-Degree Bend. (Not available with Heavy Duty Frame.).....NCO
XSB18 ☐ Folding Adjustable-Angle and Height Seat Back (Not available with Heavy Duty Frame or Ride Designs® Corbac™.).....\$215
XSB16 ☐ Depth Adjustable Back/Frame (Not Available with Hemi Height Frame, Heavy Duty Frame,
19" or 20" Seat Depths B, Swing Away Tubular or Swing Away Flip-Back Armrest. Seat Depths B of 12"-17" allow 3" of rearward
adjustment. Seat Depth B of 18" allow 2" of rearward adjustment. Requires Tension Adjustable Seat Upholstery.).....NCO

PUSH HANDLES

XSB6 ☐ Integrated Push Handles.....NCO
XSB7 ☐ Integrated Fold-Down Push Handles (Requires Tension Adjustable by Straps Back Upholstery. Not available with
Heavy Duty Frame.).....\$105
XSB9 ☐ Bolt-On Height Adjustable Push Handles (Requires Tension Adjustable by Straps Back Upholstery. Not available with
Heavy Duty Frame, Folding Adjustable-Angle and Height Seat Back, Swing Away Tubular Armrests, Swing Away - Flip Back
Tubular Armrests or Ride Designs® Corbac™.).....\$315

XSB20/XSB14 **FIXED ANGLE ADJUSTABLE SEAT BACK HEIGHT E** (11"- 15" Fixed Seat Back Height not available with Heavy Duty Frame, Swing Away Tubular Arms 12"-15", Bolt-On Height Adjustable Push Handles or Folding Stabilizer Bar.)

Push Handles (Selected on Pg. 2)	Short	Tall	Set at: (Use 1" increments)
None	<input type="checkbox"/> 11"- 15"	<input type="checkbox"/> 16"- 20"	Set at: _____
Integrated	<input type="checkbox"/> 11"- 15"	<input type="checkbox"/> 16"- 20"	Set at: _____
Integrated Fold-Down	<input type="checkbox"/> 11"- 15"	<input type="checkbox"/> 16"- 20"	Set at: _____
Bolt-On Height Adj.	-----	<input type="checkbox"/> 16"- 20"	Set at: _____

TiLite XSERIES Order Form

Effective July 5, 2016 (New Items in Red)

FOLDING FRAME

Page 3 of 6

FRAME MEASUREMENTS

XS18 FOLDING ADJUSTABLE- ANGLE AND HEIGHT SEAT BACK ③ (Not available with a Folding Stabilizer Bar or Heavy Duty Frame))

Push Handles (Selected on Pg. 2)	Short	Medium	Tall	Set at: (Use 1/2" increments)
None	<input type="checkbox"/> 11"- 13"	<input type="checkbox"/> 13.5"- 16.5"	<input type="checkbox"/> 16.5"- 19.5"	Set at: _____
Integrated	<input type="checkbox"/> 12.5"- 14.5"	<input type="checkbox"/> 13.5"- 16.5"	<input type="checkbox"/> 17"- 20"	Set at: _____
Integrated Fold-Down	<input type="checkbox"/> 15"- 17"	<input type="checkbox"/> 16"- 19"	<input type="checkbox"/> 18.5"- 20.5"	Set at: _____

CAMBER ⑦

XC8R8 ☐ 0° ☐ 2° ☐ 4° ☐ 6°..... STD

FRONT ANGLE ① Front Angle ① is measured from inside of front frame tube to top of bottom frame tube.

Aero X2 Swing Away Front End

XFA3 ☐ 70° ☐ 80° ☐ 85° (Not available with 6" Front Wheels.) ☐ 90° (Requires aluminum footplate. Not available with Heavy Duty Frame, 5" or 6" Front Wheels.)..... STD

XFA7 ☐ Elevating Leg Rest Pair (Not available with Heavy Duty Frame, 14" or narrower Seat Width A or Neoprene Impact Guards.)..... \$315 E0990

XFA0 ☐ Omit Hangers..... NCO

Aero X2 Fixed Front End

XFA6 ☐ 70° ☐ 80° ☐ 85° (85° Not available with 6" Front Wheels)..... STD

SEAT TO FOOTREST ①/FOOTRESTS Measure from front edge of seat sling to top rear of footrest.

If your Front Angle ① is 70°, 80°, OR 85° select one Seat to Footrest ① length and one Footrest type:

Standard Footrests Measurement must be at least 2" shorter than Front Seat Height ③.

☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20"

XFTR5 ☐ Flip-up Footrest..... STD

XFTR6 1-Piece Flip-up Footrest - Pivot on ☐ L OR ☐ R (Not available with Heel Loops.)..... \$160 K0045

XFTR7 ☐ 1-Piece, Over-Center Folding Footrest (Aero X2 Fixed Front Only. Not available with Heel Loops.)..... \$160 K0045

XFTR9 ☐ Depth Adjustable Flip-up Footrest..... \$100 K0045

XFTR10 ☐ Depth & Angle Adjustable Flip-up Footrest (Available for 15"- 20" Seat to Footrest ① only.)..... \$250 K0040, K0045

XFTR11 ☐ Angle Adjustable Flip-up Footrest (Available for 15"- 20" Seat to Footrest ① only.)..... \$160 K0040, K0045

High Mount Footrests (Not Available with Luggage Carriers, or Impact Guards.)

☐ 7" ☐ 7.5" ☐ 8" ☐ 8.5" ☐ 9" ☐ 9.5" ☐ 10" ☐ 10.5" ☐ 11" ☐ 11.5" ☐ 12" ☐ 12.5"

☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5"

XFTR20 ☐ High Mount Flip-up..... \$175 K0037

XFTR15 ☐ Angle Adjustable High Mount Flip-up Footrest..... \$305 K0040, K0037

XFTR16 ☐ Depth & Angle Adjustable High Mount Flip-up Footrest..... \$350 K0040, K0037

Elevating Leg Rest

☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" ☐ 20" ☐ 21"

XFTR17 ☐ Depth Adjustable Flip-up Footrest..... STD

If your Front Angle ① is 90°, select one Seat to Footrest ① length and one Footrest type:

☐ 7" ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"

XFTR5 ☐ Flip-up Footrest..... STD

XFTR9 ☐ Depth Adjustable Flip-up Footrest..... \$100 K0037

XFTR10 ☐ Depth & Angle Adjustable Flip-up Footrest..... \$250 K0040, K0037

XHLP1 ☐ Heel Loops (Not available with 13" or narrower Seat Width ① or 1-Piece Footrests.)..... \$45 E0951

CUSTOM CONFIGURATION OPTIONS

FRONT WHEELS (See Front Seat Height Ranges Chart on Page 2 for Front Seat Height ③ limitations.)

XFW24 ☐ 3"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub..... \$175

XFW14 ☐ 4"x.75" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or Glide™ Single-Sided

Suspension Fork.)..... NCO

XFW3 ☐ 4"x.75" Light-Up Micro-Caster (Not available with Power Adaptable Frame, Heavy Duty Frame or Glide™ Single-Sided

Suspension Fork.)..... \$75

XFW25 ☐ 4"x.75" LiteSpeed Billet Aluminum Wheel w/Poly Tire ☐ Silver Hub OR ☐ Black Hub..... \$135

XFW26 ☐ 4"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Assist, Power Add-on System or Heavy Duty Frame.)..... \$85

XFW27 ☐ 4"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub..... \$185

XFW6 ☐ 5"x1" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or with 90° Front Angle ①

or Glide™ Single-Sided Suspension Fork.)..... NCO

XFW28 ☐ 5"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or

90° Front Angle ①.)..... \$85

XFW29 ☐ 5"x.75" LiteSpeed Billet Aluminum Wheel w/Poly Tire ☐ Silver Hub OR ☐ Black Hub (Not available with

90° Front Angle ①.)..... \$155

XFW30 ☐ 5"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub (Not available with

90° Front Angle ①.)..... \$200

XFW9 ☐ 6"x1" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or 85° or

90° Front Angle ① or Glide™ Single-Sided Suspension Fork.)..... NCO

XFW34 ☐ 6"x1" LiteSpeed Billet Aluminum Wheel w/Poly Tire ☐ Silver Hub OR ☐ Black Hub (Not available with

85° or 90° Front Angle ① or Glide™ Single-Sided Suspension Fork.)..... \$175

XFW32 ☐ 6"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Adaptable Frame, Heavy Duty Frame,

85° or 90° Front Angle ①, or Glide™ Single-Sided Fork.)..... \$90

XFW33 ☐ 6"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub (Not available with

85° or 90° Front Angle ① or Glide™ Single-Sided Suspension Fork.)..... \$210

X2FW16 ☐ 8"x1" Plastic Wheel w/Poly Tire (Not available with Power Assist, Power Add-on, Heavy Duty Frame or Glide™ Single-Sided

Available with 70° Front Angle ① only.)..... (Not available with 12"-14" Seat Depths)..... NCO

TiLite X^{SERIES} Order Form

Effective July 5, 2016 (New Items in Red)

FOLDING FRAME

Page 4 of 6

CUSTOM CONFIGURATION OPTIONS

X2FW17 <input type="checkbox"/> 8"x1.25" Plastic Wheel w/Pneumatic Tire (Not available with Power Assist or Power Add-on, Heavy Duty Frame or Glide™ Single-Sided Suspension Fork. Available with 70° Front Angle 1 only.)- (Not available with 12"-14" Seat Depths).....	\$90	
XGL1 <input type="checkbox"/> Out-Front Glide™ Suspension Fork (Available in Black anodize only. Not available with 8" Front Wheels or Heavy Duty Frame).....	\$455	E1015
Rider's Weight: <input type="checkbox"/> 125 lbs. <input type="checkbox"/> 125-180 lbs. <input type="checkbox"/> 180-265 lbs.		
XGL2 <input type="checkbox"/> Out-Front Glide™ Single-Sided Suspension Fork (Available in Black anodize only. Not available with 6" or 8" Front Wheels or Heavy Duty Frame.).....	\$585	E1015
Rider's Weight: <input type="checkbox"/> 125 lbs. <input type="checkbox"/> 125-180 lbs. <input type="checkbox"/> 180-265 lbs.		

STANDARD REAR WHEELS (See Rear Seat Height Ranges Chart on Page 2 for Rear Seat Height **1** limitations.)

XRW14 TiLite SHADOW™.....	STD	
<input type="checkbox"/> 22"(501)		
<input type="checkbox"/> 24"(540)		
<input type="checkbox"/> 25"(559) (Requires 15" minimum Seat Depth 1 .)		
<input type="checkbox"/> 26"(590) (Requires 15" minimum Seat Depth 1 .)		

PERFORMANCE REAR WHEELS (See Rear Seat Height Ranges Chart on Page 2 for Rear Seat Height **1** limitations.)

XRW7 Spinergy SPOX <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Yellow <input type="checkbox"/> Blue <input type="checkbox"/> White	\$795	
<input type="checkbox"/> 22"(501) - Std. Hub		
<input type="checkbox"/> 24"(540) - Std. Hub		
<input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" minimum Seat Depth 1 .)		
<input type="checkbox"/> 26"(590) - Sport Hub (Requires 15" minimum Seat Depth 1 . Not available with Quad Release Axles.)		
XRW8 Spinergy LX (Available with Heavy Duty Frame only when 0 or 2 degree camber is selected.) <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Yellow <input type="checkbox"/> Blue		
<input type="checkbox"/> White <input type="checkbox"/> Go Green <input type="checkbox"/> Bright Orange <input type="checkbox"/> Flamingo Pink.....	\$995	
<input type="checkbox"/> 22"(501) - Std. Hub		
<input type="checkbox"/> 24"(540) - Std. Hub		
<input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" minimum Seat Depth 1 .)		
<input type="checkbox"/> 26"(590) - Sport Hub (Requires 15" minimum Seat Depth 1 . Not available with Quad Release Axles.)		
XRW16 Spinergy Carbon Blade (Not available with Heavy Duty Frame.).....	\$1,395	
<input type="checkbox"/> 24"(540)		
<input type="checkbox"/> 25"(559) (Requires 15" minimum Seat Depth 1 .)		
XRW17 X-Core Mag (Not available with Heavy Duty Frame, Short Tabs, Quad Release Axles, Spoke Guards or Omit Handrims Option. Requires Aluminum Silver Anodized Handrims.).....	NCO	
<input type="checkbox"/> 22"(501)		
<input type="checkbox"/> 24"(540)		

AXLES

XAXL1 <input type="checkbox"/> Stainless Quick Release.....	STD	
XAXL2 <input type="checkbox"/> Titanium Quick Release.....	\$240	
XAXL3 <input type="checkbox"/> Quad Release (Not available with Mag or 26"(590) Spinergy Spox or Spinergy LX Rear Wheels.).....	\$145	

PNEUMATIC REAR TIRES

XRTR1 <input type="checkbox"/> Treaded (Not available with 25"(559) Rear Wheels.).....	STD	
XRTR2 <input type="checkbox"/> Airless Inserts.....	\$85	E2213
XRTR3 <input type="checkbox"/> 25"(559) Primo Treaded (Standard with 25"(559) Rear Wheels.).....	STD	
XRTR4 <input type="checkbox"/> Primo (Not available with Mag Rear Wheels.).....	NCO	
XRTR8 <input type="checkbox"/> 24"(540), 25"(559) and 26"(590) Schwalbe® Marathon Plus Evolution (Not available with Mag Rear Wheels.).....	\$130	
XRTR9 <input type="checkbox"/> 24"(540) and 25"(559) Schwalbe® Slick Tires (Not available with Mag Rear Wheels.).....	\$135	
XRTR24 <input type="checkbox"/> 24"(540) Kenda Kobra Knobby Tires.....	\$80	

SOLID REAR TIRES

XRTR22 <input type="checkbox"/> 22"(501) and 24"(540) Primo Express	\$100	
XRTR23 <input type="checkbox"/> 24"(540) and 25"(559) Primo Active	\$125	

HANDRIMS

Handrims by OUT-FRONT®

XRIM1 <input type="checkbox"/> Aluminum – Silver Anodized.....	STD	
XRIM21 <input type="checkbox"/> Aluminum – Black Anodized (Not available with Mag Rear Wheels.).....	NCO	
XRIM12 <input type="checkbox"/> Titanium – Satin Finish (Not available with Mag, 20"(451) or 22"(501) Rear Wheels, or 24"(540) Knobby tires.).....	\$430	
The Natural-Fit® (Original, Larger Profile) (Not available with Mag Rear Wheels.)		
XRIM14 <input type="checkbox"/> Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$360	
XRIM15 <input type="checkbox"/> Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$405	
XRIM16 <input type="checkbox"/> Without Thumb Grip.....	\$335	
The Natural-Fit® LT (Reduced Weight, Smaller Profile) (Not available with Mag Rear Wheels.)		
XRIM17 <input type="checkbox"/> Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$390	
XRIM18 <input type="checkbox"/> Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$440	
XRIM19 <input type="checkbox"/> Without Thumb Grip.....	\$365	
XRIM10 <input type="checkbox"/> The Surge (All-In-One Oval with Gription™ Strip, Larger Profile) (Not available with Mag Rear Wheels.).....	\$375	
XRIM11 <input type="checkbox"/> The Surge LT (All-In-One Oval with Gription™ Strip, Smaller Profile) (Not available with Mag Rear Wheels.).....	\$405	
XRIM20 <input type="checkbox"/> Q-Grip (Exclusive High-Friction, Anti-Slip Coating) (Not available with Mag Rear Wheels.).....	\$315	

Additional Handrim Options

XRIM5 <input type="checkbox"/> Plastic Coated – Black (Not available with Mag Rear Wheels.).....	\$155	
XRIM6 <input type="checkbox"/> Projections - 8 Vertical Projections (Not available with 22"(501), 26"(590), Carbon Blade or Mag Rear Wheels.).....	\$345	
XRIM22 <input type="checkbox"/> Omit Handrims.....	NCO	

TiLite XSERIES Order Form

Effective July 5, 2016 (New Items in Red)

FOLDING FRAME


Page 5 of 6

CUSTOM CONFIGURATION OPTIONS

Tab Length

- XRIM7 ☐ Long Tabs.....STD
 XRIM8 ☐ Short Tabs (Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.).....\$20

WHEEL LOCKS

Wheel Locks by OUT-FRONT®  (Composite wheel locks are not always compatible with 12"-13" Seat Depth. Uni-locks will be substituted at no charge.)

- XWLK7 ☐ Compact Composite Push to Lock.....NCO
 XWLK8 ☐ Compact Composite Pull to LockNCO
 XWLK10 ☐ Extension Handles for Compact Composite Push/Pull to Lock.....\$80 E0961
 XWLK9 ☐ Composite Scissor Lock (Not available on Aero X2 Swing Away.).....\$135

Additional Wheel Locks

- XWLK3 Uni-Lock Set At: ☐ Push **OR** ☐ Pull.....\$155
 XWLK4 ☐ Extension Handles for Uni-Lock.....\$80 E0961

BACK UPHOLSTERY (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

- XBUP1 ☐ Padded Nylon (Black Only. Not available with Heavy Duty Frame, Folding Back, Folding Stabilizer Bar, Depth Adjustable Frame, Fold-Down Push Handles or Ultrasuede Backrest Stripe.).....STD
 Tension Adjustable by Straps (Included with Heavy Duty Frame. Required with Aluminum Fendered Side Guards and Fold-Down Push Handles.)
 XBUP3 ☐ Padded Nylon.....\$320 E2611
 XBUP5 ☐ AIR Breathable Upholstery (Not available with Ultrasuede or NiteLite Backrest Stripe.).....\$320 E2611
 XUPC3 Ultrasuede Backrest Stripe ☐ Black ☐ Sapphire ☐ Red Pepper ☐ Burnt Orange ☐ Ultraviolet
☐ NiteLite Reflective (Not available with Padded Nylon or AIR Breathable Back Upholstery.).....\$95
 XBUP0 ☐ Omit Back Upholstery.....NCO

SEAT UPHOLSTERY

- XSUP4 ☐ Tension Adjustable (Black only).....STD
 XSUP0 ☐ Omit Seat Upholstery.....NCO

SEAT CUSHION (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

- XCSH1 Foam Cushion ☐ 2" ☐ 3" ☐ 4".....\$105

SIDE GUARDS

- XSDG1 ☐ Cloth (Not available with Ride Designs® Corbac™ Back Support.).....\$90

Aluminum

- XSDG2 ☐ Rigid Removable ☐ Small **OR** ☐ Std.....\$185
 XSDG3 ☐ Quad Removable w/Hole ☐ Small **OR** ☐ Std.....\$200
 XSDG6 ☐ Fendered Adjustable (Requires Tension Adjustable by Straps Back Upholstery. Not available with Angle-Adjustable Fold-Down Seat Backs.).....\$315


Carbon Fiber

- XSDG7 ☐ Rigid Removable ☐ Small **OR** ☐ Std.....\$290
 XSDG8 ☐ Quad Removable w/Hole ☐ Small **OR** ☐ Std.....\$300
 XSDG11 ☐ Fendered Adjustable (Requires Tension Adjustable by Straps Back Upholstery. Not available with Angle-Adjustable Fold-Down Seat Backs.).....\$530

ARMRESTS (Desk Arms come with built-in Side Guards.)


- XARM1 Swing Away Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar with Short or Medium Seat Back Height ranges or Depth Adjustable Back.).....\$210 E0973
☐ 9" - 12" **OR** ☐ 12" - 15"

- XARM4 ☐ Swing Away - Flip Back Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar or Depth Adjustable Back.).....\$295 E0973

- XARM2 Removable Height Adjustable Desk Arm with Rigid Side Guard (Requires 14" minimum Seat Depth .).....\$315 E0973


☐ 9" to 12" **OR** ☐ 12" to 15"

☐ Standard Length Pad (10") **OR** ☐ Full Length Pad (14")

- XARM3 Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop (Requires 16" minimum Seat Depth .).....\$350 E0973


☐ 9" to 12" **OR** ☐ 12" to 15"

☐ Standard Length Pad (10") **OR** ☐ Full Length Pad (14")

- XARM8 2 Point Flip Back Desk Arm with Rigid Side Guard (Requires 15" minimum Seat Depth . Not Available with Heavy Duty, Transit Tie-Down, Transport Brackets, Bolt On Height-Adjustable Push Handles, Depth Adjustable Backrest or Black Ultrasuede Cover.).....\$315 E0973

☐ 9" to 12" **OR** ☐ 12" to 15"

☐ Standard Length Pad (10") **OR** ☐ Full Length Pad (14")

- XARM9 2 Point Flip Back Desk Arm with Rigid Side Guard and Transfer Loop (Requires 15" minimum Seat Depth . Not Available with Heavy Duty, Transit Tie-Down, Transport Brackets, Bolt On Height-Adjustable Push Handles, Depth Adjustable Backrest or Black Ultrasuede Cover.).....\$350 E0973

☐ 9" to 12" **OR** ☐ 12" to 15"

☐ Standard Length Pad (10") **OR** ☐ Full Length Pad (14")

Available Armrest Option

- XARM5 ☐ Covered in Black Ultrasuede (Not available with Desk Arms).....\$100

TiLite X^{SERIES} Order Form

Effective July 5, 2016 (New Items in Red)

FOLDING FRAME

Page 6 of 6

CUSTOM CONFIGURATION OPTIONS

SEAT BELTS

XBLT1 <input type="checkbox"/> Velcro Adjustable, Black Nylon Webbing.....	\$50	E0978
XBLT2 <input type="checkbox"/> Bodypoint® 2" Push-Button Buckle (Auto-Style.).....	\$90	E0978
XBLT4 <input type="checkbox"/> Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths A .).....	\$175	E0978
XBLT5 <input type="checkbox"/> Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths A .).....	\$165	E0978
XBL10 <input type="checkbox"/> Bodypoint® Optional Band Clamp Mounting Hardware.....	\$20	E0978

CALF STRAP

XCLF2 Velcro Adjustable.....	\$35	K0038
<input type="checkbox"/> Small (Recommended for 9"- 12" Footrest Widths)		
<input type="checkbox"/> Medium (Recommended for 11"- 15" Footrest Widths)		
<input type="checkbox"/> Large (Recommended for 14"- 17" Footrest Widths)		
<input type="checkbox"/> Xtra Large (Recommended for 16"- 19" Footrest Widths)		
XCLF3 Bodypoint® Padded Velcro Adjustable.....	\$95	K0038
<input type="checkbox"/> Medium (Recommended for 12"- 15" Footrest Widths)		
<input type="checkbox"/> Large (Recommended for 16"- 20" Footrest Widths)		
XCLF1 <input type="checkbox"/> Buckle Adjustable, Black Nylon Webbing.....	\$25	K0038

BACK PACK

XBCK3 <input type="checkbox"/> Black TiLite Back Pack.....	\$85	
--	------	--

SEAT POUCH

XPCH3 <input type="checkbox"/> Black TiLite Seat Pouch.....	\$55	
---	------	--

SPOKE GUARDS by SPOZEE.COM

XGRD1 <input type="checkbox"/> Clear with TiMotion Logo (Not available with Mag Rear Wheels.).....	\$125	K0065
--	-------	-------

ANTI-TIPS

XTIP1 <input type="checkbox"/> Rear Aluminum - Black.....	\$150	E0971
XTIP2 <input type="checkbox"/> User Friendly Flip-up.....	\$315	E0971
XTIP9 <input type="checkbox"/> Tip Assist	\$95	E0971

ACCESSORIES

XUEST1 Upper Extremity Support Tray <input type="checkbox"/> 12" to 14" OR <input type="checkbox"/> 13" to 15" Seat Widths A (Requires Armrests.).....	\$270	E0950
XLUG1 <input type="checkbox"/> Luggage Carriers (Not available with High Mount Swing Away Footrests or 1-Piece Flip-up Footrest.).....	\$135	
XGS1 <input type="checkbox"/> GRIP Solutions Lap Board.....	\$50	
XMPCT2 Neoprene Impact Guards with Ultrasuede Stripe (Not available with High Mount Footrests, 90° Hangers or Elevating Leg Rest.)		
<input type="checkbox"/> Black <input type="checkbox"/> Sapphire <input type="checkbox"/> Red Pepper <input type="checkbox"/> Burnt Orange <input type="checkbox"/> Ultraviolet <input type="checkbox"/> TiLite NiteLite Reflective.....	\$95	
XCH1 <input type="checkbox"/> Crutch Holder.....	\$210	
XHS1 <input type="checkbox"/> TiLite 9 Piece SAE Round Handle Folding Hex Key Set.....	\$25	
XAFP1 4" Logo Patch <input type="checkbox"/> Army <input type="checkbox"/> Navy <input type="checkbox"/> Air Force <input type="checkbox"/> Marines	\$95	

GENERAL

TiLite manufactures a wide variety of wheelchairs to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. TiLite wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at tilite.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, The Titanium Experts, and TiShaft are registered trademarks of TiSport, LLC and are licensed to TiLite, LLC. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. Ride Designs® and Corbac® are registered trademarks of Aspen Seating, LLC. ROHO®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. Supracor®, Stimulite® and Corbee® are registered trademarks of Supracor, Inc. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. Varilite® is a registered trademark of Cascade Designs, Inc.

QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to TiLite dealers and distributors). A TiLite Quote or a TiLite Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR