



USB
Charger
Port



NEW Oasis

INFINITY
COLLECTION

AMERICA'S #1
LIFT CHAIR BRAND

Pride[®]
Lift Chairs

The **Oasis** brings a new level of comfort to Pride's Infinity Collection. Overstuffed padding provides elegant comfort and design, while contoured seating and pocketed coils deliver individualized support and relaxation. Features include infinite-position sleep recline, zero gravity and Trendelenburg positioning, a unique quad pocket design for extra storage and an easy-to-use hand control with USB charger port.

1-800-800-8586 (U.S.) • 1-888-570-1113 (Canada)
pridemobility.com • prideliftchairs.com

Pride[®]
Mobility Products Corp.
— LIVE YOUR BEST —

NEW Oasis



Shown in DuraSoft Oat



**NEW HAND
CONTROL
WITH USB
CHARGER PORT**

DURASOFT FABRICS



EMBER



FIELD



OAT



DEEP
SKY



TIMBER

ULTRAFABRICS[®] ULTRALEATHER[™] COLORS



BUFF



FUDGE



GARNET



NILE



PECAN



RAVEN
WING



LC-580 SPECIFICATIONS:

	MEDIUM	LARGE
	Infinite-Position Zero gravity Trendelenburg Chaise Lounger 5'4"-6'0"	Infinite-Position Zero gravity Trendelenburg Chaise Lounger 5'9"-6'2"
Weight Capacity	375 lbs.	375 lbs.
Seat to Floor	20"	21"
Seat Depth	20"	22"
Seat Width	20"	22"
Top of Back to Seat	27.5"	28"
Back Style	Blown Fiber Biscuit Back	
Fully Padded Chaise	Yes	Yes
Heat and Massage	No	No
Deluxe Heat and Massage	Yes	Yes
Total Comfort Seating	No	No
Distance from the Wall	29"	30"
Overall Chair Width	37"	39"
Chair Weight	130lbs.	132lbs.
Warranty	7-Year Prorated	7-Year Prorated

FEATURES AND BENEFITS

- Contoured seating and pocketed coils offer individualized comfort and support
- Sinuous no-sag back springs combine durability with just the right amount of flexibility for personalized comfort
- Overstuffed biscuit back design provides body-conforming comfort
- Specially designed arms afford a luxurious feel and firm support
- 2" foam spans throughout the ottoman and footrest for exceptional comfort and support
- Infinite-position sleep recline, zero gravity and Trendelenburg positioning allow for endless array of comfort choices
- Industry-exclusive quad pocket design offers extra storage – easy access top pockets provide convenient stowage for your hand control and other devices
- Easy-to-use hand control with USB charger port
- Removable filling allows for custom lumbar support
- 375 lb. weight capacity
- Choice of Durasoft and optional Ultra leather[™] fabrics in an array of color choices
- Headrest and armrest covers included
- Available in medium and large sizes
- Optional deluxe heat / massage and footrest extension available

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.
The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

1-800-800-8586 (U.S.) • 1-888-570-1113 (Canada)
pridemobility.com • prideliftchairs.com

Pride[®]
Mobility Products Corp.
— LIVE YOUR BEST —