



Item #	Seat	UOM
ODYSSEYGT18CS	18" Captains Seat	1/ea



## Odyssey GT

### Full Size Scooters

- Class-leading 9.5 mph top speed and 22 mile range.
- "High/Low" switch for instant speed control.
- The Sport Style Captain's Seat in 18" width and depth features a reclining, fold-down backrest, height adjustability, swivels, and adjusts forward or backward on auto-style seat slides. (Figure A)
- Padded armrests feature easy turn cylinders for angle adjustment.
- Large, air-filled, low-profile tires on Blade style alloy rims. (Figure D)
- Full suspension on articulating frame.
- One-touch, infinite tiller angle adjustment. (Figure C)
- Convenient cup holders built into the tiller. (Figure C)
- Ergonomic delta hand controls.
- Efficient, extra bright LED lighting.
- State-of-the-art control panel with backlighting. (Figure B)
- Side marker lights, hazard lights, turn signals and reflectors for superb visibility.
- Dual rear view mirrors.
- Positioning belt included.
- Anti-tip wheels and Free-Wheel mode standard.

#### specifications

Weight Capacity	325 lbs.	Charger	5A Offboard
Top Speed	9.5 mph	Brakes	Electromagnetic
Maximum Range	22 miles	Front Wheels	11" x 3.5"
Turning Radius	47.25"	Rear Wheels	11" x 4"
Climbing Angle	6°	Base Weight	133 lbs.
Ground Clearance	3.95"	Battery Weight (pair)	66 lbs.
Dimensions	52" (L) x 26.5" (W)	Seat Weight	55 lbs.
Floor-to-Seat-Height	21.5" - 23"	Total Weight	254 lbs.
Seat Type	Sport Style Captain's	Warranty on Frame	Lifetime
Controller	PG S DRIVE 120A	Warranty on Electronics	14 months
Motors	24V x 500W x 5100 RPM	Warranty on Batteries	6 months
Batteries (pair)	12V x 50AH		