



Features

- Performance and construction **perfect for outdoor use.**
- Standard comfortable and durable **Captain's Seats** feature deep-contoured cushions and backs for maximum comfort and support.
- Captain's Seats swivel and feature a **fold-down backrest** with 25 degrees of recline and an adjustable height headrest.
- Seats have four **height adjustments**, from 27" to 30".
- Seat can be moved forward or backward using the auto-style seat slides.
- **Padded, flip-up armrests** allow for easy transfers and user access. Armrests are also removable and height and width adjustable.
- Adjustable angle, **delta tiller** allows the user to position the controls at the most comfortable angle and depth.
- **Ergonomic controls** make it easier for those with limited strength or dexterity to handle the scooter.
- **Anti-tip wheels** for greater safety.
- Additional **hand-operated drum brake** for use when extra-quick stops are necessary.
- Standard **safety package** includes headlight, horn, turn signals and dual rear view mirrors.



Dasprey™

Large Mobility Scooter

drive ActiveCare™ Mobility

Osprey™ Large Mobility Scooter

Specifications

Weight Capacity	500 lbs.
Maximum Speed	9 mph
Maximum Range	38 miles
Turning Radius	73"
Ground Clearance	7"
Overall Length	54"
Overall Width	29"
Seat Dimensions	
	20"(w) x 19"(d) x 21"(h) Captain's Seat
	22"(w) x 20"(d) x 21"(h) Captain's Seat
Floor to Seat Height	27"-30"
Motor	24V x 920W x 3700rpm
Batteries	Two 12V x 75AH (Group 24)
Battery Charger	5 Amp, Offboard
Brakes	Electromagnetic & Drum
Wheels	13" x 5" Pneumatic
Base Weight	180 lbs.
Battery Weight (pair)	94 Lbs.
Seat Weights	
20" Capt. Seat Weight	55 lbs.
22" Capt. Seat Weight	55 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Inventory Numbers

OSPREY4410MG20CS 20" Capt. Seat, Gray
 OSPREY4410MG22CS 22" Capt. Seat, Gray

Accessories (sold separately)

Carriers & Holders

- Crutch/Cane Holder
- Cylinder Tank Carrier
- Rear Basket

Carry Bags

- Armrest Bag
- Backpack

Other Accessories

- Scooter Trailer
- Scooter Cover
- Cup Holder

Warranty

Lifetime on frame, 14 months on electronic controller and drive train components, six months on batteries.

Color Options



Metallic Gray

