

S I G N A T U R E S E R I E S

# Rise



**New!**

Savannah  
PR-930



Cashmere

Indigo

Raisin

## Savannah

Blending classic styling with our unique frame design, the Savannah can be placed only 2" from the wall and still recline to the full napper position!



Rise.

Roll.

Ride.



YOUR LIFE IN MOTION



# Savannah



Savannah

Model Number	PR-930
Size	Small
Weight Capacity	350 Lb.
Recline Positions	Space Saver
Standard Fabrics	3
Fabric Options	None
Chaise Seat	NO
Back Type	Button
Head Rest Cover	NO
Arm Covers	NO
Battery Back Up	YES
Overall Width	33"
Overall Height	39"
Width Between Arms	17"
Floor to Top of Seat	19"
Seat Depth	18.5"
Seat to Top of Back	24.5"
Seat to Top of Arm	6"
Seat Angle Raised	10°
Floor to Top of Foot Rest Reclined	19.5"
Floor to Top of Back Reclined	33.5"
Back Angle Closed	105°
Back Angle Reclined	130°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	180°
Distance Required from Wall Reclined	2"
Extension reclined from Back of Seat to end of Foot Rest	39.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	1 Yr.
Warranty on Electrical Parts	1 Yr.
Pro Rated Years 4 Thru 7	NO
Warranty on Mechanical Labor	1 Yr.
Warranty on Electrical Labor	1 Yr.
Removable Back Option Available	NO
Smart Tek System Standard	YES
Customizable*	NO
Left Hand Control Option Available	NO
Extra Pocket Option Available	NO
Heat & Massage Option Available	YES
Foot Rest Extension Option Available	NO
Straight Lift Option Available	NO
User Height (suggested)	5'0" - 5'3"



501 M shown

## Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

## Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

## Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.

\* See Customize Chart for Dimension limits  
 Sizes may vary depending on fabric, filling material, upholstery and carpet thickness.  
 All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

## Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



## Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68  
 Listing # E112417  
 FNAE 80-214A

