

PX4

MORE ROOM TO ROAM



SECURE,
NFC-ENABLED
KEY FOB



FDA Class II
Medical Device*



9.3 MPH

CRUISING
SPEED



FULL

SUSPENSION



500
LBS

WEIGHT
CAPACITY



23.4 MILE
RANGE



NFC
KEY FOB

BIG PLANS REQUIRE THE BEST RIDE.

With so much to see and do, it's time to level up your outdoor mobility experience. The PX4 is here to serve as your ultimate traveling companion as it has a weight capacity of 500 pounds, a top speed of up to 9.3 mph, and an impressive 23.4-mile battery range. It's your ticket to a full day of travel or just fun in the sun.

Durable and reliable, this mobility scooter offers supreme comfort thanks to full suspension and memory foam seating in several sizes. The innovative PX4 provides room to breathe – no matter your destination. Don't be afraid to think BIG!



pridemobility.com


PRIDE[®]
MOBILITY
LIVE YOUR BEST[™]

PX4 SCOOTER

FDA Class II Medical Device*



AVAILABLE VINYL ACCENT COLORS



FEATURES

- **NEW** NFC key fob for enhanced security. Simply place or tap near the receiver to start
- 500 lb. weight capacity
- Impressive top speed of up to 9.3 mph
- Cruises up to 23.4 miles on a single battery charge
- Durable color-through shrouds
- Depth adjustable 18x18-20" memory foam seat with 6" slider with additional seat size options that include 20x18-20" and 22x20-22"
- Delta steering wheel with LCD console
- Standard mirrors and front basket to hold whatever you need for the road ahead
- Water-resistant user interface
- Full LED lighting package for maximum visibility
- Pneumatic tires mounted on split rims for excellent handling and traction
- Supplemental hand brake for additional stopping power
- Choice of five unique vinyl accent color options with a 7-year outdoor durability rating

ACCESSORIES

Only one rear accessory can be used at a time.

- 3" Single/Double Cane or Crutch Holder
- Wishbone Crutch Holder
- Cane Crutch Holder
- Forearm Crutch Holder
- Walker Holder
- Cell Phone Holder
- Oxygen Tank Holder
- Cup Holder
- Saddlebag
- Lap Belt
- Weather Cover
- Safety Flag
- Rear Basket
- Military Patches

Cup holder, rear basket and saddlebag can be ordered as a bundle.

SPECIFICATIONS

WEIGHT CAPACITY 500 lbs.

MAXIMUM SPEED¹ Up to 9.3 mph

DIMENSIONS Width: 27.05 in.³ Length: 55.9 in.

TURNING RADIUS³ 82.28 in.

RANGE PER CHARGE^{2,3} Up to 23.4 miles at 220 lbs.

TOTAL WEIGHT⁴ 201.2 lbs. heaviest piece
253.8 lbs. with seat
298.8 lbs. with batteries
351.4 lbs. with seat and batteries

WARRANTY

FRAME: Lifetime limited
ELECTRONICS: 2-year limited

DRIVETRAIN: 2-year limited
BATTERY: 6 months

*Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

1. Maximum speed varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) +10%, -5%.
2. Range per charge is established by testing conducted at maximum weight capacity, in accordance with ANSI/RESNA, WC Vol2, Section 4 & ISO 7176-4 standards. These results occur from theoretical calculations based on battery specifications and drive system performance.
3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
4. The unit's total weight may vary based on differences in battery weight by battery manufacturers.

K0899 Power Mobility Device not coded by DME PDAC. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

AUTHORIZED PROVIDER



pridemobility.com