


PHOENIX 3 WHEEL



Phoenix 3 Wheel, Compact Scooter

S35010

Red/Blue,

(Has replaced the Phantom model #'s S35001/S35002), 1/cs

FEATURES:

- For indoor and outdoor use.
- **New interchangeable color panels (red & blue) included.**
- 3 wheel model allows for 32.3" turning radius.
- Convenient, compact design makes assembly/disassembly easy.
- Easy to assemble and disassemble with 4 pieces.
- Comes with easy to remove wireless battery that can be charged separately from the scooter.
- Comes with headlight.
- **New "Finger-Pull" throttle.**
- Comfortable height adjustable, upgraded seat.
- Flip back adjustable width arms.
- Speed ranges from .5 mph to 4 mph and a cruising range of 15 miles.
- Stackable battery box will allow an optional battery for increased range up to 30 miles.
- Comes standard with anti tipppers.
- Convenient free-wheel release knob.
- **Flat free tires.**
- Large carry basket.
- Easy to adjust tiller.
- For added safety the tiller has an auto straightening function.
- Has been tested and surpassed the ANSI RESNA testing requirement.
- 300 lb. Weight Capacity.



NEW

INTERCHANGEABLE
COLOR PANELS
(RED & BLUE)
INCLUDED



STACKABLE BATTERY
DOUBLES YOUR RANGE
UP TO 20 MILES.

ADDITIONAL BATTERY
ITEM# S3BATT



EASY TO ASSEMBLE
AND DISASSEMBLE
INTO 4 PIECES



Royal Blue




Red

PHOENIX 3 Wheel

PARTS DIAGRAM

243



4. T Handle Frame Lock (S350169716)
6. Handgrips (S31003)
12. Key Switch Assembly (S31006)
Key (S3KEY)
13. Red (S31209R or B) Blue
14. Seat Post (S31049)
15. Seat Pin (S350169717)
16. Wire Harness
Front Tiller (S31048)
Rear Section (S31058)
17. Tow Handle (S31011)
24. Front Wheel Lock (S31143)
26. Front Shroud Decoration:
Red (S3501R) Right (L) Left
Blue (S3501R) Right (L) Left
33. Rear Frame (S31130)
Front Frame (S31135)
34. Motor (S31031)
Motor/Gear Assembly (S31125)
35. Rear Wheel (S31132)
38. Controller (S31008)
39. Frame Contact-Front (S31056)
Rear (S31057)

SPECIFICATIONS

ITEM #	S35010
L X W	37.8" X 20"
SADMERC CODE GROUP	GROUP 1/KO800
REAR WEIGHT	24 LBS
FRONT WEIGHT	34 LBS
SEAT WEIGHT	26 LBS
BATTERY WEIGHT	20 LBS
GROSS WEIGHT	104 LBS
DYNAMIC STABILITY INCLINE IN DEGREES	10°
MINIMUM OBSTACLE CLIMBING ABILITY	1"
SEAT DIMENSIONS	17"(D) X 17.5"(W)
SEAT TO FOOT DECK HEIGHT	17-19
DRIVE WHEEL	8" SOLID
CASTER WHEEL	8" SOLID
BATTERY (INCLUDED)	12AH, 12V X 2
CHARGER	2 AMP OFF BOARD
MAX. SPEED	4 MPH
MAX. WEIGHT LOAD	300 LBS
CRUISING RANGE	UP TO 15 MILES (UP TO 30 MILES WITH STACKABLE BATTERY)
SAFE CLIMBING ANGLE	8°
TURNING RADIUS	32.3"
GROUND CLEARANCE	2"
MOTOR	270W X 1
ANTI-TIP WHEEL	2"
CHARGING TIME, HOURS	6 HR
CONTROLLER DRIVE	70A
WARRANTY	LIFETIME ON FRAME ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)

powered mobility

drive
WWW.DRIVEMEDICAL.COM

To order call toll free: 877.224.0946