

PORTABLE RAMP SYSTEMS



Harmar Mobility produces lift and ramp products that enhance people's mobility, independence and quality of life. Products include Inside, Outside and Hybrid Auto Lifts, along with Residential and Portable Access Ramps.

HYBRID LIFT SYSTEMS

Hybrid Lifts provide the ultimate in convenience and simplicity of use by combining the best of both inside and outside lifts into unique solutions.



OUTSIDE LIFT SYSTEMS

Harmar is the leader and benchmark all others look to for ease of use, quality and performance. Models capable of fitting Class I, II or III hitches.



INSIDE LIFT SYSTEMS

The universal design of Harmar's inside lifts with lifting capacities from 200 to 400 lbs. allows installation in virtually any vehicle.



MICRO LIFT SYSTEMS

Designed to combine quality with budget, these user-friendly models feature manual, powered or combination capabilities.



Harmar Mobility
2075 47th St., Sarasota, FL 34234
Phone: 800.833.0478
Toll-Free Fax: 866.234.5680



www.harmar.com

PORTABLE RAMP SYSTEMS

SINGLE-FOLD, MULTI-FOLD, EXTENDED REACH & VEHICLE-MOUNTED RAMPS



INTRODUCING...

America's Leader in Lifts & Ramps

Harmar Mobility designs, engineers and produces American made lift and ramp products that enhance people's mobility, independence and quality of life. Our products include a complete line of Auto Lifts and Access Ramps.

Harmar Mobility products are available exclusively through a nationwide network of authorized dealers. Customers receive care from the nation's top mobility dealers. Harmar's dealers will assist in determining the best product to fit your needs and budget. It is imperative to have an authorized Harmar dealer demonstrate, install and service your Harmar lift.

Harmar's lifts and ramps are designed and built with the "Harmar difference" - they are simpler, stronger, lighter and easier to use.

Chad C. Williams
President, Harmar Mobility



Multi-Fold Ramp Shown

HARMAR UNIVERSAL SLOPE CHART

Rise (in Inches)	Solid & Solid-Fold Ramps			Multi-Fold & Extended Reach		
	Threshold Ramps	Single-Fold Ramps	Acceptable Ramp Length	Single-Fold Ramps	Multi-Fold & Extended Reach	Acceptable Ramp Length
0.5"	X	X	2'	X	X	12'
0.75"	X	X	2'	X	X	10'
1"	X	X	2'	X	X	9'
1.5"	X	X	2'	X	X	8'
2"	X	X	2'	X	X	7'
3"	X	X	2'	X	X	6'
4"	X	X	2'	X	X	5'
5"	X	X	2'	X	X	4'
6"	X	X	2'	X	X	3'
7"	X	X	2'	X	X	2'
8"	X	X	2'	X	X	16"
10"	X	X	2'	X	X	12"
11"	X	X	2'	X	X	10"
12"	X	X	2'	X	X	8"

Unoccupied chairs may be loaded up to three times this slope. Divide the recommended ramp length by three to determine minimum required ramp length for an unoccupied chair.

www.harmar.com

All Harmar Portable Ramps utilize a unique non-skid metal surface with holes designed for excellent wet or dry grip. This superior surface offers improved traction; resists dirt and wear; reduces weight; and allows water and other elements to drain, making Harmar Ramps safer and more sure-footed for both indoor and outdoor use.

RAMP SPECIFICATIONS

Width: 28.5"	Capacity: 700 lbs. (all sizes)
Length: (in feet)	Weight: (lbs.)
2'	10
3'	15
4'	20
5'	25
6'	30

Width: 28.5"	Capacity: 700 lbs. (all sizes)	
Length: (in feet)	Weight: (lbs.)	
	Per Piece	Total
6'	20	40
7'	23	45
8'	25	50
10'	35	69
12'	40	79

Width: 28.5"	Capacity: 700 lbs. (all sizes)	
Length: (in feet)	Weight: (lbs.)	
	Per Piece	Total
6'	20	40
7'	23	45
8'	25	50
10'	35	69
12'	40	79

SINGLE-FOLD RAMPS

Portable solution for smaller elevation changes.

- Available in 2' to 6' Lengths
- Simple, length-wise fold for easy transportation
- One-piece extrusion for maximum stability
- Heavy-duty, 700 Pound Capacity

MULTI-FOLD RAMPS

Portable solution for large elevation changes.

- Available in 6' to 12' Lengths
- Length and width-wise fold for easy transportation
- Patent-pending hinge allows 2-piece portability
- One-piece extrusion for maximum stability
- Heavy-duty, 700 Pound Capacity

EXTENDED REACH RAMPS

Portable solution for reaching over large vehicle bumpers.

- Available in 6' to 12' Lengths
- Similar to multi-fold ramp, only with an extended lip that bridges the bumper of an SUV or van
- Patent-pending hinge allows 2-piece portability
- One-piece extrusion for maximum stability
- Heavy-duty, 700 Pound Capacity



HITCH MOUNTED Vehicle Ramp

The AR208 Hitch-Mounted Vehicle Ramp offers the ultimate portable ramp solution.

- Compatible with all lift-gate vehicles
- No loss of cargo space
- Spring-assisted folding mechanism
- Only requires 15 lbs of force to operate
- Retain visibility through holes in ramp
- Installs in minutes with no drilling
- Heavy-duty, 700 Pound Capacity

Ramp Specifications:

Capacity:	700 lbs.
Length:	7.5'
Width:	30"
Weight:	88 lbs
Hitch Class:	I, II or III (III Standard)

