

# POSTURE-MATE™ HB

## Immersion By Design

*A unique ergonomically designed full support system for the reclining highback wheelchair*

Provides postural alignment, comfort and pressure redistribution



High Back  
Reclining  
Wheelchair  
**Not Included**

### *Indications for use:*

poor head and torso control,  
leaning, slumping,  
sliding or general  
discomfort when confined  
to the wheelchair

Velcro calf guard strap  
system for support and  
access to foot-rest functions



© 2014/1015 I-Heal Company, LLC • Patent information on file.

# POSTURE-MATE™ HB

## Immersion By Design

### KEY FEATURES:

- Accommodates postural alignment in sling-seat wheelchair styles
- Custom fit to individual wheelchair widths; 16", 18", 20"
- Contoured back with lumbar support for spinal alignment and centering
- High-specification flexible foam for pelvic immersion
- Structured seat base designed for hip and thigh support
- Removable calf guard engineered to allow access to foot-rest functions
- Contoured headrest mounts independently for comfort & alignment
- Full arm covers to protect fragile skin
- Secure straps; maintain seat, calf and headrest segments in proper position on the wheelchair base
- Velcro seams; easy detachment of system segments, secure seat section to back and calf guard, facilitate cleaning of all segments



Independently mounted contoured headrest




### STRETCH COVER DESIGN

*Bidirectional stretch cover for optimal accommodation to body contours.*

*Impervious to fluids but breathable for microclimate temperature control.*

*Wipes clean with mild soap and water or EPA approved disinfectant.*



© 2014/1015 I-Heal Company, LLC • Patent information on file.