

Backpacker PLUS

Base weight lift 155 lbs.¹

Scooter/power chair overall weight 325 lbs.

Maximum product width 27"

Maximum travel distance from the ground 26"

Motor Sealed 12 volt DC

Overall platform length 44" (41" usable surface)

Maximum product Wheelbase 35"

Warranty 3-year transferable

Optional battery pack 12 volt, on-board

Optional charger 110 volt

¹ Weight does not include any accessories

² Vehicle dependent

Can easily lift the following current Pride® Power Chairs²

- Jazzy® 600
- Jazzy® 614/614HD
- Jazzy® 1103 Ultra
- Jazzy® 1113 ATS
- Jazzy® Select™
- Jazzy® Select™ 14/14XL
- Jazzy® Select™ 6
- Jazzy® Select™ 6 Ultra
- Jazzy® Select™ Traveller
- Jet™ 3 Ultra

Can easily lift the following current Quantum® Power Chairs²

- Quantum® 600/600 XL
- Quantum® 610
- Quantum® 614
- Quantum® 1103 Ultra
- Quantum® 1121
- Quantum® R4000

Can easily lift the following current Pride® Scooters²

- Celebrity® X 3 & 4
- Go-Go® Ultra X 3 & 4
- Go-Go Elite Traveller® 3 & 4
- Go-Go Elite Traveller® Plus 3 & 4
- Go-Go Elite Traveller® Plus HD 3 & 4
- Legend® 3 & 4
- Maxima 3 & 4
- Rally® 3 & 4
- Revo® 3 & 4
- Sundancer®
- Victory® 10 3 & 4

¹ Weight does not include any accessories

² Vehicle dependent

Battery Pack

NEW, OPTIONAL MULTI-PURPOSE BATTERY BOX



- Industry first
- Industry exclusive
- Charges the mobility device from the lift as you drive

Backpacker PLUS



Backpacker PLUS

Easy to install and use. The Pride® Silver Star™ Backpacker® Plus Interior Platform Lift features a pre-assembled single piece platform that you can drive your scooter or power chair on or off from either direction. Enhanced construction and a fresh new look give the Backpacker Plus the quality you expect and deserve from Pride Silver Star Lifts.

- The quality construction and design you expect from Pride® Mobility Products
- 325 lbs. lifting capacity
- Easy-to-use 2-button hand control
- 3-year transferable warranty
- Can lift a variety of Pride® Scooters and Power Chairs

Note: The Backpacker® Plus can lift all of the listed power chairs with contoured, Synergy® or Specialty seating. Please call Customer Service when requesting information about lifting with a power positioning system. Can fit in most mini-vans with little or no vehicle modification, full-size vans require modification.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

