

NEW!



Small enough to
fit in your closet

GO-GO™

4 - WHEEL TRAVEL VEHICLE

Specifications

Model number	SC44
Front tires	2" x 7.5" (solid)
Rear tires	2" x 7.5" (solid)
Maximum speed	Up to 4 mph*
Ground clearance	1.75"
Turning radius	43"
Length	39"
Width	19"
Seat type	Foldable vinyl-covered molded foam
Seat weight	14 lbs.
Total weight	96.5 lbs.
Battery requirements	(2) 12V deep cycle (AGM or gel cell)
Battery size	12 AH
Battery weight	9 lbs. (each)
Per-charge range	Up to 10 miles*
Battery charger	Off-board, 2 amp
Drive system	24V DC motor, sealed mini-transaxle
Dual Braking System	Electronic, regenerative and electromechanical
Weight capacity	250 lbs.
Standard body colors	Red or Blue
Optional body colors	Available by special order
Warranty	1 year limited

* Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.

THE NEW 4-WHEEL GO-GO™ TRAVEL VEHICLE

The new Go-Go 4-wheel is the ideal choice for travel even if your destination includes paths that prove to be a little rough from time to time. Plus, the Go-Go 4-wheel is lightweight and easily disassembles for transport in even the smallest automobiles—no need for a lift! Go where you want... with the Go-Go 4-wheel!

Pride
Mobility Products Corp.

800-800-8586 (US)
888-570-1113 (Canada)
pridemobility.com

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

