




	Jazzy® 1121	1121HD Package
--	-------------	----------------

Weight capacity:	300 lbs.	400 lbs.
Maximum speed:	Up to 5 mph	Up to 4 mph
Ground clearance:	3.5"	3.5"
Turning radius:	20.5"	23"
Overall length:	38.5"	41"
Overall width:	23.8"	23.8"
Wheels:	Drive: 14" pneumatic	Drive: 14" solid (required)
Front Anti-tips:	6" solid	6" solid
Rear Casters:	8" solid	8" solid
Suspension:	Active-Trac®	Active-Trac®
Drive train:	Two-motor, mid-wheel drive	Two-motor, mid-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	50 amp, PG VSI controller	70 amp, PG VSI controller
Specialty controls:	N/A	N/A
Battery charger:	On-board, 5 amp (standard)	On-board, 5 amp (standard)
Per-charge range:	Up to 25 miles	Up to 25 miles
Battery requirements:	(2) 12 volt, deep cycle	(2) 12 volt, deep cycle
Battery size:	NF-22	NF-22
Transport weight		
Battery weight:	38 lbs. each	38 lbs. each
Base weight:	135.5 lbs.	144.5 lbs.
Std. seat weight:	37 lbs. (medium-back)	37 lbs. (medium-back)
Std. Seat dimensions:	18" x 18"	18" x 18"
Warranty:	Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-months limited warranty on drive motors	

Colors:



FEATURES

- Active-Trac® Suspension for enhanced performance on rough and uneven terrain, creating a smooth and stable ride
- Rear articulating caster beam for great performance and durability
- One-piece, 2" diameter tubular frame for maximum strength and enhanced performance
- Easy front access to batteries
- Easy accessible side charger port
- PG VSI integral controller

OPTIONS

- Contoured high-back
- Solid seat pan
- Solid tires
- Lap belt
- Cup holder
- Heavy-duty dust cover
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Rear basket
- Angle-adjustable footplates
- 5 amp, on-board/5 amp, off-board charger combination

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure.



- Full Platform
- Power chair lift with docking station
- Power chair lift with straps
- Backpacker



Backpacker

These lifts are compatible with the Jazzy® 1121.



Jazzy®
1121

THE MOST COMPACT AND MANEUVERABLE FULL SIZE POWER CHAIR

Jazzy® 1121

Tight-space indoor maneuverability is made easy with the Jazzy® 1121's compact overall size and mid-wheel drive design. Yet, the 1121 doesn't stop there. Outdoor performance—climbing and terrain transitions—is ensured through the use of Active-Trac® Suspension and the responsive action of the front anti-tips and rear casters. The Jazzy 1121 is the perfect power chair for your diverse and demanding lifestyle.



JAZZY® 1121
shown in Viper Blue with
medium-back seat, PG VSI
controller and foot platform

Pride
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

JZ1121SS9-05