



Quickie QXi December 2015

EIF6

Mark For:

Submitting for:

☐ Quote

☐ Order

Date: _____

PO#: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HPCPS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.

The use of **HPCPS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HPCPS	Definition	Price	✓	Option #	HPCPS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME STYLE & SETUP

Base Model

☐	196WM78	K0005	Quickie QXi (265 lbs) For 14" sizes use E1236	\$1,845.00
--------------------------	---------	-------	--	------------

Frame Style

See Seat Height Matrix on page 2

☐	196FS9		Standard Frame	NC
☐	196FS30		Hemi Frame	NC

Transit

☐	196TO38		Transit System	\$315.00
--------------------------	---------	--	----------------	----------

Frame Width

☐	196W5		14"	NC
☐	196W6		15"	NC
☐	196W7		16"	NC
☐	196W8		17"	NC
☐	196W9		18"	NC
☐	196W10		19"	NC
☐	196W11	E2201	20"	\$500.00

Seat Depth

☐	196SD7		14"	NC
☐	196SD8		15"	NC
☐	196SD9		16"	NC
☐	196SD10		17"	NC
☐	196SD11		18"	NC
☐	196SD71		19"	NC
☐	196SD13	E2203	20"	\$550.00



Quickie QXi, Ultra Lightweight Adjustable Folding Wheelchair



Quickie QXi featuring the Jay ION™ cushion

2. FRONT SEAT HEIGHT SET UP

Front Seat Height

See Front Seat Height Matrix. 12.5 & 13.0 only available w/ 60° hanger

12.5" & 13" only available w/ Alum Angle Adj Footplate

☐	196SH34	12.5"	NC	☐	196SH14	16.0"	NC
☐	196SH35	13.0"	NC	☐	196SH15	16.5"	NC
☐	196SH40	13.5"	NC	☐	196SH16	17.0"	NC
☐	196SH41	14.0"	NC	☐	196SH17	17.5"	NC
☐	196SH42	14.5"	NC	☐	196SH18	18.0"	NC
☐	196SH12	15.0"	NC	☐	196SH20	18.5"	NC
☐	196SH13	15.5"	NC	☐	196SH21	19.0"	NC
				☐	196SH23	19.5"	NC

Caster Wheels

☐	196CW3	4" Low Profile Polyurethane	NC
☐	196CW44	4" x 1.25" Semi Pneumatic	\$55.00
☐	196CW5	5" Low Profile Polyurethane	NC
☐	196CW7	5" x 1.5" Semi Pneumatic	\$85.00
☐	196CW34	6" Polyurethane	NC
☐	196CW35	6" x 1.5" Semi Pneumatic	\$85.00
☐	196CW32	8" Polyurethane	NC
☐	196CW39	8" x 1.5" Semi Pneumatic	\$85.00

Caster Fork

☐	196CF71	3.5" (standard stem bolt only)	NC
☐	196CF12	5.25"	NC
☐	196CF17	6"	NC
☐	196CF23	7"	NC

Caster Stem Bolt

☐	196CS1	Standard	NC
☐	196CS5	3/4"	NC
☐	196CS9	1 1/2"	NC

Front Seat Height Matrix

Standard Stem Bolt									
Caster	3.5" FORK				5 1/4" FORK				
	STD	HEMI	STD	HEMI	STD	HEMI	STD	HEMI	
4 Low Poly	14.5	15.0	12.5	13.0	16.0	17.0	14.0	15.0	
4 x 1.25 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5	
5 Low Poly	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5	
5 x 1.5 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5	
6 Poly	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0	
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0	
8 Poly	N/A	N/A	N/A	N/A	19.0	19.0	17.0	17.0	
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.0	16.5	17.0	
Caster	6" FORK				7" FORK				
	STD	HEMI	STD	HEMI	STD	HEMI	STD	HEMI	
4 Low Poly	17.5	18.0	15.5	16.0	18.5	19.0	16.5	17.0	
4 x 1.25 Semi Pneu	17.5	18.0	15.5	16.0	18.5	19.0	16.5	17.0	
5 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5	
5 x 1.5 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5	
6 Poly	17.5	19.0	15.5	17.0	18.5	N/A	16.5	18.0	
6 x 1.5 Semi Pneu	17.5	19.0	15.5	17.0	18.5	N/A	16.5	18.0	
8 Poly	19.5	N/A	17.5	18.0	19.5	N/A	17.5	19.0	
8 x 1.5 Semi Pneu	18.5	19.5	16.5	17.5	19.0	N/A	17.0	18.5	
Seat squeeze and seat depth can affect FSH availability									
For 3/4" and 1 1/2" Stem bolt matrix see page 5									

3. HANGER & FOOTPLATE OPTIONS

Front Hangers

Footplate mounts with Extension Tube using clamping hardware

☐	196H145	60° Swing In/Out Hanger	NC
☐	196H139	70° Swing In/Out Hanger	NC
☐	196H141	80° Swing In/Out Hanger	NC
☐	196H159	K0053 Articulating Hanger NA w/ Platform Flip Up and Alum Angle Adj Front Mount footplates.	\$350.00
☐	196H160	E0990 Elevating Hanger - Ext Mount NA w/ Platform Flip Up and Alum Angle Adj Front Mount footplates.	\$295.00
☐	196H161	E0990 Elevating Hanger - Front Mount Only available w/ front mount alum angle adj footplates and 14"-22" width	\$295.00
☐	196H144	Omit Hanger & Footplate	NC

Extension Tubes

See Lower Leg Length Matrix for range

☐	196E1	Short	NC
☐	196E2	Medium	NC
☐	196E5	Long	NC

Footplate

Heel loop not included, see Footplate Option section below

☐	196F1	Composite Footplate	NC
☐	196F3	K0040 Aluminum Angle Adj. Footplate	\$185.00
☐	196F5	Aluminum Only available w/Elevating Legrest	NC

Hanger & Footplate Options

☐	196HL12	E0951 Adjustable Heel Loop	\$35.00
☐	1968A02	Gel Cover For Hangers (Pair)	\$110.00

4. REAR SEAT HEIGHT SET UP

Rear Seat Height Matrix

Rear Wheel Size	Std Frame	Hemi Frame
20" Rear Wheel	13.0" - 16.5"	12.0" - 15.0"
22" Rear Wheel	14.0" - 17.5"	12.5" - 16.0"
24" Rear Wheel	15.0" - 18.5"	13.5" - 17.0"
25" Rear Wheel	15.5" - 19.0"	14.0" - 17.5"
26" Rear Wheel	16.0" - 19.5"	14.5" - 18.0"

Rear Seat Height

See Rear Seat Height Matrix

☐	196RH9	12"	NC
☐	196RH10	12.5"	NC
☐	196RH11	13"	NC
☐	196RH12	13.5"	NC
☐	196RH13	14"	NC
☐	196RH14	14.5"	NC
☐	196RH15	15"	NC
☐	196RH16	15.5"	NC
☐	196RH17	16"	NC
☐	196RH18	16.5"	NC
☐	196RH19	17"	NC
☐	196RH26	17.5"	NC
☐	196RH20	18"	NC
☐	196RH27	18.5"	NC
☐	196RH21	19"	NC

Seat Squeeze Matrix

Frame Depth	Max Seat Squeeze
14"	2.0"
15" - 16"	2.5"
17" - 18"	2.5"
19" - 20"	3.0"

Rear Tire & Wheel Matrix

Tire	Wheel	Mag	Spoke
Pneumatic		20,22,24,26	20,22,24,26
Pneumatic w/ Airless Insert		20,22,24,26	20,22,24,26
Full Poly Tire		20,22,24	20,22,24
High Pressure Clincher		N/A	22,24,25,26
22" = 489 mm 24" = 540 mm			

Rear Wheel Size

See Rear Tire & Wheel Matrix for availability

☐	196RW3	20" Mag Wheel	NC
☐	196RW8	22" Mag Wheel	NC
☐	196RW1	24" Mag Wheel	NC
☐	196RW90	20" Lite Spoke Wheel	\$200.00
☐	196RW91	22" Lite Spoke Wheel	\$200.00
☐	196RW92	24" Lite Spoke Wheel	NC
☐	196RW100	25" Lite Spoke Wheel	\$200.00
☐	196RW101	26" Lite Spoke Wheel	\$200.00
☐	196RW25	Omit	NC

Rear Wheel Tire

See Rear Tire & Wheel Matrix for availability

☐	196RT1	E2211 E2212	Pneumatic	NC
☐	196RT2	E2211 E2213	Pneumatic w/ Airless Insert	\$80.00
☐	196RT3	E2218	Full Poly	NC
☐	196RT6	E2211 E2212	High-Pressure Clincher	\$95.00

Handrim

Natural Fit only available w/ 22, 24, 25 Lite Spoke. Plastic NA w/ 25, 26"

☐	196HR1	Aluminum Anodized	NC
☐	196HR4	Plastic	\$95.00
☐	196HR22	Natural Fit - Standard Grip	\$375.00
☐	196HR23	Natural Fit - Super Grip	\$410.00

Wheel Lock

☐	196WL1	High Mount Push	NC
☐	196WL2	High Mount Pull	NC
☐	196WL30	Compact	\$125.00
☐	196WL10	Omit Wheel Lock	NC
☐	196WL12	Do Not Mount	NC

Extension Handle

Ext. handles only available with high mount push or pull wheel locks

☐	196EH7	E0961	Extension Handle (Pair)	\$65.00
--------------------------	--------	-------	-------------------------	---------

Axle

☐	196AX1	Quick Release Steel	NC
☐	196AX11	Quad Quick Release Steel NA on Spoke	\$125.00

5. ARMREST

Armrest - Padded Swing Away

Some limitations when used with the Angle Adjustable Back in the rear most setting

Please circle a setting, will default to 11" if no selection is made

☐	196AA18	E0973	Padded Swing-Away STD (9" to 13")	\$210.00
Setting ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"				

Armrest - Height Adjustable

☐	196AA104	E0973	Single Post Ht. Adj - Std (9.75"-14.75")	\$285.00
☐	196AA105	E0973	Single Post Ht Adj - Low (8.25"-12.25")	\$285.00
☐	196AA65	E0973	Flip Back Ht Adj LITE (10"-15")	\$295.00
☐	196AA14		Omit Armrest	NC
☐	196AA103	E0973	Dual Post Flip Back Ht Adj	\$295.00

Armpad - For Height Adjustable

☐	196AM23	Desk Length (10")	NC
☐	196AM24	Full Length - Full (14")	NC

Armrest Options

☐	1968A01	Gel Cover For Swing Away Arms (pair)	\$110.00
--------------------------	---------	--------------------------------------	----------

6. BACKREST

Back Type			
☐	196BT10	Standard	NC
☐	196BT100	Angle Adjustable Back	\$165.00
Backrest Height			
☐	196B5	Med (13"-17")	NC
☐	196B9	Tall (16"-20")	NC
Push Handles			
☐	196PH6	Integral Push Handles 8° for Std Backrest, Straight for Angle Adj. Backrest	NC
☐	196PH5	Fold-Down Push Handles Only available with Adjustable Upholstery. Used in combination w/Angle Adjust only avail in 16"-20"	\$100.00
☐	196PH2	Omit Push Handles 8° for Std Backrest, Straight for Angle Adj. Backrest	NC

Back Upholstery			
☐	196BU1	Standard Upholstery	NC
☐	196BU70	E2611 Backrest Tension Adj Ballistic Nylon *Reference 186BU70 for billing	\$315.00
☐	196BU71	E2611 Backrest Tension Adj 3DX Vented *Reference 186BU71 for billing	\$330.00
☐	196BU14	Omit Back Upholstery	-\$15.00

JAY J3 Back			
☐	196JB01	JAY J3 Back	See OF

Military Service Patches			
Only available on tension adjustable upholstery and Jay 3 Back			
☐	196VU1	US Marines Patch	\$90.00
☐	196VU2	US Army Patch	\$90.00
☐	196VU3	US Navy Patch	\$90.00
☐	196VU4	US Air Force Patch	\$90.00
☐	196VU12	US Coast Guard	\$90.00
☐	196VU5	Purple Heart	\$90.00
☐	196VU8	Purple Heart with Marines Patch	\$145.00
☐	196VU9	Purple Heart with Army Patch	\$145.00
☐	196VU10	Purple Heart with Navy Patch	\$145.00
☐	196VU11	Purple Heart with Air Force Patch	\$145.00
☐	196VU13	Purple Heart with Coast Guard	\$145.00

7. SEATING

Seat Sling			
☐	196SL23	Seat Sling Depth Adjustable	NC
☐	196SL6	Omit Seat Sling	-\$15.00

JAY Cushions (Optional)			
To order any JAY cushions on-chair, please see the JAY price list. Cushion name or part number must be filled out with desired width and depth			
☐	196JC20	JAY Bundled Cushion Cushion Name or Part # Width: _____ Depth: _____	Price List

8. Key Options & Accessories

Anti Tips			
☐	196AT1	E0971 Anti Tips (Pair)	\$140.00

Side Guards			
☐	196SG33	K0108 Plastic - Standard (8.5" High)	\$125.00
☐	196SG29	K0108 Plastic - Low (6.5" High)	\$125.00

Pack Accessories			
☐	196A1	Backpack - Black/Charcoal Adult	\$55.00
☐	1961A01	Seat Pouch Black	\$45.00

Positioning Belt			
☐	196A141	E0978 1.5" Auto Buckle w/ Clamps	\$70.00
☐	196A804	K0108 Positioning Belt Slide Clamps	\$50.00

Spoke Guards			
☐	1964A01	K0065 Quickie Black Spoke Guard	\$115.00

Miscellaneous			
☐	1968A03	Caddy	\$45.00
☐	196A2100	E2207 Crutch Holder	\$195.00

10. COLOR

Standard Frame Colors			
☐	196C140	Aztec Gold	NC
☐	196C2	Black <i>(Not avail as base w/Zebra)</i>	NC
☐	196C141	Black Cherry	NC
☐	196C10	Black Opal	NC
☐	196C38	Blue Opal	NC
☐	196C36	Candy Blue	NC
☐	196C14	Candy Purple	NC
☐	196C35	Candy Red	NC
☐	196C142	Electric Blue	NC
☐	196C143	Evergreen	NC
☐	196C126	Glow	NC
☐	196C144	Green Apple	NC
☐	196C147	Hot Sparkle Pink	NC
☐	196C148	Mauve Pink	NC
☐	196C103	Rootbeer	NC
☐	196C114	Sparkle Silver	NC
☐	196C145	Sunrise Orange	NC
☐	196C113	Titanium Color Paint	NC
☐	196C8	Yellow	NC

Matte Frame Colors			
☐	196C116	Matte Black	NC
☐	196C151	Matte Black Cherry	NC
☐	196C150	Matte Electric Blue	NC
☐	196C149	Matte Evergreen	NC
☐	196C152	Matte Purple	NC

Custom Frame Color			
☐	196C112	Desert Camouflage	\$295.00
☐	196C43	Mossy Oak Camouflage	\$295.00
☐	196C42	Stars & Stripes	\$295.00
☐	196C108	Zebra <i>(Choose above base color) Std with White unless above base paint is chosen NA w/ Blk, Matte Black, Blk Opal, Titanium)</i>	\$295.00

Additional Options Through Built for Me

To inquire about any unique customization for Quickie products, contact: Phone: (800) 333 - 4000, Fax: (800) 333 9011




Front Seat Height Matrix Cont.

3/4" STEM								
Caster	3.5" FORK				5 1/4" FORK			
	STD	HEMI	STD	HEMI	STD	HEMI	STD	HEMI
4 Low Poly	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
5 Low Poly	N/A	N/A	N/A	N/A	17.0	18.0	14.5	15.5
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	15.5	16.5
6 Poly	N/A	N/A	N/A	N/A	17.5	18.5	15.5	16.5
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	15.5	16.5
8 Poly	N/A	N/A	N/A	N/A	19.5	19.5	17.5	17.5
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	19.0	19.5	17.0	17.5
Caster	6" FORK				7" FORK			
	STD	HEMI	STD	HEMI	STD	HEMI	STD	HEMI
4 Low Poly	18.0	18.5	16.0	16.5	19.0	19.5	17.0	17.5
4 x 1.25 Semi Pneu	18.0	18.5	16.0	16.5	19.0	19.5	17.0	17.5
5 Low Poly	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0
5 x 1.5 Semi Pneu	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0
6 Poly	18.0	19.5	16.0	17.5	19.0	19.5	17.0	18.5
6 x 1.5 Semi Pneu	18.0	19.5	16.0	17.5	19.0	19.5	17.0	18.5
8 Poly	N/A	N/A	18.0	18.5	19.5	19.5	18.0	19.5
8 x 1.5 Semi Pneu	19.0	19.5	17.0	18.0	19.5	19.5	17.5	19.0
Seat squeeze and seat depth can affect FSH availability								

1 1/2" STEM								
Caster	3.5" FORK				5 1/4" FORK			
	STD	HEMI	STD	HEMI	STD	HEMI	STD	HEMI
4 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	15.5	16.5
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	15.5	16.5
5 Low Poly	N/A	N/A	N/A	N/A	18.0	19.0	16.0	17.0
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.0	19.0	16.0	17.0
6 Poly	N/A	N/A	N/A	N/A	18.5	19.5	16.5	17.5
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.5	16.5	17.5
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	18.5	18.5
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	19.5	19.5	17.5	18.0
Caster	6" FORK				7" FORK			
	STD	HEMI	STD	HEMI	STD	HEMI	STD	HEMI
4 Low Poly	19.0	19.5	17.0	17.5	N/A	N/A	18.0	18.5
4 x 1.25 Semi Pneu	19.0	19.5	17.0	17.5	N/A	N/A	17.5	18.5
5 Low Poly	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5
5 x 1.5 Semi Pneu	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0
6 Poly	18.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0
6 x 1.5 Semi Pneu	18.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0
8 Poly	N/A	N/A	18.5	19.0	N/A	N/A	18.5	19.5
8 x 1.5 Semi Pneu	N/A	N/A	18.0	19.0	N/A	N/A	18.5	19.5
Seat squeeze and seat depth can affect FSH availability								

Lower Leg Length Range Matrix

Extension Mount				
Hanger	Ext Tube	Footplate		
		Composite	Alum Angle Adj	Aluminum
60° Swing In/Out	Short	11.25" - 14.75"	9.25" - 13.25"	NA
	Medium	11.25" - 16.00"	9.25" - 14.75"	NA
	Long	12.50" - 16.00"	10.75" - 16.25"	NA
70° Swing In/Out	Short	12.00" - 15.75"	10.00" - 14.50"	NA
	Medium	12.00" - 16.50"	10.75" - 16.50"	NA
	Long	13.75" - 16.50"	12.25" - 16.50"	NA
80° Swing In/Out	Short	10.50" - 15.25"	10.50" - 15.25"	NA
	Medium	11.75" - 16.50"	11.75" - 16.50"	NA
	Long	13.50" - 16.50"	11.75" - 16.50"	NA
All values are ± 0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights.				



Hanger	Extension	Seat Sling to Top of Footrest		Min Allowable FSH
		Min	Max	
ALR	Short	9 1/2	15 1/2	11 1/2
	Med	11	17 1/2	13
	Long	13	19 1/2	15
ELR Ext. Mount	Short	9	15	11
	Med	9	17	11
	Long	10	19	12
ELR Front Mount	No Ext (top hole)	4 1/2	8 1/2	6 1/2
	2" ext tube	8	11	10
	4" ext tube	10	13	12
	6" ext tube	12	15	14

