



Quickie 2 July 2014

EIQ2N

Mark For:**Submitting for:**☐ Quote☐ Order

Date: _____

PO#: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____

ST: _____

ZIP: _____

Dealer Phone: () _____

Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email**ADDITIONAL SHIPPING INFORMATION**

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____

ST: _____

ZIP: _____

Ship To Phone: () _____

Fax: () _____

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.

The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME STYLE

E1236 coding still applies for EIQ2N with certain sizes

☐	184WM32	K0005	Quickie 2	\$2,050.00
--------------------------	---------	-------	-----------	------------

2. XTENDER POWER ASSIST OPTION

Rear axle position available 15.5"-17.5". Front seat height must be greater than rear and NA w/ Rounded frame. NA with 3/4" & 1 1/2" caster stems - OR - Depth Adj. Back, HD Package or Transit Option. Rear seat heights of 19"-21" available through Built 4 Me.

Xtender Power Assist

☐	184XT2	E0986	Xtender Power Assist 4X Multiplier w/ Pneumatic Tires - NICKEL	\$6,595.00
☐	184XT3	E0986	Xtender Power Assist 4X Multiplier w/ Pneumatic Tires - LITHIUM	\$7,595.00

Xtender Handrim Options

☐	184XH3		Plastic Coated	NC
☐	184XH4		Stainless Steel	NC

Xtender Accessory

☐	184XA1		Spare NiMH Power Pack - NICKEL	\$1,395.00
☐	184XA7		Spare NiMH Power Pack - LITHIUM	\$1,795.00

Xtender Wheel Options

☐	184XA2	E2211 E2212	SPARE 24" Lite Spoke Pneu Alum	NC
☐	184XA4	E2211 E2213	SPARE 24" Lite Spoke Pneu Plastic	\$90.00

3. FRAME SET UP**Rear Frame Type (265 lbs. Weight Capacity)**

☐	184RS1		Rear Side Frame Straight	NC
☐	184RS2		Rear Side Frame Rounded NA w/ Xtender Power Assist	NC
☐	184RS3		Depth Adjustable Back (DAB) NA w/ swing-away, flip-back armrests, or Xtender. Backrest set 2" rear of seat tube w/ 1" adjustable forward.	\$195.00

Heavy Duty Package (350 lbs. weight capacity)

☐	184HV1		Heavy Duty (18"-22" w. only) NA w/ Xtender or Small Frame. Includes Reinforced Back Upholstery & Stabilizer Bar.	\$500.00
--------------------------	--------	--	---	----------

Frame Width

☐	184W3		12"	NC
☐	184W4		13"	NC
☐	184W5		14"	NC
☐	184W6		15"	NC
☐	184W7		16"	NC
☐	184W8		17"	NC
☐	184W9		18"	NC
☐	184W10		19"	NC
☐	184W11	E2201	20"	\$500.00

Frame Width (Continued)

☐	184W12	E2201	21" (Only Avail. w/ Heavy Duty)	\$500.00
☐	184W13	E2201	22" (Only Avail. w/ Heavy Duty)	\$500.00

Seat Depth

☐	184SD5		12"	NC
☐	184SD6		13"	NC
☐	184SD7		14"	NC
☐	184SD67		15"	NC
☐	184SD9		16"	NC
☐	184SD69		17"	NC
☐	184SD70		18"	NC
☐	184SD71		19"	\$75.00
☐	184SD13	E2203	20"	\$550.00

Frame Length

Hemi frame lengths NA w/ Xtender Power Assist

☐	184FL83		Small (12"-14")	NC
☐	184FL9		Regular (13"-16")	NC
☐	184FL18		Long (15"-18")	NC
☐	184FL25		X-Long (17"-20")(Steel Reinforced)	NC
☐	184FL17		Hemi (13"-16")	NC
☐	184FL24		Long Hemi (15"-18")	NC
☐	184FL29		X-Long Hemi (17"-20")(Steel Reinforced)	NC

Customer Service: 800-333-4000 Fax: 800-300-7502 Email: orders@sunmed.com www.sunrisemedical.com

Specifications are subject to change without notice. Copyright © 2012 Sunrise Medical (US) LLC. (11.05) Quickie 2 MK-100059 REV B

4. FRONT SEAT HEIGHT SET UP

Front Seat Height

Please see matrix on Page 9 for additional information

☐	184SH40	13.5" (Only available on Hemi Frame)	NEW!	NC
☐	184SH41	14" (Only available on Hemi Frame)		NC
☐	184SH42	14.5"		NC
☐	184SH12	15"		NC
☐	184SH13	15.5"		NC
☐	184SH14	16"		NC
☐	184SH15	16.5"		NC
☐	184SH16	17"		NC
☐	184SH17	17.5"		NC
☐	184SH18	18"		NC
☐	184SH20	18.5"		NC
☐	184SH21	19"		NC
☐	184SH23	19.5"		NC
☐	184SH24	20"		NC
☐	184SH26	20.5"		NC
☐	184SH27	21"		NC
☐	184SH29	21.5" (NA with Hemi Frame)		NC
☐	184SH36	22" (NA with Hemi Frame)		NC

Caster Wheels

All 8" casters NA w/ 90° hangers, ELRs or small frame length. Also NA w/ combination of regular front frames, composite footplates, 24" wheels, and 70 and 80 degree hangers. See matrix on page 9.

☐	184CW1	3" Micro	NC
☐	184CW3	4" Low Profile Polyurethane NA w/ HD package	NC
☐	184CW31	4" Micro Lighted	\$70.00
☐	184CW72	4" x 1" Aluminum Polyurethane	\$100.00
☐	184CW44	4" x 1.25" Semi Pneumatic	\$55.00
☐	184CW64	4" x 1.5" Aluminum Soft Roll	\$175.00
☐	184CW5	5" Low Profile Polyurethane	NC
☐	184CW51	5" Micro Lighted	\$85.00
☐	184CW65	5" x 1" Aluminum Polyurethane	\$110.00
☐	184CW7	5" x 1.5" Semi-Pneumatic	\$85.00
☐	184CW6	5" Soft Roll (Grey) NA w/ Frog Legs	\$85.00
☐	184CW66	5" x 1.5" Aluminum Soft Roll	\$185.00
☐	184CW34	6" Polyurethane	NC
☐	184CW33	E2214 E2215 6" Pneumatic	\$75.00
☐	184CW35	6" x 1.5" Semi Pneumatic	\$85.00
☐	184CW68	6" x 1.5" Aluminum Soft Roll	\$195.00
☐	184CW16	7" Polyurethane	NC
☐	184CW32	8" Polyurethane	NC
☐	184CW21	E2214 E2215 8" Pneumatic	\$75.00
☐	184CW39	8" x 1.5" Semi Pneumatic	\$85.00
☐	184CW36	E2214 E2215 8" x 2" Pneumatic	\$145.00
☐	184CW22	E2214 K0108 8" x 2" Pneumatic w/Airless Insert	\$225.00


Caster Options

☐	184OC1	K0073	Caster Pin Locks (N/A w/8"x2") NA w/ Xtender, Frog Legs or Transit	\$195.00
--------------------------	--------	-------	---	----------

Caster Fork

Please see matrix on Page 9 for additional information

☐	184CF46		Multi Position Fork-Short	NC
☐	184CF47		Multi Position Fork- Tall	NC
☐	184CF62		3" Fork	NC
☐	184CF12		5 1/4" Fork	NC
☐	184CF17		6 " Fork	NC
☐	184CF23		7" Fork	NC
☐	184CF29		8" x 2" Fork	NC
☐	184CF53	E1015	Frog Legs Suspension (0-100 lbs)	\$425.00
☐	184CF54	E1015	Frog Legs Suspension (101-190 lbs)	\$425.00
☐	184CF55	E1015	Frog Legs Suspension (191-265 lbs)	\$425.00



5. HANGER & FOOTPLATE OPTIONS

NEW!

Hangers

Front mount hangers only available w/ aluminum angle adj footplates.
ALR/ELR Hangers NA w/ Push To Lock wheel locks.

☐	184H145	60° Swing In/Out Hanger	NC
☐	184H139	70° Swing In/Out Hanger	NC
☐	184H141	80° Swing In/Out Hanger	NC
☐	184H154	60° Swing In/Out Hanger Hemi	NC
☐	184H155	70° Swing In/Out Hanger Hemi	NC
☐	184H140	70° Swing In/Out Front Mount	NC
☐	184H142	80° Swing In/Out Front Mount	NC
☐	184H158	90° Swing In/Out Front Mount	NC
☐	184H159	Articulating Hanger NA w/ Platform Flip Up and Alum Angle Adj Front Mount footplates.	\$350.00
☐	184H160	Elevating Hanger - Ext Mount NA w/ Platform Flip Up and Alum Angle Adj Front Mount footplates.	\$295.00
☐	184H161	Elevating Hanger - Front Mount Only available w/ front mount alum angle adj footplates and 14"-22" width	\$295.00
☐	184H32	Omit Hanger	-\$15.00

Footplates

Heel loop or leg strap not included, See Footplate Options Section

☐	184F1	Composite (13"-22"w)	NC
☐	184F271	K0040 E0951 Composite Angle Adj-(16"-22"w)	\$150.00
☐	184F3	K0040 Alum Angle Adj - Ext Mt (12"-22"w)	\$185.00
☐	184F8	K0040 Alum Angle Adj - Front Mt(12"-22"w)	\$185.00
☐	184F5	Aluminum Footplate - Fixed Available w/ Ext Mount ELR or ALR	NC
☐	184F48	K0040 Platform Flip Up (12"-22"w) NA w/ Hemi Hangers, ELR/ALR or Front Mount	\$190.00
☐	184F7	Omit Footplate & Ext Tube	-\$15.00

Footrest Extension Tubes

Required w/ extension tube mounted hangers. See leg length matrix

☐	184E1	Short	NC
☐	184E2	Medium	NC
☐	184E5	Long	NC
☐	184E16	High Mount NA w/ one-arm drive option or ELR/ALR or 90° front mount. Only available w/angle adj footplates	\$55.00
☐	184E23	2" - Front Mount Only NA w/ 6", 7" or 8" casters with 90° hangers, or ELR's	\$55.00
☐	184E24	4" - Front Mount Only NA w/ 90° hangers or ELR/ALR	\$55.00

Footplate Options (Ext Tube & Front Mount)

☐	184HL1	E0951 Heel Loop	\$30.00
☐	184HL5	K0038 Leg Strap	\$35.00
☐	184HL4	Omit Heel Loop	NC
☐	184HL6	Omit Leg Strap	NC

Lower Leg Length Setting (Extension Mount Only)

Recommend 2" less than FSTFH. (*)denotes frame

☐	184LO01	9.0" Hemi Only*	NC	☐	184LO13	15.0"	NC
☐	184LO02	9.5" Hemi Only*	NC	☐	184LO14	15.5"	NC
☐	184LO03	10.0"	NC	☐	184LO15	16.0"	NC
☐	184LO04	10.5"	NC	☐	184LO16	16.5"	NC
☐	184LO05	11.0"	NC	☐	184LO17	17.0"	NC
☐	184LO06	11.5"	NC	☐	184LO18	17.5"	NC
☐	184LO07	12.0"	NC	☐	184LO19	18.0"	NC
☐	184LO08	12.5"	NC	☐	184LO20	18.5"	NC
☐	184LO09	13.0"	NC	☐	184LO21	19.0"	NC
☐	184LO10	13.5"	NC	☐	184LO22	19.5"	NC
☐	184LO11	14.0"	NC	☐	184LO23	20.0"	NC
☐	184LO12	14.5"	NC	☐	184LO24	20.5" Std. Only*	NC

Lower Leg Length Matrix

2" & 4" extensions for front mount hangers shipped separate, not mounted. Front seat height must be at least 2" taller than lower leg length on all swing in/out hangers. Front seat height must be at least 2" taller than maximum lower leg length on all front mount hangers.

Hanger	Extension	Seat Sling to Top of Footrest in inches		*Minimum Allowable FSH
		Min	Max	
60° Hemi Swing In/Out	Short Hemi	10	13	14
	Short	-	-	-
	Medium	12.5	15.5	15.5
	Long	14.5	17.5	17.5
70° Hemi Swing In/Out	Short Hemi	13	17	14
	Short	-	-	-
	Medium	13.5	18.5	16
	Long	15.5	20.5	18
60° Swing In/Out	Short	11.5	15	15
	Medium	11.5	17	15
	Long	14	18.5	17
70° Swing In/Out	Short	13	18	15
	Medium	13	20	15.5
	Long	16	20.5	17.5
80° Swing In/Out	Short	13.5	17	15.5
	Medium	13.5	19	15.5
	Long	15.5	20.5	17.5
60° Swing In/Out w/Platform Flip Up	Short	11.5	14.5	14
	Long	14	18.5	16
70° Swing In/Out w/Platform Flip Up	Short	13	15.5	15
	Long	15.5	20	17.5
80° Swing In/Out w/Platform Flip Up	Short	13.5	14.5	15.5
	Long	15	19	17
70° Front Mount	No Ext	6	11.5	14
	2"	6	14	16
	4"	6	16	18
80° Front Mount	No Ext	6	12	14
	2"	6	14.5	16.5
	4"	6	16.5	18.5
90° Front Mount	No Ext	5.5	11.5	14
	Short	16.5	20	15
	Medium	17.5	20.5	16
ALR	Long	19	20.5	17
	Short	15.5	19.5	15
	Medium	15.5	20.5	16
ELR Ext. Mount	Long	17	20.5	17
	No Ext	5.5	9.5	13

6. REAR SEAT HEIGHT SET UP

Rear Seat Height

21.5" & 22" only available with amputee or straight axle plate & mountain wheel

☐ 184RH11	13.0"	NC	☐ 184RH26	17.5"	NC
☐ 184RH12	13.5"	NC	☐ 184RH20	18.0"	NC
☐ 184RH13	14.0"	NC	☐ 184RH27	18.5"	NC
☐ 184RH14	14.5"	NC	☐ 184RH21	19.0"	NC
☐ 184RH15	15.0"	NC	☐ 184RH37	19.5"	NC
☐ 184RH16	15.5"	NC	☐ 184RH22	20.0"	NC
☐ 184RH17	16.0"	NC	☐ 184RH38	20.5"	NC
☐ 184RH18	16.5"	NC	☐ 184RH23	21.0"	NC
☐ 184RH19	17.0"	NC	☐ 184RH39	21.5"	NC
			☐ 184RH40	22.0"	NC

Axle Plate

☐ 184AL33		Standard Adjustable	NC
☐ 184AL5		Curved Adjustable	\$85.00
☐ 184AL6	E0959	Amputee Adjustable NA w/ Padded Armrest	\$175.00
☐ 184AL29		Offset Adjustable Rear seat height will be 1/4" lower than standard RSTFH	NC

One Arm Drive

NA w/ Xtender, 12", 21" & 22" wide, HD package, curved & amputee axle plates, and small tray table. NA w/ combination of Padded SA arms and Anti tips. Only available w/ 20", 22", 24" mag wheels.

☐ 184AL16	E0958	One Arm Drive Adj. Axle Plate	\$850.00
☐ 184OA1		Small Handrim (NA w/ 20" wheels)	NC
☐ 184OA2		Standard Handrim	NC
184HM1	☐	Right Hand Mount	NC
184HM2	☐	Left Hand Mount	NC

Rear Wheel Type

Not applicable w/ Xtender Power Assist Option. Ultra Lightweight Wheels only available w/ Alum Anodized Handrim. Handrims cannot be omitted on Lite Spoke Wheels. Contact Built 4 Me for Spinergy color options.

☐ 184RW3	20" Mag Wheel	NC
☐ 184RW8	22" Mag Wheel	NC
☐ 184RW1	24" Mag Wheel	NC
☐ 184RW102	24" Mag Wheel 5 Spoke	NC NEW!
☐ 184RW90	20" Lite Spoke Wheel	\$200.00
☐ 184RW91	22" Lite Spoke Wheel	\$200.00
☐ 184RW92	24" Lite Spoke Wheel	NC
☐ 184RW93	25" Lite Spoke Wheel	NC
☐ 184RW94	26" Lite Spoke Wheel	\$80.00
☐ 184RW77	22" Ultra Lightweight Spoke Wheel	\$250.00
☐ 184RW78	24" Ultra Lightweight Spoke Wheel	\$250.00
☐ 184RW28	22" Spinergy SPOX Black	\$795.00
☐ 184RW95	24" Spinergy SPOX Black	\$795.00
☐ 184RW96	25" Spinergy SPOX Black	\$795.00
☐ 184RW97	26" Spinergy SPOX Black	\$795.00
☐ 184RW32	22" Spinergy LX Black	\$995.00
☐ 184RW98	24" Spinergy LX Black	\$995.00
☐ 184RW99	25" Spinergy LX Black	\$995.00
☐ 184RW40	26" Spinergy LX Black	\$995.00
☐ 184RW50	24" Mountain	\$250.00
☐ 184RW25	Omit Rear Wheel and Axle (Intended Size for Set Up _____)	-\$55.00
☐ 184RW49	Omit Rear Wheel (Use with Xtender)	NC

Rear Wheel Tire

☐ 184RT1	E2211 E2212	Pneumatic	NC
☐ 184RT2	E2211 E2213	Pneumatic Airless Insert	\$80.00
☐ 184RT3	E2218	Full Poly	NC
☐ 184RT6	E2211 E2212	High-Pressure Clincher	\$95.00
☐ 184RT27	E2211 E2212	Knobby V-TRAK	\$95.00
☐ 184RT21	E2211 E2212	Schwalbe Marathon Plus Evolution	\$125.00
☐ 184RT22	E2218	Kik® Mako Black	\$95.00
☐ 184RT12	E2211 E2212	Mountain Tire Only available on 24" Mountain Rear Wheel	\$95.00


Rear Seat Height Availability Matrices

Rear tire, side frame, and axle plate selections dictate seat height availability.
*Not all combinations are possible within acceptable seat height ranges

Rear Tire	Min RSH	Max RSH
20" Pneumatic	13"	19.0"
20" Poly	13"	19.0"
22" Pneumatic	14.0"	20.0"
22" Poly	14.0"	20.0"
22" HP Clincher	13.5"	19.5"
22" Schwalbe	13.5"	19.5"
24" Pneumatic	15"	21.0"
24" Poly	15"	21.0"
24" HP Clincher	14.5"	20.5"
24" Knobby	15"	21.0"
24" Schwalbe	14.5"	20.5"
24" Kik	14.5"	20.5"
24" Mountain	15.0"	21.0"
25" HP Clincher	15.0"	21.0"
25" Schwalbe	15.0"	21.0"
25" Kik	14.5"	20.5"
26" Pneumatic	16.0"	22.0"
26" HP Clincher	15.5"	21.5"
26" Schwalbe	15.5"	21.5"

Rear Side Frame	Min RSH	Max RSH
Curved	13"	21.0"
Straight	13.5"	22.0"
DAB	14.0"	22.0"

Axle Plate	Min RSH	Max RSH
Curved	13.5"	22.0"
Standard	14"	22.0"
Offset	13.5"	21.5"
Amputee	14"	22.0"



	Wheel	5 Spoke Mag	Mag	Lite Spoke	Ultra Lt wt	Spinergy
Tire						
Pneumatic		24	20,22,24	20,22,24,26		24,26
Pneumatic w/Airless Insert		24	20,22,24	20,22,24,26		24,26
Full Poly Tire		24	20,22,24	20,22,24		
High Pressure Clincher				22,24,25,26	22,24	22,24,25,26
Knobby V-TRAK				24		24
Schwalbe Marathon				24,25,26	22,24	24,25,26
Kik Mako Black		24	24	24,25		24,25

22" = 489 mm 24" = 540 mm

Customer Service: 800-333-4000 Fax: 800-300-7502 Email: orders@sunmed.com www.sunrisemedical.com

Specifications are subject to change without notice. Copyright © 2012 Sunrise Medical (US) LLC. (11.05) Quickie 2 MK-100059 REV B

6. REAR SEAT HEIGHT SET UP (CONT.)

Handrim

*Natural Fits available in 22", 24", and 25". Handrim selection NA w/ Xtender.
Ultra lightweight wheels only available w/ alum anodized handrim. Plastic coated limited to certain rims and sizes*

☐	184HR1	Aluminum Anodized	NC
☐	184HR4	Plastic Coated	\$95.00
☐	184HR22	Natural Fit - Standard Grip	\$375.00
☐	184HR23	Natural Fit - Super Grip	\$410.00
☐	184HR30	Natural Fit LT- Standard Grip	\$375.00
☐	184HR31	Natural Fit LT- Super Grip	\$410.00
☐	184HR24	Spinergy Flex Rim Available on 24" & 25" Spinergy LX wheel	\$525.00

Handrim (Cont)

☐	184HR7	Omit Handrim Option may not be available on all handrims	NC
☐	184HR8	8 Oblique NA w/ Xtender, plastic coated handrims. Only available with 24" or 26" Std Mag or Lite Spoke rims.	\$325.00

Axle

☐	184AX1	Quick Release Steel	NC
☐	184AX11	Quad Quick Release Std Mag and Mountain Wheel only	\$125.00

7. BACKREST

Back Height

☐	184B1	Low (9"-12")	NC
☐	184B5	Medium (12"-15")	NC
☐	184B9	Tall (16"-19")	NC
☐	184B66	Depth Adjustable 12"	NC
☐	184B67	Depth Adjustable 14"	NC
☐	184B68	Depth Adjustable 16"	NC
☐	184B69	Depth Adjustable 18"	NC
☐	184B70	Depth Adjustable 20"	NC

Back Height Setting (Std or Angle Adjustable)

☐	184UX4	9"	NC	☐	184UX18	16"	NC
☐	184UX6	10"	NC	☐	184UX20	17"	NC
☐	184UX8	11"	NC	☐	184UX22	18"	NC
☐	184UX10	12"	NC	☐	184UX24	19"	NC
☐	184UX12	13"	NC	☐	184UX26	20"	NC
☐	184UX14	14"	NC	☐	184UX28	21"	NC
☐	184UX16	15"	NC	20" & 21" only w/ Ang Adj Back			

Back Type

☐	184BT10	Standard	NC
☐	184BT20	Angle Adjustable Back (3-5")	\$185.00
☐	184BT22	Angle Adjustable Back (5"-7")	\$185.00
☐	184BT6	Depth Adjustable (Required with DAB)	NC

Back Post Options (Std or Angle Adjustable)

Omit push handles NA w/ stabilizer bar or transit

☐	184BP1	Push Handle Straight	NC
☐	184BP3	Push Handle 8° Med or Tall only. NA w/ Angle Adjustable Backrest	\$25.00
☐	184BP2	Omit Push Handle Straight (NA w/HD)	NC
☐	184BP4	Omit Push Handle 8° (NA w/HD)	\$25.00
☐	184BP8	Omit Back Post & Upholstery Use w/ DAB solid back & Hardware	-\$35.00
☐	184BP5	Angle Adj. Straight Push Handles	NC
☐	184BP6	Angle Adj. Omit Push Handles	NC

Back Post Options (Depth Adjustable Back)

Diagram of DAB in Appendix on page 9. SA Stroller Handle height range is 21.25"-28.5"

☐	184BP33	DAB - Swing Away Stroller Handles	\$100.00
☐	184BP27	DAB - Post Standard (19")	NC
☐	184BP28	DAB - Post Medium (22")	NC
☐	184BP29	DAB - Post Tall (25")	NC

Stabilizer Bar

☐	184SZ1	Stabilizer Bar Standard with Heavy Duty Package	\$105.00
--------------------------	--------	---	----------

Backrest Upholstery (Standard or Angle Adjustable)

☐	184BU1	Standard Upholstery	NC
☐	184BU70	E2611 Backrest Tension Adjust Ballistic NA on Heavy Duty NEW!	\$315.00
☐	184BU71	E2611 Backrest Tension Adjust 3DX Vented NA on Heavy Duty NEW!	\$330.00
☐	184BU82	E2611 Jay GO Backrest NA on Heavy Duty NEW!	\$400.00
☐	184BU28	Omit DAB Solid Back and Hardware	NC
☐	184BU29	Omit DAB Solid Back Include Hardware	\$80.00
☐	184BU23	Depth Adjustable Solid Back	\$195.00
☐	184BU14	Omit Back Upholstery	-\$15.00

Backrest Option

☐	184OD4	Reinforced Upholstery Only avail w/ Std Back Upholstery. Std w/HD package	\$50.00
--------------------------	--------	---	---------

JAY Back (Optional)

☐	184JB01	J3 Back (See J3 Order Form)	Price List
--------------------------	---------	-----------------------------	------------

Stroller Handles (Removable Push Handles)

Only available w/ omit push handles. NA w/Xtender, DAB, amputee axle plate, push handles when angle adjustable back is selected, flip-back or flip-back height adjustable armrests, and JAY backs.

☐	184SR2	Regular Stroller Handles	\$125.00
☐	184SR3	Tall Stroller Handles	\$125.00

Angle Adjustable Backrest Chart


Back Tube Height	Low	Medium	Tall
Pivot Regular set at 3"	12" 13" 14"	15" 16" 17"	17" 18" 19"
Pivot Regular set at 4"	13" 14" 15"	16" 17" 18"	18" 19" 20"
Pivot Regular set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot Tall set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot Tall set at 6"	15" 16" 17"	18" 19" 20"	N/A
Pivot Tall set at 7"	16" 17" 18"	19" 20" 21"	N/A
Adjustable Pivot Point Measurement from Seat Sling			
If back height option is not chosen, middle option will be present			

7. BACKREST (CONT.)

Angle Adjustable Backrest Pivot Height

Standard push handles only. Some interference w/ Padded SA arms. NA w/ flip-back or flip-back height adjustable armrests (classic).

☐	184AH1	Angle Adjustable Pivot 3"	NC
☐	184AH2	Angle Adjustable Pivot 4"	NC
☐	184AH3	Angle Adjustable Pivot 5"	NC
☐	184AH4	Angle Adjustable Pivot 6"	NC
☐	184AH5	Angle Adjustable Pivot 7"	NC



Preset Angle Adjustable Backrest (-15° to +3°) (Optional)

☐	184BD5	Back Angle - 15°	NC
☐	184BD4	Back Angle - 12°	NC
☐	184BD3	Back Angle - 6°	NC
☐	184BD2	Back Angle - 0°	NC
☐	184BD1	Back Angle +3°	NC

Military Service Patches (Optional)

NA with standard Upholstery

☐	184VU1	Marines Patch	\$90.00
☐	184VU2	Army Patch	\$90.00
☐	184VU3	Navy Patch	\$90.00
☐	184VU4	Air Force Patch	\$90.00
☐	184VU12	Coast Guard Patch	\$90.00
☐	184VU5	Purple Heart Patch	\$90.00
☐	184VU8	Purple Heart with Marines Patch	\$145.00
☐	184VU9	Purple Heart with Army Patch	\$145.00
☐	184VU10	Purple Heart with Navy Patch	\$145.00
☐	184VU11	Purple Heart with Air Force Patch	\$145.00
☐	184VU13	Purple Heart with Coast Guard Patch	\$145.00

8. ARMREST

Padded Swing Away

NA w/ DAB, Xtender or Amputee axle. Some limitations may exist with integral push handles, some RSTFHs and wheel sizes. **Now available w/Transit.

☐	184AA18	E0973	Padded Swing-Away STD (9" to 13")	\$210.00
		Setting	☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"	

Height Adjustable

☐	184AA11	E0973	Single Post Hgt. Adj. Std (8.5"-13.5")	\$285.00
☐	184AA29	E0973	Single Post Hgt. Adj. Low (7"-11")	\$285.00

Flip Back

NA w/ DAB, Xtender, small frame length, amputee axle, highest axle on rounded rear frame, or angle adjustable back.

☐	184AA6	E0973	Flip Back Height Adj (Classic 8.5"-13.5")	\$295.00
☐	184AA65	E0973	Flip Back Height Adj (LITE 8.5"-13.5")	\$295.00
☐	184AA5	K0015	Flip Back (incl. Desk Length Pad)	NC

Flip Up (For use with Depth Adjustable Back)

Requires an Armpad or Foam from the Arm Pad section below

☐	184AA7	E0973	Length Adjustable Locking Flip-up	\$265.00
☐	184AA1	E0973	Adj. Locking Flip-up (Std. w/ foam)	\$255.00

Omit Armrest

☐	184AA14		Omit Armrest	NC
--------------------------	---------	--	--------------	----

Armpad for Height Adjustable

For flip back, single post, and flip up- contour pad not recommended for tray tables. Contour only available with single post height adjustable.

☐	184AM23		Desk Length - Classic (10")	NC
☐	184AM21		Desk Length - Contour (10")	NC
☐	184AM24		Full Length - Classic (14")	NC
☐	184AM22		Full Length - Contour (14")	NC
☐	184AM9		Armpad Foam - Short (Flip Up armrests)	NC
☐	184AM10		Armpad Foam - Long (Flip Up armrests)	NC

9. KEY OPTIONS & ACCESSORIES

Wheel Lock

Ergo & compact NA w/ small frame or amputee axle plate. Additional No Gos may exist depending on wheel size and position

☐	184WL1	High Mount Push	NC
☐	184WL2	High Mount Pull	NC
☐	184WL27	Ergo Scissor	\$110.00
☐	184WL30	Compact	\$125.00
☐	184WL10	Omit	NC

Wheel Lock Option

☐	184WL24	Unilateral Wheel Lock Right	\$200.00
☐	184WL23	Unilateral Wheel Lock Left	\$200.00
☐	184WL12	Do Not Mount	NC

Extension Handle

NA w/ ergo scissor or compact wheel locks

☐	184EH7	E0961	Extension Handle (Pair)	\$65.00
--------------------------	--------	-------	-------------------------	---------

Grade Aid

Only available w/ high mount push or pull to lock wheel locks on pneumatic or pneumatic w/ airless tire options

☐	184GA1	E0974	Grade Aid	\$150.00
--------------------------	--------	-------	-----------	----------

Seat Sling

☐	184SL1	Standard Seat Sling	NC
☐	184SL6	Omit Seat Sling	-\$15.00

Cushion

☐	184CU1	E2601	2" Cushion	\$120.00
--------------------------	--------	-------	------------	----------

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

☐	184JC20	JAY Bundled Cushion		
		Cushion Name or Part		Price List
		Width: _____	Depth: _____	

Side Guard

Plastic side guards NA w/ height adj. arms. Side guard interference w/ Padded SA arms at some heights may exist.

☐	184SG33	K0108	Plastic Standard (8.5" High)	\$125.00
☐	184SG29	K0108	Plastic Low (6.5" High)	\$125.00

Anti-Tip

☐	184AT1	E0971	Rear Anti-Tip (NC w/ Xtender)	\$140.00
--------------------------	--------	-------	-------------------------------	----------

Transit Option

NA on Stroller Handles (Removable Push Handles) or Xtender

☐	184TO46		Transit System	\$315.00
--------------------------	---------	--	----------------	----------

Backpack

☐	184A1		Quickie Backpack Black	\$55.00
☐	184A2		Quickie Kids Backpack Black	\$50.00

Seat Pouch

☐	184A8		Quickie Seat Pouch Black	\$45.00
--------------------------	-------	--	--------------------------	---------

Additional Cover & Pack Accessories

☐	184A500		Gel Cover Hanger/Frame	\$105.00
☐	184A994		Gel Cover Swing Away Arm	\$105.00

Positioning Belts

☐	184A141	E0978	Belt 1.5" Auto Buckle w/Clamps	\$70.00
☐	184A117	E0978	Belt 2" Aircraft Buckle w/Clamps	\$110.00
☐	184A118	E0978	Belt 2" Aircraft Buckle Pad w/Clamp	\$120.00
☐	184A120	K0108	Positioning Belt Slide Clamps	\$50.00

Spoke Guards

Only available w/ Std Mag, Mountain, and Lite Spoke wheels

☐	184A302	K0065	Quickie Clear Spoke Guards	\$115.00
☐	184A23	K0065	Black Spoke Guards	\$115.00

Miscellaneous

☐	184A82		Caddy (NA w/ front mount hangers)	\$45.00
☐	184A49	E2207	Crutch Holder	\$195.00
☐	184A55	K0038	Extra Leg Strap Black	\$30.00
☐	184IG2		Impact Guard Neoprene	\$85.00
☐	184A73	K0105	IV Holder	\$275.00
☐	184A145	E2208	O2 Holder (NA w/ transit)	\$310.00
☐	184A83		Tool Kit	\$45.00

Tray Tables

☐	184A30	E0950	Small (13"-14")	\$250.00
☐	184A31	E0950	Medium (15"-17")	\$250.00
☐	184A32	E0950	Large (18"-20")	\$250.00
☐	184A33	E0950	X-Large (21"-22")	\$250.00

Color selection on next page

9. KEY OPTIONS & ACCESSORIES (CONT.)

Standard Frame Colors			
☐	184C140	Aztec Gold	NC
☐	184C2	Black <i>(Not avail as base w/Zebra)</i>	NC
☐	184C141	Black Cherry	NC
☐	184C10	Black Opal	NC
☐	184C38	Blue Opal	NC
☐	184C36	Candy Blue	NC
☐	184C14	Candy Purple	NC
☐	184C35	Candy Red	NC
☐	184C142	Electric Blue	NC
☐	184C143	Evergreen	NC
☐	184C126	Glow	NC
☐	184C144	Green Apple	NC
☐	184C147	Hot Sparkle Pink	NC
☐	184C148	Mauve Pink	NC
☐	184C103	Rootbeer	NC
☐	184C114	Sparkle Silver	NC
☐	184C145	Sunrise Orange	NC
☐	184C113	Titanium Color Paint	NC
☐	184C8	Yellow	NC

Matte Frame Colors			
☐	184C116	Matte Black	NC
☐	184C151	Matte Black Cherry	NC
☐	184C150	Matte Electric Blue	NC
☐	184C149	Matte Evergreen	NC
☐	184C152	Matte Purple	NC

Custom Frame Color			
☐	184C112	Desert Camouflage	\$295.00
☐	184C43	Mossy Oak Camouflage	\$295.00
☐	184C42	Stars & Stripes	\$295.00
☐	184C108	Zebra <i>(Choose above base color) Std with White unless above base paint is chosen NA w/ Black</i>	\$295.00

Additional Options Through Built for Me

To inquire about any unique customization for Quickie products, contact: Phone: (800) 333 - 4000, Fax: (800) 333 9011



Commonly Requested Options Include:

- ✓ Up to 24" wide
- ✓ Reclining Back
- ✓ Custom Anodizing
- ✓ Dual Mount Wheel Locks - Attendant Foot Lock
- ✓ Fold-Down Push Handle
- ✓ Dynamic Back
- ✓ Hub Wheel Locks on Mag or Spoke Wheels



10. APPENDIX

Quickie 2 Seat To Floor Height Matrix - Standard Frame *Hemi options will result in a reduction of approx 2 1/2"
 * Values Given are Estimates. Actual seat to floor height may vary based on seat depth and front seat to floor height selection.

All 3" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	16.5	17.0
	3/4"	-	-	-	17.0	17.5
	1 1/2"	-	-	-	18.0	18.5
Multi Position Tall	Std	-	-	-	18.5	19.0
	3/4"	-	-	-	19.0	19.5
	1 1/2"	-	-	-	20.0	20.5
3"	Std	-	-	-	-	15.0
	3/4"	-	-	-	-	15.5
	1 1/2"	-	-	-	-	16.5
5 1/4"	Std	-	-	-	-	17.0
	3/4"	-	-	-	-	18.0
	1 1/2"	-	-	-	-	18.5
6"	Std	-	-	-	-	18.0
	3/4"	-	-	-	-	18.5
	1 1/2"	-	-	-	-	19.5
7"	Std	-	-	-	-	19.0
	3/4"	-	-	-	-	19.5
	1 1/2"	-	-	-	-	20.5
8" X 2"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	16.5	17.0
	1 1/2"	-	-	-	18.0	18.5

All 4" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	16.5	17.5
	3/4"	-	-	-	17.0	18.0
	1 1/2"	-	-	-	18.5	19.0
Multi Position Tall	Std	-	-	-	18.5	19.5
	3/4"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.5	21.0
3"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	-	16.5	17.5
	3/4"	-	-	-	18.0	18.5
	1 1/2"	-	-	-	18.5	19.0
6"	Std	-	-	-	18.0	18.5
	3/4"	-	-	-	18.5	19.0
	1 1/2"	-	-	-	19.5	20.0
7"	Std	-	-	-	19.0	19.5
	3/4"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.5	21.0
8" X 2"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	16.5	17.5
	1 1/2"	-	-	-	18.0	19.0

All 5" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	16.5	17.0
	3/4"	-	-	-	17.5	18.0
	1 1/2"	-	-	-	18.0	18.5
Multi Position Tall	Std	-	-	-	18.5	19.0
	3/4"	-	-	-	19.0	19.5
	1 1/2"	-	-	-	20.0	20.5
3"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	-	17.0	17.5
	3/4"	-	-	-	18.5	19.0
	1 1/2"	-	-	-	19.0	19.5
6"	Std	-	-	-	18.0	18.5
	3/4"	-	-	-	19.0	19.5
	1 1/2"	-	-	-	20.0	20.5
7"	Std	-	-	-	19.0	19.5
	3/4"	-	-	-	20.0	20.5
	1 1/2"	-	-	-	21.0	21.5
8" X 2"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	17.5	18.0
	1 1/2"	-	-	-	19.0	19.5

All 6" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	17.0	17.5
	3/4"	-	-	-	18.5	19.0
	1 1/2"	-	-	-	19.0	19.5
Multi Position Tall	Std	18.0	18.5	19.0	19.5	20.0
	3/4"	-	-	-	20.5	21.0
	1 1/2"	-	-	-	21.0	21.5
3"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	-	17.5	18.0
	3/4"	-	-	-	18.5	19.0
	1 1/2"	-	-	-	19.5	20.0
6"	Std	-	-	-	17.5	18.0
	3/4"	-	-	-	18.5	19.0
	1 1/2"	-	-	-	20.0	20.5
7"	Std	-	-	-	18.5	19.0
	3/4"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.5	21.0
8" X 2"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	-	18.5
	1 1/2"	-	-	-	-	20.0

All 7" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	19.0
	3/4"	-	-	-	-	19.5
	1 1/2"	-	-	-	-	20.5
Multi Position Tall	Std	-	-	-	19.5	20.0
	3/4"	-	-	-	20.0	20.5
	1 1/2"	-	-	-	21.0	21.5
3"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	-	18.5	19.0
	3/4"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.0	20.5
6"	Std	-	-	-	19.0	19.5
	3/4"	-	-	-	20.0	20.5
	1 1/2"	-	-	-	21.0	21.5
7"	Std	-	-	-	19.5	20.0
	3/4"	-	-	-	20.0	20.5
	1 1/2"	-	-	-	21.0	21.5
8" X 2"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	-	-
	1 1/2"	-	-	-	-	-

All 8" x 1.5" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Multi Position Tall	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
3"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	-	18.5	19.0
	3/4"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.0	20.5
6"	Std	-	-	-	19.0	19.5
	3/4"	-	-	-	20.0	20.5
	1 1/2"	-	-	-	21.0	21.5
7"	Std	-	-	-	19.5	20.0
	3/4"	-	-	-	20.0	20.5
	1 1/2"	-	-	-	21.0	21.5
8" X 2"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	-	-
	1 1/2"	-	-	-	-	-

All 8" x 2" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Multi Position Tall	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
3"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
6"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
7"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
8" X 2"	Std	-	-	-	19.5	20.5
	3/4"	-	-	-	20.0	21.0
	1 1/2"	-	-	-	22.0	-
Frog Legs	Std	-	-	-	-	-
	1 1/2"	-	-	-	-	-

All Other 8" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Multi Position Tall	Std	19.0	19.5	20.0	20.5	21.0
	3/4"	-	-	-	21.5	22.0
	1 1/2"	-	-	-	22.0	-
3"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	-	-	19.5
	3/4"	-	-	-	-	20.0
	1 1/2"	-	-	-	-	21.0
6"	Std	-	-	-	19.5	20.0
	3/4"	-	-	-	20.5	21.0
	1 1/2"	-	-	-	21.0	21.5
7"	Std	-	-	-	19.5	20.0
	3/4"	-	-	-	20.5	21.0
	1 1/2"	-	-	-	21.5	22.0
8" X 2"	Std	-	-	-	-	-
	3/4"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	-	-
	1 1/2"	-	-	-	-	-

Depth Adjustable Back Frame

