



Quickie GT December 2014

EIR3

Mark For:

Submitting for:

☐ Quote

☐ Order

Date: _____

PO#: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.

The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME STYLE & SETUP

Base Model

☐	145WM13	K0005	GT (265 lbs)	\$2,325.00
--------------------------	---------	-------	--------------	------------

Frame Angle (see Figure 1)

☐	145FA1		75°	NC
☐	145FA8		80°	NC
☐	145FA2		85°	NC

Frame Inset Per Side (See Figure 2)

☐	145FF1		0"	NC
☐	145FF2		1" (NA w/ 12" wide)	NC
☐	145FF3		2" (NA w/ 12", 13", or 14" wide)	NC

Frame Inset

Footrest Width

0	Subtract 0 from Frame Width
1"	Subtract 2" from Frame Width
2"	Subtract 4" from Frame Width


Frame Style (See Figure 3)

☐	145FS34		Classic Frame	NC
☐	145FS35		Performance Frame	NC
☐	145FS36		XTR Suspension Frame Only available in Performance Frame Set Up	\$895.00
☐	145SK1		Light Shock (up to 180 lbs.)	NC
☐	145SK3		Heavy Shock (up to 265 lbs.)	NC


Depth Adjustable Backrest (Optional)

See **Figure 5** Depth Adjust Backrest extends the frame rear of the seat sling, select a factory setting. NA on 19" and 20" Seat Depths.


☐	145DP02	0 Setting Backrest at standard position, field adjustable to -1" or -2"	\$165.00
☐	145DP03	-1 Setting Backrest -1" offset, field adjustable to flush or -2"	\$165.00
☐	145DP04	-2 Setting Backrest -2" offset, field adjustable to flush or -1"	\$165.00



0 POSITION



-1" POSITION



-2" POSITION

Figure 1




Figure 2





Figure 3



Classic Frame - Traditional GT seat sling to front frame bend angle dimensions.

Performance Frame - Seat Sling to front frame bend angle is 1.5" shorter than the classic position.

Front of back cane post to front edge of seat sling will remain the same for both Classic and Performance.

1. FRAME STYLE & SETUP (CONT)

Seat Width

See Figure 2. Measure from outside of seat tube to outside of opposite seat tube. XTR and J3 Back NA w/ 12" and 13" widths

☐	145S3	12"	NC
☐	145S4	13"	NC
☐	145S5	14"	NC
☐	145S6	15"	NC
☐	145S7	16"	NC
☐	145S8	17"	NC
☐	145S9	18"	NC
☐	145S10	19"	NC
☐	145S11	E2201 20"	\$500.00

Seat Depth

See Figure 1. Measures from front of back post to front edge of seat sling. Seat sling will be 1" shorter than depth ordered, but actual seat depth will remain the same as ordered.

☐	145SD5	12"	NC
☐	145SD6	13"	NC
☐	145SD7	14"	NC
☐	145SD8	15"	NC
☐	145SD9	16"	NC
☐	145SD10	17"	NC
☐	145SD11	18"	NC
☐	145SD12	19"	NC
☐	145SD13	E2203 20"	\$550.00

2. FRONT SEAT HEIGHT SET UP

Front Seat Height

See Figures 4 and 5

☐	145SH14	16"	NC	☐	145SH20	18.5"	NC
☐	145SH15	16.5"	NC	☐	145SH21	19"	NC
☐	145SH16	17"	NC	☐	145SH23	19.5"	NC
☐	145SH17	17.5"	NC	☐	145SH24	20"	NC
☐	145SH18	18"	NC	☐	145SH26	20.5"	NC
☐				☐	145SH27	21"	NC

Front to Rear Seat Height Differences

Frame Depth	Max Seat Squeeze
12"/13"	2.5"
14"/15"	3.0"
16"/17"	3.5"
18"/19"	4.0"
20"	4.5"

Caster Wheels

☐	145CW1	3" Micro	NC
☐	145CW3	4" Low Profile Polyurethane	NC
☐	145CW31	4" Micro Lighted	\$70.00
☐	145CW72	4" x 1" Aluminum Polyurethane	\$100.00
☐	145CW44	4" x 1.25" Semi Pneumatic	\$55.00
☐	145CW64	4" x 1.5" Aluminum Soft Roll	\$175.00
☐	145CW5	5" Low Profile Polyurethane	NC
☐	145CW51	5" Micro Lighted	\$85.00
☐	145CW6	5" Soft Roll (Grey) NA w/ Frog Legs	\$85.00
☐	145CW65	5"x 1" Aluminum Polyurethane	\$110.00
☐	145CW7	5"x 1.5" Semi Pneumatic	\$85.00
☐	145CW66	5"x 1.5" Aluminum Soft Roll	\$185.00
☐	145CW34	6" Polyurethane	NC
☐	145CW35	6"x 1.5" Semi Pneumatic	\$85.00
☐	145CW68	6"x 1.5" Aluminum Soft Roll	\$195.00

Caster Fork

See Figures 5 and 6

☐	145CF27	Short Fork	NC
☐	145CF45	Short Wide Fork (w/ 1.25" or 1.5" caster)	NC
☐	145CF32	Tall Wide Fork	NC

Caster Fork Suspension

☐	145CF53	E1015 Frog Legs 0 weight to 100 lbs	\$425.00
☐	145CF54	E1015 Frog Legs 101 lbs weight to 190 lbs	\$425.00
☐	145CF55	E1015 Frog Legs 191 lbs weight to 265 lbs	\$425.00

Figure 4




Figure 5

Front		Fork			
STFH		Short/Single	Short Wide	Tall Wide	Frog Legs
Caster	3"	16"-20"	16"-20.5"	16.5"-21"	17.5"-21"
	4"	16"-20.5"	16.5"-21"	16.5"-21"	17.5"-21"
	5"	17",18",19",20",21"	17"-21"	17"-21"	18.5"-21"
	6"	NA	NA	18"-21"	19",20",21"

Figure 6

Caster Size		Fork	Tubular/Flat			Flared			Platform		
			Inset			Inset			Inset		
			0"	1"	2"	0"	1"	2"	0"	1"	2"
75° Frame Angle	3"	Short, Short W, Tall W				NA					
		Single Sided				NA					
		Froglegs				NA					
		Short, Short W, Tall W				NA					
		Single Sided				NA					
		Froglegs				NA	2"				
	4"	Short, Short W, Tall W	NA			NA	NA				
		Single Sided	NA			NA	NA				
		Froglegs	NA			NA	NA		.5"		
		Tall Wide	NA			NA	NA	2"			
		Froglegs	NA	NA		NA	NA	1.5"	NA	.5"	
	5"	Short, Short W, Tall W				NA					
		Single Sided				NA					
		Froglegs				NA					
	6"	Short, Short W, Tall W				NA					
		Single Sided				NA					
		Froglegs				NA					
80° Frame Angle	3"	Short, Short W, Tall W				NA					
		Single Sided				NA					
		Froglegs				NA					
	4"	Short, Short W, Tall W				NA					
		Single Sided				NA					
		Froglegs				NA					
	5"	Short, Short W, Tall W				NA					
		Single Sided				NA					
		Froglegs				NA					
	6"	Short, Short W, Tall W				NA					
		Single Sided				NA					
		Froglegs				NA					
85° Frame Angle	3"	Short, Short W, Tall W				NA					
		Single Sided				NA					
		Froglegs				NA					
	4"	Short, Short W, Tall W				NA					
		Single Sided				NA					
		Froglegs				NA					
	5"	Short, Short W, Tall W				NA					
		Single Sided				NA					
		Froglegs				NA					
	6"	Short, Short W, Tall W				NA					
		Single Sided				NA					
		Froglegs				NA					

Go

No Go

= "Soft" Go with limited seat angle
Min seat angle listed in cell

3. FOOTPLATE OPTIONS

Footrest			
<i>Heel loop or leg strap not included, See Footplate Options Section</i>			
☐	145F228	Aluminum Tubular	
☐	145F229	Aluminum Tubular Long	NC
☐	145F230	Aluminum Tubular w/ Flat Cover	NC
☐	145F231	Aluminum Tubular Long w/ Flat Cover	nc
☐	145F234	K0040 Aluminum Platform (Ang. Adj.)	\$95.00
☐	145F235	K0040 Aluminum Platform Long (Ang. Adj.)	\$95.00
☐	145F278	K0040 High Mount Platform NA w/ 12", 13" W with 1" or 2" Inset or 14" or 15" W with 2" insets.	\$95.00

Adaptive Platform Spacer			
<i>NA with Tubular Footrest</i>			
☐	145AD1	Adaptive Platform 2" Spacer	\$95.00
☐	145AD2	Adaptive Platform 4" Spacer	\$95.00
☐	145AD3	Adaptive Platform 6" Spacer	\$95.00

Footrest Options			
☐	145HL5	K0038 Leg Strap	\$35.00
☐	145DI01	Removable Snap On Cover NA w/ Platform Footrests	\$60.00
☐	145DI02	Wear Pads NA w/ Platform Footrests	\$50.00
☐	145A179	K0108 Positioning Plates for Platform Footrest	\$50.00
☐	145A54	K0038 Extra Leg Strap	\$30.00

4. REAR SEAT HEIGHT SET UP

Rear Seat Height			
<i>14" & 14.5" only available with 20" rear wheel. NA with Amputee Axle Plate.</i>			
<i>13" & 13.5" only available with 20" or 22" rear wheel. NA with Amputee Axle Plate.</i>			
☐	145RH11	13"	NC
☐	145RH12	13.5"	NC
☐	145RH13	14"	NC
☐	145RH14	14.5"	NC
☐	145RH15	15"	NC
☐	145RH16	15.5"	NC
☐	145RH17	16"	NC
☐	145RH18	16.5"	NC
☐	145RH19	17"	NC
☐	145RH26	17.5"	NC
☐	145RH20	18"	NC
☐	145RH27	18.5"	NC
☐	145RH21	19"	NC
☐	145RH37	19.5"	NC
☐	145RH22	20"	NC

Rear Wheel Tire			
<i>See Availability Matrix</i>			
☐	145RT1	E2211 E2212 Pneumatic	NC
☐	145RT2	E2211 E2213 Pneumatic w/ Airless Insert	\$80.00
☐	145RT3	E2218 Full Poly	NC
☐	145RT6	E2211 E2212 High-Pressure Clincher	\$95.00
☐	145RT27	E2211 E2212 Knobby V- Trak	\$95.00
☐	145RT21	E2211 E2212 Schwalbe Marathon Plus Evolution	\$125.00
☐	145RT38	E2218 SHOX Performance Poly	\$95.00
☐	145RT12	E2211 E2212 Mountain Tire	\$95.00

Rear Wheel Size			
<i>Ultra Lightweight Wheels only available with aluminum anodized handrim</i>			
<i>Contact Built 4 Me for Spinergy color options</i>			
☐	145RW3	20" Mag Wheel	NC
☐	145RW8	22" Mag Wheel	NC
☐	145RW1	24" Mag Wheel	NC
☐	145RW90	20" Lite Spoke Wheel	\$200.00
☐	145RW91	22" Lite Spoke Wheel	\$200.00
☐	145RW92	24" Lite Spoke Wheel	NC
☐	145RW93	25" Lite Spoke Wheel	NC
☐	145RW94	26" Lite Spoke Wheel	\$80.00
☐	145RW77	22" Ultra Lightweight Spoke Wheel	\$250.00
☐	145RW78	24" Ultra Lightweight Spoke Wheel	\$250.00
☐	145RW28	22" Spinergy SPOX Black	\$795.00
☐	145RW42	24" Spinergy SPOX Black	\$795.00
☐	145RW96	25" Spinergy SPOX Black	\$795.00
☐	145RW44	26" Spinergy SPOX Black	\$795.00
☐	145RW32	22" Spinergy LX Black	\$995.00
☐	145RW88	24" Spinergy LX Black	\$995.00
☐	145RW99	25" Spinergy LX Black	\$995.00
☐	145RW40	26" Spinergy LX Black	\$995.00
☐	145RW50	24" Mountain	\$250.00

Rear Tire & Wheel Matrix				
Tire	Wheel	Mag	Lite Spoke	Ultra Lightweight
Pneumatic		20,22,24	20,22,24,26	
Pneumatic w/ Airless Insert		20,22,24	20,22,24	
Full Poly		20,22,24	20,22,24	
High-Pressure Clincher			22,24,25,26	22,24
Knobby V Trak by Kenda			24	24
Schwalbe Marathon			24,25,26	
SHOX Performance Poly		24	24,25	
22" = 489 mm 24" = 540 mm 25" = 559 26" = 590				

Handrim			
<i>Aluminum and Plastic are not available on all wheels in all sizes. Natural Fits are NA w/ 26".</i>			
<i>Handrim Selection not applicable with Xtender Power Assist</i>			
☐	145HR1	Aluminum Anodized	NC
☐	145HR4	Plastic Coated NA w/ 25" & 26" Lite Spoke Wheels	\$95.00
☐	145HR22	Natural Fit - Standard Grip	\$375.00
☐	145HR23	Natural Fit - Super Grip	\$410.00
☐	145HR30	Natural Fit LT - Standard Grip	\$375.00
☐	145HR31	Natural Fit LT - Super Grip	\$410.00
☐	145HR24	Spinergy Flex Rim Only available on Spinergy LX wheels 24" & 25"	\$525.00
☐	145HR8	8 Oblique Projections NA w/ Kik tires, only available with Mag or Spoke	\$325.00
☐	145HR7	Omit Handrim NA w/ Ultra Lightweight Wheels	NC

4. REAR SEAT HEIGHT SET UP (CONT)

Xtender Power Assist

Rear seat height is 18.5" to 20" & 0 camber, NA w/ Amputee Axle, Padded Swing Arms.
Front To Rear Seat Height Range 0 to 1" difference, 24" wheels only. Not available with XTR or standard rear suspension. Anti-Tips no charge with Xtender

☐	145XT2	E0986	Xtender Power Assist 4X Multiplier NICKEL	\$6,595.00
--------------------------	--------	-------	---	------------

☐	145XT3	E0986	Xtender Power Assist 4X Multiplier LITHIUM	\$7,595.00
--------------------------	--------	-------	--	------------

Please select one of the following handrim options

☐	145XH3		Plastic Coated	NC
☐	145XH4		Stainless Steel	NC

Xtender - Other Options

☐	145XA1		Spare NiMH Power Pack - NICKEL	\$1,395.00
☐	145XA7		Spare NiMH Power Pack - LITHIUM	\$1,795.00

Axle

☐	145AX1		Quick Release Steel	NC
☐	145AX13		Quad Quick Release Only available with Mag wheels	\$125.00

Axle Plates

Selection not required for XTR Suspension

☐	145AL40		Standard Short (COG 1" to 4.5") Some limitations with STFH	NC
☐	145AL43		Standard Long (COG 1" to 4.5")	NC
☐	145AL37	E0959	Amputee (COG -1" to 2.5")	\$175.00
☐	145AL38	E1015	Suspension (COG-.5 to +3) NA with XTR Suspension Frame	\$395.00
☐	145AL39		Xtender (COG 1")	NC

		Short Axle Plate- Wheel Size:				
		20	22	24	24 (559)	26
Rear Seat to Floor Height	13	14	15	15.5	16	
	13.5	14.5	15.5	16	16.5	
	15	16	17	17.5	18	
	15.5	16.5	17.5	18	18.5	
	16	17	18	18.5	19	

Rear Suspension Rating

Selection not required for XTR Suspension

☐	145SV04		From 0 weight to 100 lbs (Yellow)	NC
☐	145SV05		From 101 lbs weight to 190 lbs (Red)	NC
☐	145SV06		From 191 lbs weight to 265 lbs (Blue)	NC

Camber

0,3,6 available with std frame. 0,4,8 available with XTR frame


☐	145CI1		0°	NC
☐	145CI2		3°	NC
☐	145CI3		4°	NC
☐	145CI4		6°	NC
☐	145CI5		8°	NC

Center of Gravity

Measure from front of back post to center of rear axle. XTR Suspension Frame will be factory set at 1.5", Xtender 1"

☐	145CC11		-1"	NC
☐	145CC10		-.5"	NC
☐	145CC12		0"	NC
☐	145CC1		.5"	NC
☐	145CC2		1"	NC
☐	145CC3		1.5"	NC
☐	145CC4		2"	NC
☐	145CC5		2.5"	NC
☐	145CC6		3"	NC
☐	145CC7		3.5"	NC
☐	145CC8		4"	NC
☐	145CC9		4.5"	NC

Figure 7



5. BACKREST

Back Type

☐	145BT3		Non-Folding Backrest	NC
☐	145BT82		Folding Lock Down Angle Adj	NC

Height Adjustable Backrest

☐	145B110		10" - 14" Back Height	NC
☐	145B127		12" - 16" Back Height	NC
☐	145B111		14" - 18" Back Height	NC
☐	145B137		16" - 20" Back Height	NC

Back Upholstery

☐	145BU1		Standard Upholstery	NC
☐	145BU8		Adjustable	\$175.00
☐	145BU57		Vented	\$200.00
☐	145BU14		Omit Upholstery (Use w/ JAY Backs)	-\$15.00

JAY J3 Back

☐	145JB01		JAY J3 Back	See OF
--------------------------	---------	--	-------------	--------

Military Service Patches

Only available on tension adjustable upholstery and Jay 3 Back

☐	145VU1		US Marines Patch	\$90.00
☐	145VU2		US Army Patch	\$90.00
☐	145VU3		US Navy Patch	\$90.00
☐	145VU4		US Air Force Patch	\$90.00
☐	145VU12		US Coast Guard	\$90.00
☐	145VU5		Purple Heart	\$90.00
☐	145VU8		Purple Heart with Marines Patch	\$145.00
☐	145VU9		Purple Heart with Army Patch	\$145.00
☐	145VU10		Purple Heart with Navy Patch	\$145.00
☐	145VU11		Purple Heart with Air Force Patch	\$145.00
☐	145VU13		Purple Heart with Coast Guard	\$145.00

Push Handles

☐	145PH6		Integral Push Handles	NC
☐	145PH5		Fold-Down Push Handles NA w/ 10-14" Back Height or Standard Upholstery.	\$100.00
☐	145PH2		Omit Push Handles	NC

5. ARMREST

Armrest - Padded Swing Away

Includes angle adjustable receiver hardware. NA w/ Xtender. Some limitations may exist w/ integral push handles, some RSHs and wheel sizes

☐	145AA18	E0973	Padded Swing-Away - Standard	\$210.00
--------------------------	---------	-------	------------------------------	----------

Armrest - Height Adjustable

☐	145AA11	E0973	Single Post Ht. Adjust. Std (8.5"-13.5")	\$285.00
☐	145AA29	E0973	Single Post Ht. Adjust. Low (7"-11")	\$285.00

Armpad - For Height Adjustable

Suggest using classic model for tray tables

☐	145AM23		Desk Length (10") Classic	NC
☐	145AM21		Desk Length (10") Contour	NC
☐	145AM24		Full Length (14") Classic	NC
☐	145AM22		Full Length (14") Contour	NC

6. Key Options & Accessories

Wheel Lock

Not all wheel locks available in all seat depths, wheel sizes and COG settings. XTR suspension only available w/ high mount push

☐	145WL1		High Mount Push	NC
☐	145WL2		High Mount Pull	NC
☐	145WL27		Ergo Scissor	NC
☐	145WL30		Compact	\$125.00
☐	145WL10		Omit Wheel Lock	NC
☐	145WL12		Do Not Mount	NC

Extension Handle

Only available w/ High Mount Push or Pull Wheel Locks

☐	145EH7	E0961	Extension handles (Pair)	\$65.00
--------------------------	--------	-------	--------------------------	---------

Grade Aid

Only available w/ high mount push or pull. NA w/ clincher, poly, Kik, or Schwalbe, or Xtender

☐	145GA1	E0974	Grade Aid	\$150.00
--------------------------	--------	-------	-----------	----------

Anti Tips

☐	145AT1	E0971	Anti Tips (Pair)	\$140.00
☐	145AT3		Omit Rear Anti Tip	NC

Seat Sling

☐	145SL1		Standard Seat Sling	NC
☐	145SL6		Omit Seat Sling	-\$15.00

Cushion

☐	145CU1	K0669	2" Cushion	\$120.00
--------------------------	--------	-------	------------	----------

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

☐	145JC20	JAY Bundled Cushion		
		Cushion Name or Part #		Price List
		Width:	Depth:	

Side Guards

All side guards NA w/ Height Adj. Arms. Some interference may exist w/ padded SA arms at lowest settings

☐	145SG33	K0108	Plastic - Standard	\$125.00
☐	145SG29	K0108	Plastic - Low	\$125.00
☐	145SG12	K0108	Carbon Fiber NA w/ push or pull wheel locks	\$450.00

Pack Accessories

☐	145A1		Quickie Backpack Black	\$55.00
☐	145A2		Quickie Kids Backpack Black	\$50.00
☐	145A8		Seat Pouch Black	\$45.00
☐	145A500		Gel Cover Hanger/Frame	\$105.00
☐	145A994		Gel Cover Swing Away Arm	\$105.00

Positioning Belt

☐	145A142	E0978	Belt 1.5" Auto Buckle Universal w/clamps	\$70.00
☐	145A12	E0978	Belt 2" Aircraft Buckle	\$110.00
☐	145A13	E0978	Belt 2" Aircraft Buckle Padded	\$120.00

Spoke Guards

☐	145A23	K0065	Quickie Black Spoke Guard	\$115.00
☐	145A302	K0065	Quickie Clear Spoke Guard	\$115.00

Miscellaneous

☐	145A32	E0950	Wheelchair Tray	\$250.00
☐	145IG6		Neoprene Impact Guards	\$85.00
☐	145A51	E2207	Crutch Holder Only available with Integral or Fold Down push handles	\$195.00
☐	145A85		Bolt-On Push Handles	\$115.00
☐	145A96		Presta Valve Adapter	\$5.00
☐	145A83		Tool Kit	\$45.00

6. Key Options & Accessories (Cont)

Standard Frame Colors

☐	145C140	Aztec Gold	NC
☐	145C2	Black <i>(Not avail as base w/Zebra)</i>	NC
☐	145C141	Black Cherry	NC
☐	145C10	Black Opal	NC
☐	145C38	Blue Opal	NC
☐	145C36	Candy Blue	NC
☐	145C14	Candy Purple	NC
☐	145C35	Candy Red	NC
☐	145C142	Electric Blue	NC
☐	145C143	Evergreen	NC
☐	145C126	Glow	NC
☐	145C144	Green Apple	NC
☐	145C147	Hot Sparkle Pink	NC
☐	145C148	Mauve Pink	NC
☐	145C103	Rootbeer	NC
☐	145C114	Sparkle Silver	NC
☐	145C145	Sunrise Orange	NC
☐	145C113	Titanium Color Paint	NC
☐	145C8	Yellow	NC

Matte Frame Colors

☐	145C116	Matte Black	NC
☐	145C151	Matte Black Cherry	NC
☐	145C150	Matte Electric Blue	NC
☐	145C149	Matte Evergreen	NC
☐	145C152	Matte Purple	NC

Custom Frame Color

☐	145C112	Desert Camouflage	\$295.00
☐	145C43	Mossy Oak Camouflage	\$295.00
☐	145C42	Stars & Stripes	\$295.00
☐	145C108	Zebra <i>(Choose above base color) Std with White unless above base paint is chosen NA w/ Black</i>	\$295.00

Anodized Color Package

Camber tubes are titanium grey except for the black option. Frog Legs black only. Black only w/ amputee, Xtender, suspension axle plates.

☐	145EB02	Black	NC
☐	145EB06	Blue	\$135.00
☐	145EB05	Green	\$135.00
☐	145EB04	Purple	\$135.00
☐	145EB03	Red	\$135.00
☐	145EB01	Ti Grey	\$135.00

