

The Quickie 2 is already the most popular ultra-light folding frame wheelchair available today, and we've found a way to make it even better. With twelve new features and options, the Quickie 2 now offers more fine tune adjustments for improved versatility, customization and ease of use. With the Quickie 2 Family the possibilities are endless.



Quickie 2 Highlights:

- Fully adjustable axle plate with camber, center of gravity, wheel base width and seat height adjustments
- Even more caster adjustability to accommodate any modifications of the frame that affect alignment
- Angle-adjustable backrest
- Fully modular design for complete customization
- The Quickie True Fit Program allows the chair to grow
- High quality Quickie components and hardware ensure durability and performance
- Widest variety of hanger, footplate, frame, caster, wheel, backrest and armrest options
- New Jay Precision™ Back option
- Travel wheel option
- One-arm drive option (see Quickie Kids picture)
- Uni-lateral (one-hand) wheel lock system option



FOLDING FAMILY

Quickie Folding Chairs – A complete folding line to suit your mobility needs.

Ultra-Light



Lightweight



Positioning



Pediatrics



Quickie 2 & Quickie 2 Kids & Quickie 2 HD

Front Seat Height (For Standard Frames Only)

CASTER FORK MOUNTING HOLE	5" LP POLY 5 1/4" FORK	5" LP POLY 7" FORK	6" PNEU 6" FORK	6" POLY 6" FORK	8" PN 5 1/4" FORK
Top	16 15/16"	18 3/4"	18 5/16"	18 3/16"	18 1/2"
Bottom	17 15/16"	19 3/4"	19 5/16"	19 3/16"	19 1/2"

Rear Seat Height (For Standard Frames and Standard Axle Plate Position Only)

REAR AXLE/HOLE POSITION	PNEU	FP	PNEU	LP	FP	HP	PNEU/ KEVLAR	HP CLNCHR
1 AND 2	18 7/8"	19 1/2"	19 3/4"	20"	20 1/4"	20 3/8"	20 1/2"	20 1/4"
2 AND 3	17 3/4"	18 3/8"	18 5/8"	18 7/8"	19 1/8"	19 1/4"	19 3/8"	20 1/8"
3 AND 4	16 5/8"	17 1/4"	17 1/2"	17 3/4"	18"	18 1/8"	18 1/4"	19"
4 AND 5	15 1/2"	16 1/8"	16 3/8"	16 5/8"	16 7/8"	17"	17 1/8"	17 7/8"

Overall Width (With Standard Camber and Padded Swing-Away Arm)

FRAME WIDTH	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
SPOKE WHEEL W/HANDRIMS	20 1/2"	21 1/2"	22 1/2"	23 1/2"	24 1/2"	25 1/2"	26 1/2"	27 1/2"	28 1/2"	29 1/2"	30 1/2"

Quickie 2	Quickie 2 HD	Quickie 2 Kids
Weight: 27 lbs		
Weight Capacity: 250 lbs	Weight Capacity: 350 lbs	Weight Capacity: 250 lbs

Quickie 2HP

Front Seat Height (For Standard Frames Only)

CASTER FORK MOUNTING HOLE	3"MICROS 5 1/4" FORK	3"MICROS 6" FORK	3"MICROS 7" FORK	5"LP POLY 5 1/4" FORK	5"LP POLY 6" FORK	5"LP PC 7" FORK
TOP	16 5/16"	17 1/8"	18 1/8"	17 3/16"	18"	19"
BOTTOM	17 5/16"	18 1/8"	19 1/8"	18 3/16"	19"	20"

Rear Seat Height (For Standard Frames and Standard Axle Plate Position Only)

REAR AXLE/HOLE POSITION	PNEU	LATEX	HP	PNEU	TEN/B-BALL	LATEX	PNEU
1 AND 2	19 3/4"	20 1/4"	20 3/8"	20 1/2"	20 5/8"	21 1/4"	21 5/8"
2 AND 3	18 5/8"	19 1/8"	19 1/4"	19 3/8"	19 1/2"	20 1/8"	20 1/2"
3 AND 4	17 1/2"	18"	18 1/8"	18 1/4"	18 3/8"	19"	19 3/8"
4 AND 5	16 3/8"	16 7/8"	17"	17 1/8"	17 1/4"	17 1/8"	18 1/4"

Overall Width (With Standard Camber and Padded Swing-Away Arm)

FRAME WIDTH	12"	13"	14"	15"	16"	17"	18"	19"	20"
SPOKE/MAG WHEEL W/HANDRIMS	21 1/2"	22 1/2"	23 1/2"	24 1/2"	25 1/2"	26 1/2"	27 1/2"	28 1/2"	29 1/2"

Quickie 2 HP
Weight: 22 lbs
Weight Capacity: 250 lbs

Quickie M6

FRAME WIDTH:	24"	26"	28"	30"
OVERALL WIDTH:	32.5"	34.5"	36.5"	38.5"
FOLDED WIDTH:	14"			

Quickie M6

Weight:	53 lbs. w/legrests (26 x 20)
Weight Capacity:	650 lbs

Quickie Recliner

Front Seat Height

CASTER FORK MOUNTING HOLE	5" LP POLY 6" FORK	6" PNEU 6" FORK	6" POLY 6" FORK	8" PNEU 5 1/4" FORK	8" POLY 5 1/4" FORK
TOP	17 3/4"	18 5/16"	18 3/16"	N/A	N/A
BOTTOM	18"	19 5/16"	19 3/16"	19 1/2"	19 5/16"

Rear Seat Height

REAR AXLE HOLE POSITION	22" FP Poly	24" FP Poly	24" HP Clincher	24" Pneu/ Kevlar
Top	15"	15 3/4"	15 7/8"	16"
Middle	16"	16 3/4"	16 7/8"	17"
Bottom	17"	17 3/4"	17 7/8"	18"

Armrest Height

Height-Adjustable

9" - 14"

Footrest Height

SWING-AWAY

ELR 14" - 21" 16 1/2" - 19"

Overall Width (Without Armrests)

Frame Width	14"	15"	16"	17"	18"	19"	20"	Folded
SPOKE/MAG WHEEL W/HANDRIMS & W/HANDRIMS & QUICK-RELEASE AXLES	22 1/4"	23 1/4"	24 1/4"	25 1/4"	26 1/4"	27 1/4"	28 1/4"	12 1/2"
SPOKE/MAG WHEEL W/HANDRIMS & THREADED AXLES	22"	23"	24"	25"	26"	27"	28"	12 1/2"

Quickie Recliner

Weight: 42 lbs. (w/ swing-away footrests, height-adjustable armrests, headrest extension and mag. wheels)

Weight Capacity: 250 lbs

Quickie LX & LXI

Front Seat Height

CASTER FORK MOUNTING HOLE	5"LP Poly 4 3/4" Fork	5"LP Poly 5 1/4" Fork	6"Poly 6" Fork	6"Pneu 6" Fork	8"Poly 5 1/4" Fork	8"Pneu 5 1/4" Fork	8"x2" Caster 8"x2" Fork
TOP	N/A	N/A	17"	17"	N/A	N/A	N/A
BOTTOM	16"	16 1/2"	18"	18"	18"	18"	N/A

Rear Seat Height

REAR AXLE HOLE POSITION (LX ONLY)	20" LP POLY	20" PNEU	22" LP POLY	22" PNEU	24" LP POLY	24" PNEU
TOP	15"	15 1/2"	16"	16 1/2"	17"	17 1/2"
BOTTOM	16"	16 1/2"	17"	17 1/2"	18"	18 1/2"

REAR AXLE HOLE POSITION (ADJ.AXLE)	20" LP POLY	20" PNEU	22" LP POLY	22" PNEU	24" LP POLY	24" PNEU
TOP	N/A	15"	15 1/2"	16"	16"	17"
MIDDLE	15 1/2"	16"	16 1/2"	17 1/2"	17 1/2"	18"
BOTTOM	17"	17 1/2"	18"	18 1/2"	18 1/2"	19 1/2"

Armrest Height

Footrest Height

FLIP-BACK HEIGHT ADJ.	FLIP-BACK/ PADDED S/A HEIGHT ADJ.	HT. ADJ.	SWING-AWAY	ELR/ALR
9 1/4"	8 3/4"-12 5/8"	8 3/4"	7 1/4"-12 1/4"	12 3/4"-15 1/4"
			8 3/4"-13 3/4"	14"-16 1/2"
				15 1/4"-17 3/4"
				16 1/2"-19"
				17 3/4"-20 1/4"

Overall Width (Without Armrests)

FRAME WIDTH	12"	14"	16"	18"	20"	FOLDED
MAG WHEEL W/HANDRIMS	19"	21"	23"	25"	27"	12 1/2"
SPOKEWHEEL W/HANDRIMS	20"	22"	24"	26"	28"	12 1/2"

Quickie Recliner

Weight: 28 lbs
Weight Capacity: 250 lbs

SPECIFICATIONS

Quickie Revolution							
Front Seat Height (For Standard Frames Only)							
CASTER FORK	5" MICRO	5" POLY	6" PNEU	6" POLY			
MOUNTING HOLE	5 1/4"	5 1/4"	5 1/4"	5 1/4"			
TOP	N/A	17 1/4"	17 1/8"	17 3/4"			
BOTTOM	17 1/4"	18 1/4"	18 7/8"	18 3/4"			
REAR SEAT HEIGHT (FOR STANDARD FRAMES AND AXLE POSITION ONLY)							
CAMBER/REAR	22"	24"	24"	24" HT	26"	26"	26" HT
AXLE HOLE	PNEU	PNEU	LATEX	TEN/BBALL	HP	PNEU	TEN/BBALL
POSITION			TURBO TRD	CLINCHER		TURBO TRD	
0°, 4°/NO SPACER	15 1/4"	16 1/8"	15 3/4"	16 1/4"	16 7/8"	17 1/4"	17 5/8"
8°, 12°/NO SPACER	14 1/2"	15 3/8"	15"	15 1/2"	16 1/8"	16 1/2"	16 7/8"
8°, 4° + 3/4°	16"	16 7/8"	16 1/2"	17"	17 5/8"	18"	18 3/8"
8°, 12° + 3/4°	15 1/4"	16 1/8"	15 3/4"	16 1/4"	16 7/8"	17 1/4"	17 5/8"
0°, 4° + 1 1/2°	16 3/4"	17 5/8"	17 1/4"	17 3/4"	18 3/8"	18 3/4"	19 1/8"
8°, 12° + 1 1/2°	15"	16 7/8"	16 1/2"	17"	17 5/8"	18"	18 3/8"
Overall Width (With 4° of camber)							
Frame Width	14"	15"	16"	17"	18"	19"	20"
Spoke/Wheel*							
with Handrims	24"	25"	26"	27"	28"	29"	30"
Quickie M6							
Weight:	23 lbs.						
Weight Capacity:	250 lbs						

Quickie 2 True Fit Program

Physical changes can cause a chair not to fit properly. To ensure proper fit, Sunrise offers one-time replacement for growth of the seat sling, back upholstery, footplates, hangers, back tubes, frame and crossbrace at no extra charge. (For up to three years from delivery on all Quickie 2 family wheelchairs shipped after November 15, 1995.)

ACTUAL MEASUREMENTS MAY VARY ACCORDING TO CHAIR CONFIGURATION AND ACCESSORIES. PLEASE CONSULT ORDER FORM FOR COMBINATIONS AND RESTRICTIONS. SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE.

LIFETIME FRAME AND AXLE WARRANTY

We highly recommend distribution through a certified Rehab Professional when fitting and ordering any Quickie Products.

