

The Quickie Recliner has all of the features of a Quickie 2 plus a reclining backrest. The Quickie 2 can be converted to a recliner and back again. The ease of use and effectiveness of this chair make it a big hit with caregivers.

**Quickie Recliner Highlights:**

- Reclines 90°-160°
- Ease of use ensured by the dual trigger release with pin lock on the stabilizer bar
- 10" removable headrest extension



The Quickie Recliner has all of the features of a Quickie 2 plus a reclining backrest. The Quickie 2 can be converted to a recliner and back again. The ease of use and effectiveness of this chair make it a big hit with caregivers.

**Quickie Recliner Highlights:**

- Reclines 90°-160°
- Ease of use ensured by the dual trigger release with pin lock on the stabilizer bar
- 10" removable headrest extension



FOLDING FAMILY

Quickie Folding Chairs – A complete folding line to suit your mobility needs.

Ultra-Light



Lightweight



Positioning



Pediatrics



Quickie 2 & Quickie 2 Kids & Quickie 2 HD

Front Seat Height (For Standard Frames Only)

CASTER FORK MOUNTING HOLE	5" LP POLY 5 1/4" FORK	5" LP POLY 7" FORK	6" PNEU 6" FORK	6" POLY 6" FORK	8" PN 5 1/4" FORK
Top	16 15/16"	18 3/4"	18 5/16"	18 3/16"	18 1/2"
Bottom	17 15/16"	19 3/4"	19 5/16"	19 3/16"	19 1/2"

Rear Seat Height (For Standard Frames and Standard Axle Plate Position Only)

REAR AXLE/HOLE POSITION	PNEU	FP	PNEU	LP	FP	HP	PNEU/ KEVLAR	HP CLNCHR
1 AND 2	18 7/8"	19 1/2"	19 3/4"	20"	20 1/4"	20 3/8"	20 1/2"	20 1/4"
2 AND 3	17 3/4"	18 3/8"	18 5/8"	18 7/8"	19 1/8"	19 1/4"	19 3/8"	20 1/8"
3 AND 4	16 5/8"	17 1/4"	17 1/2"	17 3/4"	18"	18 1/8"	18 1/4"	19"
4 AND 5	15 1/2"	16 1/8"	16 3/8"	16 5/8"	16 7/8"	17"	17 1/8"	17 7/8"

Overall Width (With Standard Camber and Padded Swing-Away Arm)

FRAME WIDTH	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
SPOKE WHEEL W/HANDRIMS	20 1/2"	21 1/2"	22 1/2"	23 1/2"	24 1/2"	25 1/2"	26 1/2"	27 1/2"	28 1/2"	29 1/2"	30 1/2"

Quickie 2	Quickie 2 HD	Quickie 2 Kids
Weight: 27 lbs		
Weight Capacity: 250 lbs	Weight Capacity: 350 lbs	Weight Capacity: 250 lbs

Quickie 2HP

Front Seat Height (For Standard Frames Only)

CASTER FORK MOUNTING HOLE	3"MICROS 5 1/4" FORK	3"MICROS 6" FORK	3"MICROS 7" FORK	5"LP POLY 5 1/4" FORK	5"LP POLY 6" FORK	5"LP PC 7" FORK
TOP	16 5/16"	17 1/8"	18 1/8"	17 3/16"	18"	19"
BOTTOM	17 5/16"	18 1/8"	19 1/8"	18 3/16"	19"	20"

Rear Seat Height (For Standard Frames and Standard Axle Plate Position Only)

REAR AXLE/HOLE POSITION	PNEU	LATEX	HP	PNEU	TEN/B-BALL	LATEX	PNEU
1 AND 2	19 3/4"	20 1/4"	20 3/8"	20 1/2"	20 5/8"	21 1/4"	21 5/8"
2 AND 3	18 5/8"	19 1/8"	19 1/4"	19 3/8"	19 1/2"	20 1/8"	20 1/2"
3 AND 4	17 1/2"	18"	18 1/8"	18 1/4"	18 3/8"	19"	19 3/8"
4 AND 5	16 3/8"	16 7/8"	17"	17 1/8"	17 1/4"	17 1/8"	18 1/4"

Overall Width (With Standard Camber and Padded Swing-Away Arm)

FRAME WIDTH	12"	13"	14"	15"	16"	17"	18"	19"	20"
SPOKE/MAG WHEEL W/HANDRIMS	21 1/2"	22 1/2"	23 1/2"	24 1/2"	25 1/2"	26 1/2"	27 1/2"	28 1/2"	29 1/2"

Quickie 2 HP
Weight: 22 lbs
Weight Capacity: 250 lbs

Quickie M6

FRAME WIDTH:	24"	26"	28"	30"
OVERALL WIDTH:	32.5"	34.5"	36.5"	38.5"
FOLDED WIDTH:	14"			

Quickie M6

Weight:	53 lbs. w/legrests (26 x 20)
Weight Capacity:	650 lbs

Quickie Recliner

Front Seat Height

CASTER FORK MOUNTING HOLE	5" LP POLY 6" FORK	6" PNEU 6" FORK	6" POLY 6" FORK	8" PNEU 5 1/4" FORK	8" POLY 5 1/4" FORK
TOP	17 3/4"	18 5/16"	18 3/16"	N/A	N/A
BOTTOM	18"	19 5/16"	19 3/16"	19 1/2"	19 5/16"

Rear Seat Height

REAR AXLE HOLE POSITION	22" FP Poly	24" FP Poly	24" HP Clincher	24" Pneu/ Kevlar
Top	15"	15 3/4"	15 7/8"	16"
Middle	16"	16 3/4"	16 7/8"	17"
Bottom	17"	17 3/4"	17 7/8"	18"

Armrest Height

Height-Adjustable

9" - 14"

Footrest Height

SWING-AWAY	ELR
14" - 21"	16 1/2" - 19"

Overall Width (Without Armrests)

Frame Width	14"	15"	16"	17"	18"	19"	20"	Folded
SPOKE/MAG WHEEL W/HANDRIMS & W/HANDRIMS & QUICK-RELEASE AXLES	22 1/4"	23 1/4"	24 1/4"	25 1/4"	26 1/4"	27 1/4"	28 1/4"	12 1/2"
SPOKE/MAG WHEEL W/HANDRIMS & THREADED AXLES	22"	23"	24"	25"	26"	27"	28"	12 1/2"

Quickie Recliner

Weight:	42 lbs. (w/ swing-away footrests, height-adjustable armrests, headrest extension and mag. wheels)
Weight Capacity:	250 lbs

Quickie LX & LXI

Front Seat Height

CASTER FORK MOUNTING HOLE	5"LP Poly 4 3/4" Fork	5"LP Poly 5 1/4" Fork	6"Poly 6" Fork	6"Pneu 6" Fork	8"Poly 5 1/4" Fork	8"Pneu 5 1/4" Fork	8"x2" Caster 8"x2" Fork
TOP	N/A	N/A	17"	17"	N/A	N/A	N/A
BOTTOM	16"	16 1/2"	18"	18"	18"	18"	N/A

Rear Seat Height

REAR AXLE HOLE POSITION (LX ONLY)	20" LP POLY	20" PNEU	22" LP POLY	22" PNEU	24" LP POLY	24" PNEU
TOP	15"	15 1/2"	16"	16 1/2"	17"	17 1/2"
BOTTOM	16"	16 1/2"	17"	17 1/2"	18"	18 1/2"

REAR AXLE HOLE POSITION (ADJ.AXLE)	20" LP POLY	20" PNEU	22" LP POLY	22" PNEU	24" LP POLY	24" PNEU
TOP	N/A	15"	15 1/2"	16"	16"	17"
MIDDLE	15 1/2"	16"	16 1/2"	17 1/2"	17 1/2"	18"
BOTTOM	17"	17 1/2"	18"	18 1/2"	18 1/2"	19 1/2"

Armrest Height

Footrest Height

FLIP-BACK HEIGHT ADJ.	FLIP-BACK/ PADDED S/A HEIGHT ADJ.	HT. ADJ.	SWING-AWAY	ELR/ALR
9 1/4"	8 3/4"-12 5/8"	8 3/4"	7 1/4"-12 1/4"	12 3/4"-15 1/4"
			8 3/4"-13 3/4"	14"-16 1/2"
				15 1/4"-17 3/4"
				16 1/2"-19"
				17 3/4"-20 1/4"

Overall Width (Without Armrests)

FRAME WIDTH	12"	14"	16"	18"	20"	FOLDED
MAG WHEEL W/HANDRIMS	19"	21"	23"	25"	27"	12 1/2"
SPOKEWHEEL W/HANDRIMS	20"	22"	24"	26"	28"	12 1/2"

Quickie Recliner

Weight:	28 lbs
Weight Capacity:	250 lbs