



Aqua Creek Products

THE LEADERS IN RECREATION, FITNESS AND ABILITY!


Item# F-25ABLE

Bathroom Accessibility

The Aquatic Bathtub Lift Elite or ABLE is a wall mounted bath tub lift. The ABLE provides stability for those who need assistance to use their bathtubs.

The manually pivoting seat provides complete control on entering the tub and is thin enough to allow deep submersion into the water.

Site construction and product assembly are required.



FEATURES:

- 300 lb (136 kg) weight capacity
- Durable powder coat finish
- Stainless steel construction
- Rechargeable battery with wall mount charger
- Water resistant handset
- Folding seat with seatbelt
- Requires little to no assistance to operate
- Maximum tub height 22"
- Minimum ceiling height 8'
- Five year structural warranty
- Five year pro-rated electronics warranty

PROUDLY SOLD BY: