



Item #	Seat	Color	UOM
SFSCOUTDLX3	16.5" Folding Seat	Red/Blue ■ ■	1/ea



Spitfire Scout 3 - DLX

Compact Travel Scooter

- Quick and easy disassembly.
- Glossy finish interchangeable color shrouds in Red and Blue.
- Armrests are padded and adjustable.
- Quick connect batteries.
- Ergonomic throttle control.
- Height-adjustable swivel seat with fold-down backrest.
- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Delta tiller.
- Easy freewheel operation

specifications

Weight Capacity	300 lbs.	Charger	2A Offboard
Top Speed	4.25 mph	Brakes	Electromagnetic
Maximum Range	9 miles	Front Wheels	8" x 2" Flat-free
Turning Radius	51.2"	Rear Wheels	8" x 2" Flat-free
Climbing Angle	6°	Base Weight	49 lbs.
Ground Clearance	2.5"	Battery Weight (pair)	20 lbs.
Dimensions	42.5" (L) x 20.5" (W)	Seat Weight	16 lbs.
Floor-to-Seat Height	17" - 19"	Total Weight	85 lbs.
Seat Dimensions	16.5" (W) x 16" (D)	Heaviest Piece Weight (Front Section)	29 lbs.
Controller	CTE 45A	Warranty on Frame	Lifetime
Motor	24V x 270V x 4700 rpm	Warranty on Electronic Controller/Drive Train Components	14 months
Batteries (pair)	12V x 12AH	Warranty on Batteries	6 months



Item #	Seat	Color	UOM
SFSCOUTDLX4	16.5" Folding Seat	Red/Blue	1/cs



A

Spitfire Scout 4 - DLX

Compact Travel Scooter

- Quick and easy disassembly. (Figure A)
- Glossy finish interchangeable color shrouds in Red and Blue.
- Armrests are padded and adjustable.
- Quick connect batteries.
- Ergonomic throttle control.
- Height-adjustable swivel seat with fold-down backrest.
- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Delta tiller.
- Easy freewheel operation

specifications			
Weight Capacity	300 lbs.	Charger	2A Offboard
Top Speed	4 mph	Brakes	Electromagnetic
Maximum Range	9 miles	Front Wheels	8" x 2" Flat-free
Turning Radius	53.75"	Rear Wheels	8" x 2" Flat-free
Climbing Angle	6°	Base Weight	58 lbs.
Ground Clearance	2.5"	Battery Weight (pair)	20 lbs.
Dimensions	42.5" (L) x 20.5" (W)	Seat Weight	16 lbs.
Floor-to-Seat Height	17" - 19"	Total Weight	94 lbs.
Seat Dimensions	16.5" (W) x 16" (D)	Heaviest Piece Weight (Font Section)	39 lbs.
Controller	CTE 45A	Warranty on Frame	Lifetime
Motor	24V x 270W x 4700 rpm	Warranty on Electronic Controller/Drive Train Components	14 months
Batteries (pair)	12V x 12AH	Warranty on Batteries	6 months